

# Flash Chromatography



373 - 382

## Flash Chromatography

### CLARICEP™ Irregular Flash Media

Irregular CS Silica Columns.....	375
Irregular CM Silica Columns.....	377

### CLARICEP Spherical Flash Media

Spherical Silica Columns.....	378
Spherical C18 Columns.....	379
Spherical AQ C18 Columns.....	380

## Seamlessly Upgrade from Traditional Columns to CLARICEP Flash

Bonna-Agela Technologies have developed a technology that effectively deactivates the silica surface. As a result, CLARICEP Flash columns have less surface activity than ordinary silica columns and demonstrate significantly improved chromatographic performance.

### Traditional Column

- High surface activity that causes instability of certain compounds
- Unwanted tailing or overly long retention of basic compounds due to secondary ionic reactions or metal chelating effects
- Poor reproducibility
- Limited selectivity range
- Pressure limited

Vs.

### CLARICEP Column

- Deactivated silica surface promotes compound stability
- Excellent peak shape and performance for both acidic and basic compounds
- High quality and reproducibility
- Wide range of selectivities
- High pressure tolerance



# Flash Chromatography

## CLARICEP™ Irregular CS Silica Columns


- Ultra pure silica packing
- Additional acid and deionized water wash
- Narrow particle size range
- Carefully controlled water content

### Technical Specifications

Surface Area:	480 m <sup>2</sup> /g
Surface pH:	6.3-7.2
Water Content:	3.0-5.0%
Average Particle Size:	40-60 μm
Average Pore Size:	60 Å

## Improved Peak Shape and Faster Analysis

### Aniline Peak Symmetry and Retention Test



#### Flash Conditions:

**Column:** Claricep Irregular Silica CS (40-60 μm, 60 Å, 40 g)  
**Brand I:** Flash Irregular Silica (40 g)  
**Mobile Phase:** Dichloromethane/ Methanol (99 : 1)  
**Flow Rate:** 20 mL/min  
**Detector:** UV @ 254 nm  
**Temperature:** Ambient  
**Retention Time:** CLARICEP CS: 4.090 min  
Brand I: 4.373 min  
**Sample:** Aniline


Comparative separations may not be representative of all applications.

# Flash Chromatography

## CLARICEP™ Irregular CS Silica Columns (cont'd)


### The Purity of Sesamol in Sesame Oil

#### HPLC Analysis



**Column:** Fully Porous, 5 µm, C18 Column  
**Dimensions:** 4.6 x 150 mm  
**Mobile Phase:** Methanol/Water (75:25)


#### Flash Purification



#### Flash Conditions:

**Column:** Claricep™ Irregular Silica CS (40-60 µm, 60 Å, 12 g)  
**Part No.:** [CS140012-0](#)  
**Mobile Phase:** Acetic ether/ Petroleum ether (12:88)  
**Flow Rate:** 18 mL/min  
**Injection Volume:** 4 mL  
**Sample Concentration:** 400 mg/20 mL  
**Instrument:** CHEETAH™ MP 100

#### Purity Confirmation



# Flash Chromatography

## CLARICEP™ Irregular CM Silica Columns

- Significantly improved performance over regular flash columns
- Silica deactivated by proprietary process
- Alternative selectivity for complex purification requirements


### Technical Specifications

Average Particle Size:	40-60 µm
Average Pore Size:	60 Å
Surface Area:	480 m <sup>2</sup> /g

## Better Peak Shape With CM Silica

### HPLC Test:

Unmodified and Deactivated Silica were packed into individual stainless steel columns (4.6 x 150 mm) and then evaluated on a HPLC System



### HPLC Conditions:

**Dimensions:** 150 x 4.6 mm  
**Mobile Phase:** Dichloromethane/Methanol (98:2)  
**Flow Rate:** 1.8 mL/min  
**Injection Volume:** 5 µL  
**Temperature:** 30 °C  
**Detector:** UV @ 254 nm  
**Sample:** Catechol 100 µg/mL

# Flash Chromatography

## CLARICEP™ Spherical Silica Columns

### Technical Specifications

Average Particle Size:	20 µm	20-35 µm	
Surface Area:	320 m <sup>2</sup> /g	480 m <sup>2</sup> /g	320 m <sup>2</sup> /g
Water Content:	3.0 - 5.0%	3.0 - 5.0%	3.0 - 5.0%
Average Pore Size:	100 Å	60 Å	100 Å

- Higher Resolution
- Better Purification
- Lower Backpressure
- Faster Flow Rate
- Higher Loading Capacity
- Lower Backpressure
- Faster Flow Rate

### Purification of a Sample with Methacrylic Acid Ester Target Compound




#### Sample Information:

The sample is colorless liquid, with about 60% target compound by weight

Dissolve 0.2 mL of sample into 1.5 mL ethanol sonication

#### Structure:

Small molecular weight with UV absorption of methacrylic acid ester R: no UV absorption



#### Flash Conditions:

**Column A:** Claricep Spherical Silica (20-35 µm, 100 Å, 12 g, 2 columns in tandem)

**Column B:** Claricep Spherical Silica (20 µm, 100 Å, 12 g, 2 columns in tandem)

**Part Nos.:** Column A: [SS130012-0](#)  
Column B: [SS120012-0](#)

**Mobile Phase:** A: Hexane B: Ethanol

Gradient:	Time/min	B %
	0	5
	20	55

**Flow Rate:** 12 mL/min

**Detector:** UV @ 254/220 nm

**Sample Loading:** 0.2 mL

Claricep Flash silica 20 µm is a better choice for complex sample polarity. It provides higher resolution and better purification performance.

# Flash Chromatography


## CLARICEP™ Spherical C18 Columns

- Significantly improved performance over regular flash columns
- Silica deactivated by proprietary process
- Alternative selectivity for complex purification requirements

### Technical Specifications

Average Particle Size:	20 µm	20 - 35 µm		40 - 60 µm
Average Pore Size:	100Å	60Å	100Å	100Å
Carbon Loading:	14%	15%	14%	14%

### Tetrandrine Extracted from Natural Products, Formulation



#### Flash Conditions:

**Column A:** Brand X Flash Irregular C18 (40-60 µm 100Å, 12 g, 3 columns in tandem)

**Column B:** Claricep Spherical C18 (20-35 µm 100Å, 12 g, 3 columns in tandem)

**Column C:** Claricep Spherical C18 (20 µm 100Å, 12 g, 3 columns in tandem)

**Part Nos.:** Column B: [S0230012-0](#)

Column C: [C-S0120012-0](#)

**Mobile Phase:** A: Water

B: Methanol with 0.06% diethylamine

Gradient:	Time/min	% B
	0	85
	100	85

**Detector:** UV @ 254/282 nm

**Sample:** Tetrandrine

### Did You Know?



**Flash Chromatography** also known as medium pressure chromatography is:

- A pressure driven hybrid for medium and short column chromatography optimized for rapid separation
- Popularized years ago by Clark Still of Columbia University
- An alternative to slow and inefficient gravity-fed chromatography

Comparative separations may not be representative of all applications.

# Flash Chromatography


## CLARICEP™ Spherical AQ C18 Columns

- Greater polar retention under reversed phase
- Applicable for both hydrophilic and hydrophobic compounds
- Useful for mixtures of compounds with varying polarities

### Technical Specifications

Surface Area:	300 m <sup>2</sup> /g	320 m <sup>2</sup> /g	320 m <sup>2</sup> /g
Average Particle Size:	40-60 μm	20-35 μm	20 μm
Average Pore Size:	100 Å	100 Å	100 Å
Carbon Loading:	14%	15%	15%

### High Resolution Separation of Iridoids



#### Flash Conditions:

**Column:** Claricep Spherical AQ C18 (20-35 μm, 100 Å)  
**Mobile Phase:** Methanol/ Water/Formic Acid  
**Flow Rate:** 26 mL/min  
**Detector:** UV @ 231/214 nm  
**Sample:** Iridoid Compounds





# Flash Chromatography


## CLARICEP™ Screw-on Flash Columns

This first series of Flash Screw-on is a new feature of Claricep columns that allows the user to load solid samples directly on the column.

### i-Series

The i-Series features a managed column head space with a secured screw-on lid. This new design allows either loading of liquid samples directly onto the column head or loading of impregnated solid sample directly into the space. Users will benefit from:

- Choice of loading method based on sample properties
- Narrow band for liquid samples because of wide loading area
- Dry-loading of solid impregnated samples minimizes band broadening
- Customized loading method upon user preference
- Higher load tolerance and reusable consumables




### s-Series

The s-Series columns are fully packed without a head void and a secure screw-on lid. In addition, the Luer lock fittings for both inlet and outlet allows easy operation of tandem columns or the coupling of a loading cartridge.

Compatible with the following instrument:

- Biotage




### c-Series

The c-Series shares the same design, but the column outlet does not have a luer lock structure with a screw-on lid, this simplifies tubing connection across various flash systems.

Compatible with instruments from the following:

- Teledyne ISCO
- Agela Technologies
- Buchi
- Grace
- Interchim
- Other brands



## Easy to Order Screw-On Flash Columns

### Ordering Information Screw-On Flash Columns

Part No.	Description
<b>For i-series</b>	
SN	Part Number starts with "SN". As an example, to order <a href="#">CS140012-0</a> in i-series, Part Number to order is <a href="#">SN-CS140012-0</a>
<b>For s-series</b>	
S	Part Number starts with "S". As an example, to order <a href="#">CS140012-0</a> in s-series, Part Number to order is <a href="#">S-CS140012-0</a>
<b>For c-series</b>	
C	Part Number starts with "C". As an example, to order <a href="#">CS140012-0</a> in c-series, Part Number to order is <a href="#">S-CS140012-0</a>



Note: Request an external option to load any dry sample loading via empty tubes (20/40/80/120g tube equivalent size) by contacting a Phenomenex representative.

# Flash Chromatography

## Claricep™ Irregular and Spherical Silica and Spherical Bonded Phase Flash Columns

### Ordering Information

#### Irregular Silica Phase

Type	Average Particle Size (µm)	Pore Size (µm)	Part No.	Silica Amount (g)	Unit (pk)
Silica (CS) Standard Silica	40 - 60 µm	60 Å	<a href="#">CS140004-0</a>	4	20
			<a href="#">CS140012-0</a>	12	20
			<a href="#">CS140020-0</a>	20	20
			<a href="#">CS140040-0</a>	40	10
			<a href="#">CS140080-0</a>	80	5
			<a href="#">CS140120-0</a>	120	5
			<a href="#">CS140330-0</a>	330	1
			<a href="#">CS140800-0</a>	800	1
			<a href="#">CS1401500-0</a>	1500	1

Also available in i-series, s-series & c-series (12g, 20g, 40g, 80g, & 120g) per request (Contact your Sales Rep)

#### Spherical Silica Phase

Type	Average Particle Size (µm)	Pore Size (µm)	Part No.	Silica Amount (g)	Unit (pk)
Spherical Silica	20 - 35 µm	60 Å	<a href="#">SS130004-0</a>	4	20
			<a href="#">SS130012-0</a>	12	20
			<a href="#">SS130020-0</a>	20	20
			<a href="#">SS130040-0</a>	40	10
			<a href="#">SS130080-0</a>	80	5
			<a href="#">SS130120-0</a>	120	5
			<a href="#">SS130330-0</a>	330	1

Also available in i-series, s-series & c-series (12g, 20g, 40g, 80g, & 120g) per request (Contact your Sales Rep)

#### Spherical Bonded Phase

Type	Average Particle Size (µm)	Pore Size (µm)	Part No.	Silica Amount (g)	Unit (pk)
C18	40 - 60 µm	100 Å	<a href="#">S0240004-0</a>	4	20
			<a href="#">S0240012-0</a>	12	20
			<a href="#">S0240020-0</a>	20	20
			<a href="#">S0240040-0</a>	40	10
			<a href="#">S0240080-0</a>	80	5
			<a href="#">S0240120-0</a>	120	5
			<a href="#">S0240330-0</a>	330	1
			<a href="#">S0240800-0</a>	800	1
			<a href="#">S02401500-0</a>	1500	1

Also available in i-series, s-series & c-series (12g, 20g, 40g, 80g, & 120g) per request (Contact your Sales Rep)

#### Spherical Bonded Phase

Type	Average Particle Size (µm)	Pore Size (µm)	Part No.	Silica Amount (g)	Unit (pk)
C18	20 - 35 µm	100 Å	<a href="#">S0230004-0</a>	4	20
			<a href="#">S0230012-0</a>	12	20
			<a href="#">S0230020-0</a>	20	20
			<a href="#">S0230040-0</a>	40	10
			<a href="#">S0230080-0</a>	80	5
			<a href="#">S0230120-0</a>	120	5
			<a href="#">S0230330-0</a>	330	1

Also available in i-series, s-series & c-series (12g, 20g, 40g, 80g, & 120g) per request (Contact your Sales Rep)

#### Spherical Bonded Phase

Type	Average Particle Size (µm)	Pore Size (µm)	Part No.	Silica Amount (g)	Unit (pk)
AQ C18	20 - 35 µm	100 Å	<a href="#">S0230004-0</a>	4	20
			<a href="#">S0230012-0</a>	12	20
			<a href="#">S0230020-0</a>	20	20
			<a href="#">S0230040-0</a>	40	10
			<a href="#">S0230080-0</a>	80	5
			<a href="#">S0230120-0</a>	120	5
			<a href="#">S0230330-0</a>	330	1

Also available in i-series, s-series & c-series (12g, 20g, 40g, 80g, & 120g) per request (Contact your Sales Rep)



Reach out to your technical sales representative to ask about our CLARICEP flash bulk options.



NOTE: Additional CLARICEP Flash Column Phases, Formats, and Sizes Available. Contact Your Sales Rep or Our Technical Support.