

## GC Column Selection Guidelines

The Master Resolution Equation.....	88
Choosing Your Selectivity.....	88-90
Choosing Your Dimensions.....	91-92
Column Selection Charts.....	93-99

## Zebtron™ GC Columns

Meet Your GC Column Family.....	101
Zebtron Unlimited.....	102-131
Zebtron PLUS.....	132-143
Zebtron Inferno™.....	144-151
Zebtron Essentials.....	152-162
Zebtron Guard Columns.....	163-164

## GC Accessories

Inlet Supplies.....	165-180
Column Unions & Splitters.....	181
Gas Management.....	182-183
Gas & Moisture Traps.....	184
Tools & Maintenance Kits.....	184
GC Column Check Mixes.....	185

“The chromatography quality and performance are excellent [with Zebtron]. Column bleed is minimal at 320 °C. Peak quality remains good for 5 to 6 months averaging 40 injections in a 24 hour period, 6 to 7 days per week.”

**Kevin Walkup**  
Specialized Assays, Inc.


The opinions stated herein are solely those of the speaker and not necessarily those of any company or organization.

## The Basic Principle of Column Selection

How do you choose a column? Do you reach into a cabinet of mystery columns, look to your favorite 5% phenyl phase, or borrow one from a colleague? Understanding how column parameters impact key elements of the master resolution equation will help you quickly make the right column selection for successful separations.

$$R_s = \left[ \frac{\sqrt{N}}{4} \right] \times \left[ \frac{\alpha - 1}{\alpha} \right] \times \left[ \frac{k}{k + 1} \right]$$

	Efficiency Term	Selectivity Term	Retention Term
<b>Relates to:</b>	Column Length Column ID	Column Phase	Column ID Film Thickness
<b>Other Considerations:</b>	Carrier Gas Linear Velocity	Temperature	Temperature



## Selectivity Has the Biggest Impact on Resolution


Resolution between two analytes is mainly determined by the selectivity of the stationary phase. By increasing the resolution between two compounds, the total analysis time can often be reduced significantly!

### Selectivity vs. Polarity

Polarity gives a general guideline for sample capacity and separation, which can affect peak shape and resolution. However, two columns may have similar polarity but show different separation profiles due to dissimilar phase chemistries. For example, ZB-35 and ZB-1701 are close in polarity, but the cyanopropyl group makes ZB-1701 very different from ZB-35 in terms of selectivity.


#### ZB-35

Polarity: 18



#### ZB-1701

Polarity: 19



# Choosing Your Selectivity (cont'd)

## The 3 Most Prevalent GC Interactions

The following selection guidelines can be a starting point for choosing Zebron™ columns in common selectivities. Please contact your Phenomenex representative for additional assistance.

### Dispersive Forces (Van der Waals Interactions)

- Weakest of all intermolecular forces and occurs between non-polar compounds
- Separation is based on boiling point (classic example – hydrocarbon separation in SimDist analysis)




### Dipole-Dipole Interactions

- Either permanently present or induced by analyte-stationary phase interactions
- Higher dipole-dipole interaction can help separate compounds with similar boiling points, but different chemical structures

### Hydrogen Bonding (Acid-Base Interactions)


- Can cause poor peak shape or irreversible binding to the inlet liner or to the column itself
- Zebron columns are specially deactivated to minimize these interactions

## Choosing A “5” Phase




# Choosing Your Selectivity *(cont'd)*

## Choosing A "1" Phase



## Choosing A "PEG" Phase






# Choosing Your Dimensions

## Length

Longer columns can improve resolution, but they will also increase run times. Under isothermal conditions, doubling column length only increases resolution by 41 %, but doubles the run time! Choose a column length that balances efficiency with acceptable run times.

Short	Good Starting Length	Long
<b>15 m or less</b>	<b>30 m</b>	<b>60 m or more</b>
<b>Applications</b> <ul style="list-style-type: none"><li>• High boilers</li><li>• GC-MS applications</li></ul> <b>Advantages</b> <ul style="list-style-type: none"><li>• Faster run times</li><li>• Higher temp. limits</li><li>• Lower bleed</li><li>• Higher efficiency</li></ul> <b>Disadvantages</b> <ul style="list-style-type: none"><li>• Less inert</li><li>• Limited retention</li></ul>		<b>Applications</b> <ul style="list-style-type: none"><li>• Complex samples with closely eluting peaks</li><li>• Low boilers</li><li>• Less active samples</li><li>• Complex temperature ramps</li></ul> <b>Advantages</b> <ul style="list-style-type: none"><li>• Better resolution</li></ul> <b>Disadvantages</b> <ul style="list-style-type: none"><li>• Slow run times</li></ul>

## Try The GC Column Finder!

Easily select a column by part number, manufacturer, industry, application, or official method **in under 1 minute.**




[www.phenomenex.com/FindGC](http://www.phenomenex.com/FindGC)

# Choosing Your Dimensions *(cont'd)*


## Internal Diameter

Column internal diameter (ID) has a major impact on both resolution and sample capacity. Unlike column length, using smaller ID columns can actually lead to faster run times, because the column length required with a small ID is often shorter due to increased efficiency.

Narrow	Good Starting ID	Wide
0.10, 0.18, 0.20 mm	0.25 mm	0.32, 0.53 mm
<p><b>Applications</b></p> <ul style="list-style-type: none"> <li>• Complex samples</li> </ul> <p><b>Advantages</b></p> <ul style="list-style-type: none"> <li>• Faster run times</li> <li>• Better resolution</li> </ul> <p><b>Disadvantages</b></p> <ul style="list-style-type: none"> <li>• Lower sample capacity</li> <li>• Easily overloaded</li> </ul>		<p><b>Applications</b></p> <ul style="list-style-type: none"> <li>• Dirty samples</li> <li>• Highly concentrated samples</li> </ul> <p><b>Advantages</b></p> <ul style="list-style-type: none"> <li>• Increased sample capacity</li> <li>• Good for on-column injections</li> </ul> <p><b>Disadvantages</b></p> <ul style="list-style-type: none"> <li>• Decreased efficiency</li> <li>• May need higher flow rates</li> <li>• Not compatible with most GC-MS</li> </ul>

## Film Thickness

Film thickness determines solute retention and plays an important role in column sample capacity. Thin film columns are faster and provide higher resolution, but lower sample capacity. In most instances, choose the thinnest film possible that still provides adequate retention. When working with active samples, using a slightly thicker film can significantly improve peak shape.

Thin	Good Starting Film	Thick
0.10, 0.18 $\mu$ m	0.25 $\mu$ m	0.50 $\mu$ m or more
<p><b>Applications</b></p> <ul style="list-style-type: none"> <li>• High boilers</li> <li>• GC-MS applications</li> </ul> <p><b>Advantages</b></p> <ul style="list-style-type: none"> <li>• Faster run times</li> <li>• Higher temp. limits</li> <li>• Lower bleed</li> <li>• Higher efficiency</li> </ul> <p><b>Disadvantages</b></p> <ul style="list-style-type: none"> <li>• Less inert</li> <li>• Limited retention</li> </ul>		<p><b>Applications</b></p> <ul style="list-style-type: none"> <li>• Low boilers</li> <li>• Gases, solvents, purgeables, volatiles</li> <li>• Purity testing</li> </ul> <p><b>Advantages</b></p> <ul style="list-style-type: none"> <li>• Better inertness</li> <li>• Higher capacity</li> </ul> <p><b>Disadvantages</b></p> <ul style="list-style-type: none"> <li>• Slow run times</li> <li>• Lower temp. limits</li> <li>• Higher bleed</li> </ul>

# Cross-Reference by Manufacturer

## Upgrade to Zebron!

Our commitment to quality and innovation is what makes Zebron GC columns well-suited for any application. Performance is GUARANTEED.

Zebron Phase	Zebron Composition	Restek®	Agilent®	Supelco®	SGE®	OV
ZB-1	100% Dimethylpolysiloxane	Rtx®-1, Rtx-1PONA, Rtx-1 F&F	DB®-1, <a href="#">DB-2887</a> , DB-1 EVDX, HP-1, HP-101, HP-PONA,Ultra 1, CP-Sil 5 CB	SPB®-1, SPB-1 TG, SE-30, MET-1, SPB-1 Sulfur, SPB-HAP	BP1, BP1-PONA, BPX1-SimD	OV-1
ZB-DHA-PONA	100% Dimethylpolysiloxane	Rtx-DHA	HP-PONA, DB-PETRO, CP-Sil PONA CB	Petrocol®-DH		
ZB-1 PLUS™	100% Dimethylpolysiloxane	Rtx-1ms, Rxi®-1ms	DB-1ms, DB-1ms Ultra Inert, HP-1ms, HP-1ms Ultra Inert, CP-Sil 5 CB MS, VF-1ms	MDN-1, Equity®-1	SolGel-1ms™	
ZB-1HT Inferno™	100% Dimethylpolysiloxane	Rxi-1HT	DB-1ht, CP-SimDist	Petrocol 2887		
ZB-1XT SimDist	100% Dimethylpolysiloxane	MXT®-1HT SimDist, MXT-1, MXT-1 SimDist, <a href="#">MXT-2887</a>	CP-SimDist UltiMetal, CP-Sil 8 CB UltiMetal, BPX1-SimD, DB-PS1, DB-HT SimDis, <a href="#">DB-PS2887</a>			
ZB-5	5% Phenyl 95% Dimethylpolysiloxane	Rtx-5	DB-5, HP-5, Ultra 2, HP-PAS-5, CP-Sil 8 CB, HP-5ms Ultra Inert	MDN-5, SPB-5, PTE-5, SE-54, PTA-5, Equity-5, Sac-5	BP5, BPX5	OV-5
ZB-5 PLUS™	5% Phenyl 95% Dimethylpolysiloxane	Rtx-5ms, Rxi-5ms, Rtx-5Amine	DB-5, HP-5ms, HP-5msi	MDN-5S		
ZB-5HT Inferno	5% Phenyl 95% Dimethylpolysiloxane	Rxi-5HT, Rtx-5HT Stx®-5HT, XTI®-5HT	DB-5ht, VF-5ht	HT-5		
ZB-5ms	5% Phenyl-Arylene 95% Dimethylpolysiloxane	Rtx-5Sil MS, Rxi-5Sil MS	DB-5ms, DB-5.625, VF-5ms, DB-5ms EVDX, CP-Sil 8 CB MS			
ZB-5MS PLUS™	5% Phenyl-Arylene 95% Dimethylpolysiloxane	Rxi-5Sil MS	DB-5ms Ultra Inert, HP-5ms Ultra Inert, DB-5ms, VF-5ms	SLB®-5ms		
ZB-SemiVolatiles	5% Phenyl-Arylene 95% Dimethylpolysiloxane	Rxi-5Sil MS, Rxi-5ms	DB-5ms Ultra Inert	SLB-5ms		
ZB-35	35% Phenyl 65% Dimethylpolysiloxane	Rtx-35, Rtx-35ms	DB-35, DB-35ms, HP-35, HP-35ms	MDN-35, SPB-35, SPB-608	BPX35, BPX608	OV-11
ZB-35HT Inferno	35% Phenyl 65% Dimethylpolysiloxane			Phenomenex Exclusive		
ZB-50	50% Phenyl 50% Dimethylpolysiloxane	Rtx-50	DB-17, DB-17HT, DB-17ms, HP-50+, DB-17 EVDX, CP-Sil 24 CB	<a href="#">SP-2250</a> , SPB-17, SPB-50	BPX50	OV-17
ZB-624	6% Cyanopropylphenyl 94% Dimethylpolysiloxane	Rtx-1301, Rtx-624	<a href="#">DB-1301</a> , DB-624, DB-VRX, HP-VOC, <a href="#">CP-1301</a> , CP-Select 624 CB	<a href="#">SPB-1301</a> , SPB-624	BP624	OV-624
ZB-624 PLUS™	Proprietary	Rxi-624Sil MS	CP-Select 624 CB, DB-624UI Ultra Inert			
<a href="#">ZB-1701</a>	14% Cyanopropylphenyl 86% Dimethylpolysiloxane	Rtx-1701	<a href="#">DB-1701</a> , CP-Sil 19 CB	<a href="#">SPB-1701</a> , Equity-1701	BP10	<a href="#">OV-1701</a>
<a href="#">ZB-1701P</a>	14% Cyanopropylphenyl 86% Dimethylpolysiloxane		<a href="#">DB-1701P</a>			
ZB-FAME	High Cyanopropyl		CP-Sil 88, HP-88, DB-23	SP®-2560, <a href="#">SP-2380</a>		
ZB-WAX	Polyethylene Glycol	Rtx-WAX, Famewax, Stabilwax-DB	DB-WAXetr, HP-INNOWax, CP-Wax 57 CB	MET-Wax, Omegawax	SolGel-WAX™	
ZB-WAX PLUS™	Polyethylene Glycol	Stabilwax®	DB-WAX, CAM, HP-20M, Carbowax 20M, CP-Wax 52 CB	SUPELCOWAX® 10	BP20	Carbowax 20M
ZB-FFAP	Nitroterephthalic Acid Modified Polyethylene Glycol	Stabilwax-DA	DB-FFAP, HP-FFAP, CP-FFAP CB CP-Wax 58 FFAP CB	Nukol, <a href="#">SPB-1000</a>	BP21	OV-351
ZB-MultiResidue™-1	Proprietary	Rtx-CLPesticides, Stx-CLPesticides				
ZB-MultiResidue-2	Proprietary	Rtx-CLPesticides2, Stx-CLPesticides2				
ZB-CLPesticides-1	Proprietary	Rtx-CLPesticides, Stx-CLPesticides				
ZB-CLPesticides-2	Proprietary	Rtx-CLPesticides2, Stx-CLPesticides2				
ZB-XLB	Proprietary	Rtx-XLB, Rxi-XLB	DB-XLB, VF-XMS	MDN-12		
ZB-XLB-HT Inferno	Proprietary			Phenomenex Exclusive		
ZB-Drug-1	Proprietary			Phenomenex Exclusive		
ZB-BAC-1	Proprietary	Rtx-BAC1	<a href="#">DB-ALC1</a>			
ZB-BAC-2	Proprietary	Rtx-BAC2	<a href="#">DB-ALC2</a>			
ZB-Bioethanol	Proprietary			Phenomenex Exclusive		
ZB-PAH-EU	Proprietary	Rxi-PAH	DB-PAH EU			
ZB-PAH-CT	Proprietary		PAH-Select			
ZB-Dioxin	Proprietary	Rtx-Dioxin 2	DB-Dioxin, DB-225, DB-5MSUI	<a href="#">SP-2330</a>		

This section is, neither in terms of manufacturers nor in terms of their products, a complete list, and the accuracy of the data is not guaranteed. Small differences in dimensions or performance might be possible and slight adjustments to your application may be necessary.

# Environmental Selection Chart

Listed below are recommended Zebron columns for environmental and EPA methods. Other columns may also be used for these analyses — please contact Phenomenex for your specific GC column needs.


Drinking Water	Method #	Description	Primary Column	Confirmation Column	Page	
	501.3	Trihalomethanes by GC-MS with Selected Ion Monitoring (SIM)	ZB-624, ZB-624 <sup>PLUS</sup> <sup>™</sup>		157, 140	
	502.2	Volatile Halogenated Organics by Purge & Trap GC/PID/ELCD	ZB-624, ZB-624 <sup>PLUS</sup>		157, 140	
	503.1	Volatile Aromatics and Unsaturated Organics by Purge & Trap GC	ZB-624, ZB-624 <sup>PLUS</sup>		157, 140	
	504.1	1,2-Dibromoethane (EDB), 1,2-Dibromo-3-chloropropane (DBCP), and 1,2,3-Trichloropropane (123TCP) by GC	ZB-CLPesticides-1 ZB-MultiResidue <sup>™</sup> -1	ZB-CLPesticides-2 ZB-MultiResidue-2	116 118	
	505	Organohalide Pesticides & Aroclors by GC-ECD	ZB-CLPesticides-1 ZB-MultiResidue-1	ZB-CLPesticides-2 ZB-MultiResidue-2	116 118	
	507	Nitrogen & Phosphorus Containing Pesticides by GC-NPD	ZB-MultiResidue-1 ZB-CLPesticides-2	ZB-MultiResidue-2 ZB-CLPesticides-2	118 116	
	508	Chlorinated Pesticides by GC-ECD	ZB-CLPesticides-1 ZB-MultiResidue-1	ZB-CLPesticides-2 ZB-MultiResidue-2	116 118	
	509	Ethylene Thiourea (ETU) by GC-NPD	ZB-WAX <sup>PLUS</sup> <sup>™</sup>	ZB-1701	138, 158	
	513	2, 3, 7, 8-Tetrachlorodibenzo-p-dioxin by GC-HRMS	ZB-SemiVolatiles		114	
	515.3	Chlorinated Acids by Liquid-Liquid Extraction, Derivatization and GC-ECD	ZB-XLB	ZB-35	162, 155	
	521	Nitrosamines by Solid Phase Extraction (SPE) and GC-MS/MS with Large Volume Injection	ZB-SemiVolatiles		114	
	522	1,4-Dioxane by Solid Phase Extraction (SPE) and GC-MS with Selected Ion Monitoring (SIM)	ZB-SemiVolatiles		114	
	523	Triazine Pesticides and their Degradates by GC-MS	ZB-50		156	
	524.3	Purgeable Organic Compounds by GC-MS	ZB-624, ZB-624 <sup>PLUS</sup>		157, 140	
	525.2	Semi-volatile Organic Chemicals by Solid Phase Extraction (SPE) and GC-MS	ZB-SemiVolatiles		114	
	526	Selected Semi-volatile Organic Compounds by Solid Phase Extraction (SPE) and GC-MS	ZB-SemiVolatiles		114	
	527	Selected Pesticides and Flame Retardants by Solid Phase Extraction (SPE) and GC-MS	ZB-5 <sup>PLUS</sup> <sup>™</sup>		134	
	528	Phenols by Solid Phase Extraction (SPE) and GC-MS	ZB-SemiVolatiles	ZB-35	114, 155	
	529	Explosives and Related Compounds by Solid Phase Extraction (SPE) and GC-MS	ZB-5 <sup>PLUS</sup> <sup>™</sup>		134	
	548	Endothall by Aqueous Derivatization, Liquid-Solid Extraction, and GC-ECD	ZB-SemiVolatiles	ZB-35	114, 155	
	551.1	Chlorinated Solvents & Disinfection Byproducts by Liquid-Liquid Extraction and GC-ECD	ZB-35		155	
	552.3	Haloacetic Acids and Dalapon by Liquid-Liquid Extraction, Derivatization, and GC-ECD	ZB-CLPesticides-1 ZB-XLB	ZB-CLPesticides-2 ZB-35	116, 116 162, 155	
	556	Carbonyl Compounds by Pentafluorobenzylhydroxylamine Derivatization and GC-ECD	ZB-SemiVolatiles	ZB-1701	114, 158	
	Waste Water	Method #	Description	Primary Column	Confirmation Column	Page
		601	Purgeable Halocarbons by Purge & Trap GC	ZB-624, ZB-624 <sup>PLUS</sup>		157, 140
		602	Purgeable Aromatics by Purge & Trap GC	ZB-624, ZB-624 <sup>PLUS</sup>		157, 140
		603	Acrolein & Acrylonitrile Purge & Trap GC	ZB-624, ZB-624 <sup>PLUS</sup>		157, 140
604		Phenols by GC-ECD	ZB-SemiVolatiles		114	
606		Phthalate Esters by GC-ECD	ZB-5 <sup>PLUS</sup> <sup>™</sup>		134	
607		Nitrosamines by GC-NPD	ZB-SemiVolatiles		114	
608		Organochlorine Pesticides and PCBs by GC-ECD	ZB-MultiResidue-1	ZB-MultiResidue-2	118, 118	
609		Nitroaromatics & Isophorone by GC-FID and GC-ECD	ZB-SemiVolatiles		114	
610		Polynuclear Aromatic Hydrocarbons by GC-FID	ZB-PAH-EU ZB-PAH-CT		106 110	
611		Haloethers by GC-ECD	ZB-SemiVolatiles	ZB-SemiVolatiles	114	
612		Chlorinated Hydrocarbons by GC-ECD	ZB-SemiVolatiles		114	
613		2,3,7,8-Tetrachlorodibenzo-p-dioxin by GC-MS	ZB-SemiVolatiles		114	
615		Chlorinated Herbicides by GC-ECD	ZB-CLPesticides-1 ZB-XLB	ZB-CLPesticides-2 ZB-35	116, 116 162, 155	
619		Triazine Herbicides by GC-MS	ZB-50		156	
622		Organophosphorus Pesticides by GC-MS	ZB-MultiResidue-1		118	
624		Purgeable Volatiles by Purge & Trap GC-MS	ZB-624		157	
625		Base/Neutral and Acids by GC-MS	ZB-SemiVolatiles		114	
1613		Tetra- through Octa-Chlorinated Dioxins & Furans by Isotope Dilution HRGC/HRMS	ZB-Dioxin	ZB-SemiVolatiles	102, 114	
1614		Polybrominated Diphenyl Esters (PBDEs) by HRGC/HRMS	ZB-5HT Inferno <sup>™</sup> ZB-SemiVolatiles		146 114	
1618		Organohalide Pesticides, Organophosphorus Pesticides, and Phenoxy-Acid Herbicides by GC	ZB-MultiResidue-1	ZB-MultiResidue-2	118, 118	
1624		Volatile Organic Compounds by Isotope Dilution GC-MS	ZB-624, ZB-624 <sup>PLUS</sup>		157, 140	
1625		Semi-volatile Organic Compounds by Isotope Dilution GC-MS	ZB-SemiVolatiles		114	
1653		Chlorinated Phenols by In-Situ Acetylation and GC-MS	ZB-SemiVolatiles		114	
1657		Organophosphorous Pesticides by GC-FPD	ZB-MultiResidue-1	ZB-MultiResidue-2	118, 118	
1658		Phenoxy-Acid Herbicides by GC-ECD	ZB-MultiResidue-1	ZB-MultiResidue-2	118, 118	
1659		Dazomet by GC-NPD	ZB-MultiResidue-1	ZB-MultiResidue-2	118, 118	
1666		Pharmaceutical Volatile Organic Compounds by Purge & Trap GC or Isotope Dilution GC-MS	ZB-SemiVolatiles (Direct Injection) ZB-624 (Purge & Trap), ZB-624 <sup>PLUS</sup>		114 157 140	
1668	Polychlorinated Biphenyl (PCB) Congeners by HRGC/HRMS	ZB-MultiResidue-1	ZB-1	118, 152		
1671	Pharmaceutical Manufacturing Volatile Organic Compounds by GC-FID	ZB-1		152		
7850	White Phosphorus (P4) by Solvent Extraction and GC-NPD	ZB-1		152		



# Environmental Selection Chart


Listed below are recommended Zebron columns for environmental and EPA methods. Other columns may also be used for these analyses – please contact Phenomenex for your specific GC column needs.

Solid Waste	Method #	Description	Primary Column	Confirmation Column	Page
	8010B	Halogenated Volatile Organics by GC-ELCD	ZB-624, ZB-624 <sup>PLUS</sup> <sup>™</sup>		157, 140
	8015C	Nonhalogenated Organics by GC	ZB-5HT Inferno <sup>™</sup>		146
	8020A	Aromatic Volatile Organics by GC-PID	ZB-WAX, ZB-WAX <sup>PLUS</sup> <sup>™</sup>		160 138
	8021B	Aromatic and Halogenated Volatiles by GC-PID or GC-ELCD	ZB-624, ZB-624 <sup>PLUS</sup>	ZB-1 (thick phase)	157, 140, 152
	8030A	Acrolein and Acrylonitrile by GC-FID	ZB-624, ZB-624 <sup>PLUS</sup>		157, 140
	8032A	Acrylamide by GC-ECD	ZB-5HT Inferno		146
	8041	Phenols by GC-ECD or GC-FID	ZB-SemiVolatiles		114
	8061A	Phthalate Esters by GC-ECD	ZB-SemiVolatiles	ZB-1701	114, 158
	8081B	Organochlorine Pesticides by GC-ECD	ZB-MultiResidue <sup>™</sup> -1 ZB-CLPesticides-1	ZB-MultiResidue-2 ZB-CLPesticides-2	118 116
	8082A	Polychlorinated Biphenyls (PCBs) by GC-ECD	ZB-MultiResidue-1 ZB-CLPesticides-1	ZB-MultiResidue-2 ZB-CLPesticides-2	118 116
	8091	Nitroaromatics and Cyclic Ketones by GC-ECD or GC-NPD	ZB-SemiVolatiles	ZB-1701	114, 158
	8095	Explosives by GC-ECD	ZB-50		156
	8100	Polynuclear Aromatic Hydrocarbons by GC-FID	ZB-SemiVolatiles ZB-35		114 155
	8121	Chlorinated Hydrocarbons by GC-ECD	ZB-MultiResidue-1	ZB-MultiResidue-2	118
	8131	Aniline and Selected Derivatives by GC-NPD	ZB-SemiVolatiles	ZB-1	114, 152
	8141B	Organophosphorus Pesticides by GC-FPD or GC-NPD	ZB-MultiResidue-1 ZB-CLPesticides-1	ZB-MultiResidue-2 ZB-CLPesticides-2	118 116
	8151A	Chlorinated Herbicides by GC-ECD	ZB-CLPesticides-1 ZB-XLB	ZB-CLPesticides-2 ZB-35	116, 116 162, 155
	8260B	Volatile Organic Compounds by GC-MS	ZB-624, ZB-624 <sup>PLUS</sup>		157, 140
	8270D	Semi-volatile Organic Compounds by GC-MS	ZB-SemiVolatiles		114
	8272	Polynuclear Aromatic Hydrocarbons (PAHs) by SPME and GC-MS with Selected Ion Monitoring (SIM)	ZB-SemiVolatiles, ZB-35		114 155
	8280B	Polychlorinated Dibenzo-P-Dioxins (PCDDs) and Polychlorinated Dibenzofurans (PCDFs) By HRGC-LRMS	ZB-SemiVolatiles		114
	8290A	Polychlorinated Dibenzo-P-Dioxins (PCDDs) and Polychlorinated Dibenzofurans (PCDFs) By HRGC/HRMS	ZB-SemiVolatiles		114
	8410	Semi-Volatile Organic Compounds by GC-FTIR	ZB-SemiVolatiles		114
	8430	Bis(2-chloroethyl) Ether and Hydrolysis Products by Direct Aqueous Injection GC-FTIR	ZB-WAX <sup>PLUS</sup>		138

Air	Method #	Description	Primary Column	Page
	TO-1	Volatile Organic Compounds by Thermal Adsorption and GC-MS	ZB-1 <sup>PLUS</sup> <sup>™</sup>	132
	TO-2	Volatile Organic Compounds by Carbon Molecular Sieve Adsorption and GC-MS	ZB-1 <sup>PLUS</sup>	132
	TO-3	Volatile Organic Compounds by Cryogenic Preconcentration Techniques and GC-FID /ECD	ZB-1 <sup>PLUS</sup>	132
	TO-4A	Pesticides and Polychlorinated Biphenyls (PCBs) by High Volume Polyurethane Foam (PUF) Sampling and GC	ZB-MultiResidue-1	118
	TO-7	N-Nitrosodimethylamine by GC-MS	ZB-WAX <sup>PLUS</sup>	138
	TO-9A	Polychlorinated, Polybrominated, and Brominated/Chlorinated Dibenzo-p-Dioxins and Dibenzofurans by HRGC/HRMS	ZB-SemiVolatiles	114
	TO-10A	Pesticides and Polychlorinated Biphenyls (PCBs) by Low Volume Polyurethane Foam (PUF) Sampling and GC	ZB-MultiResidue-1	118
	TO-13A	Polycyclic Aromatic Hydrocarbons (PAHs) by GC-MS	ZB-SemiVolatiles	114
	TO-14A	Volatile Organic Compounds by Specially Prepared Canisters and GC	ZB-1 <sup>PLUS</sup>	132
	TO-15	Volatile Organic Compounds by Specially Prepared Canisters and GC-MS	ZB-1 <sup>PLUS</sup>	132

# Food & Flavors Selection Chart

Listed below are recommended Zebron columns for food safety, food quality, and flavor/fragrance methods. Other columns may also be used for these analyses — please contact Phenomenex for your specific GC column needs.

Food Safety	Compound Class	Analysis	Recommended Columns	Page
	<b>Pesticides &amp; Antimicrobials</b>	Multi-Residue Pesticide Screening	ZB-MultiResidue™-1 and -2	118
		Organochlorine Pesticides in Water	ZB-MultiResidue-1 and -2	118
		Organochlorine Pesticides in Foods of Plant Origin	ZB-MultiResidue-1 and -2	118
		Organophosphorus Pesticides in Foods of Plant Origin	ZB-MultiResidue-1 and -2	118
		Triazine Pesticides in Water	ZB-50	156
		Triazine Pesticides in Foods of Plant Origin	ZB-50	156
		Chloramphenicol in Foods of Animal Origin	ZB-1 <sup>PLUS</sup> ™	132
		<b>Environmental Contaminants</b>	Polybrominated Diphenyl Ethers (PBDEs) in Food	ZB-5MS <sup>PLUS</sup> ™, ZB-SemiVolatiles, ZB-35
Polychlorinated Biphenyls (PCBs) in Water	ZB-MultiResidue-1, ZB-XLB-HT Inferno™		118, 150	
Polychlorinated Dibenzo-dioxins (PCDDs) in Food	ZB-5MS <sup>PLUS</sup> , ZB-SemiVolatiles		136, 114	
Polychlorinated Dibenzo-furans (PCDFs) in Food	ZB-5MS <sup>PLUS</sup> , ZB-SemiVolatiles		136, 114	
Polycyclic Aromatic Hydrocarbons (PAHs) in Water	ZB-PAH-EU, ZB-PAH-CT, ZB-5MS <sup>PLUS</sup> , ZB-SemiVolatiles, ZB-35		106, 110, 136, 114, 155	
<b>Food Contact Materials</b>	Food Packaging Volatiles	ZB-624, ZB-624 <sup>PLUS</sup> ™	157, 140	
	Dioxins and Furans in Food	ZB-Dioxin	102	
	Melamine in Food	ZB-XLB-HT Inferno	150	
	Cyanuric Acid in Food	ZB-XLB-HT Inferno	150	
	Phthalates in Food	ZB-5MS <sup>PLUS</sup>	136	
	Residual Solvents in Food	ZB-624, ZB-624 <sup>PLUS</sup> , ZB-WAX <sup>PLUS</sup>	157, 140, 138	
	Bisphenol A & F (BPA/BPF) in Food	ZB-5MS <sup>PLUS</sup>	136	
<b>Additives &amp; Preservatives</b>	Parabens in Food	ZB-5MS <sup>PLUS</sup>	136	
	Chloropropanols (3-MCPD) in Food	ZB-5MS <sup>PLUS</sup>	136	
	Flavor Additives (Borneol)	ZB-MultiResidue-1	118	
	Phenolic Antioxidants (BHA & BHT) in Food	ZB-50	156	
	Tocopherols in Food	ZB-5MS <sup>PLUS</sup>	136	
<b>Process Contaminants</b>	Acrylamide in Foods	ZB-5HT Inferno	146	
	Acrylamide, Acrylonitrile, and Acrolein in Water	ZB-624, ZB-624 <sup>PLUS</sup>	157, 140	
	Benzene in Food	ZB-WAX <sup>PLUS</sup>	138	
	Glycols in Food	ZB-WAX <sup>PLUS</sup>	138	
<b>Hormones</b>	Steroid Hormones in Food	ZB-5MS <sup>PLUS</sup> , ZB-1 <sup>PLUS</sup>	136, 132	

## Try The GC Column Finder!


Easily select a column by part number, manufacturer, industry, application, or official method **in under 1 minute.**




[www.phenomenex.com/FindGC](http://www.phenomenex.com/FindGC)

# Food & Flavors Selection Chart


Listed below are recommended Zebron columns for food safety, food quality, and flavor/fragrance methods. Other columns may also be used for these analyses — please contact Phenomenex for your specific GC column needs.


Food Quality	Compound Class	Analysis	Recommended Columns	Page
	<b>Fatty Acids &amp; FAMES</b>	Food Industry Fatty Acid Methyl Esters (FAMES)	ZB-FAME	112
		Marine Oil Fatty Acid Methyl Esters (FAMES)	ZB-FAME	112
		Saw Palmetto Fatty Acid Methyl Esters (FAMES)	ZB-FAME	112
		Free Fatty Acids	ZB-FFAP	161
		Essential Fatty Acids (EFAs) Omega-3 and Omega-6	ZB-FAME	112
	<b>Triglycerides</b>	Butter, Canola Oil, Olive Oil, and Peanut Oil Triglycerides	ZB-5HT Inferno™	146
	<b>Alcoholic Beverages</b>	Cognac Compounds	ZB-WAX <sup>PLUS</sup> ™	138
		Distilled Liquor Screen	ZB-FFAP	161
		Ethanol in Beer	ZB-Bioethanol	122
		Sulfur in Beer	ZB-1 <sup>PLUS</sup> ™	132
		Whiskey Compounds	ZB-WAX <sup>PLUS</sup>	138
		Wine Compounds	ZB-WAX, ZB-WAX <sup>PLUS</sup>	160, 138
	<b>Other Acids</b>	Organic Acids	ZB-FFAP	160
		Amino Acids	ZB-50	156
	<b>Sterols</b>	Sterols in Lard, Margarine, Peanut Butter, or Olive Oil	ZB-5HT Inferno	146
<b>Sugars</b>	Alditol Acetates	ZB-5MS <sup>PLUS</sup> ™	136	
	Trimethylsilyl (TMS) Sugars	ZB-MultiResidue™-1	118	

Flavors & Fragrances	Compound Class	Analysis	Recommended Columns	Page
	<b>Essential Oils</b>	Cold-Pressed Orange Oil	ZB-WAX <sup>PLUS</sup>	138
		Ginkgo Biloba Oil, Lavender Oil, and Ylang Ylang Oil	ZB-1 <sup>PLUS</sup>	132
		Peppermint Oil	ZB-WAX	160
		Rose Oil	ZB-XLB	162
		Spearmint Oil	ZB-5MS <sup>PLUS</sup>	136
	<b>Flavors</b>	Flavors Screening	ZB-FFAP	161
		Flavor Allergens	ZB-5MS <sup>PLUS</sup>	136
		Flavor Volatiles	ZB-1 <sup>PLUS</sup> , ZB-WAX <sup>PLUS</sup> , ZB-624	132, 138, 157
		Alcoholic Beverage Profile	ZB-FFAP	161
		Honey Profile	ZB-WAX <sup>PLUS</sup>	138
	<b>Fragrances</b>	Fragrance Screening	ZB-WAX <sup>PLUS</sup> , ZB-624	138, 157
		Fragrance Allergens	ZB-1 <sup>PLUS</sup>	132

# Pharmaceutical Selection Chart

Listed below are recommended Zebron columns for USP and pharmaceutical methods. Other columns may also be used for these analyses – please contact Phenomenex for your specific GC column needs.

USP	Phase Composition	Recommended Columns	Page
	G1 Dimethylpolysiloxane Oil	ZB-1, ZB-1 <sup>PLUS</sup> ™, ZB-1HT Inferno™	152, 132, 144
	G2 Dimethylpolysiloxane Gum	ZB-1, ZB-1 <sup>PLUS</sup> , ZB-1HT Inferno	152, 132, 144
	G3 50 % Phenyl 50 % Methylpolysiloxane	ZB-50	156
	G5 Not less than 70 % of 3-Cyanopropylpolysiloxane	ZB-FAME	112
	G8 80 % Bis (3-Cyanopropyl-20 % 3-Cyanopropylphenylpolysiloxane)	ZB-FAME	112
	G9 Methylvinylpolysiloxane	ZB-1 <sup>PLUS</sup> , ZB-1HT Inferno, ZB-1	132, 144, 152
	G14 Polyethylene Glycol (Average MW 950-1,050)	ZB-WAX, ZB-WAX <sup>PLUS</sup> ™	160, 138
	G15 Polyethylene Glycol (Average MW 3,000-3,700)	ZB-WAX, ZB-WAX <sup>PLUS</sup>	160, 138
	G16 Polyethylene Glycol (Average MW 15,000)	ZB-WAX, ZB-WAX <sup>PLUS</sup>	160, 138
	G17 75 % Phenyl 25 % Methylpolysiloxane	ZB-50	156
	G20 Polyethylene Glycol (Average MW of 380-420)	ZB-WAX, ZB-WAX <sup>PLUS</sup>	160, 138
	G25 Polyethylene Glycol TPA (Carbowax 20M Terephthalic Acid)	ZB-FFAP	161
	G27 5 % Phenyl 95 % Methylpolysiloxane	ZB-5, ZB-5 <sup>PLUS</sup> ™, ZB-5HT Inferno	153, 134, 146
	5 % Phenyl-Arylene 95 % Methylpolysiloxane	ZB-5ms, ZB-5MS <sup>PLUS</sup> ™, ZB-SemiVolatiles	154, 136, 114
	G28 25 % Phenyl 75 % Methylpolysiloxane	ZB-35, ZB-35HT Inferno	155, 148
	G32 20 % Phenylmethyl 80 % Dimethylpolysiloxane	ZB-35, ZB-35HT Inferno	155, 148
	G35 Polyethylene Glycol & Diepoxide Esterified with Nitrotterephthalic Acid	ZB-FFAP	161
	G36 1 % Vinyl 5 % Phenylmethylpolysiloxane	ZB-5, ZB-5 <sup>PLUS</sup> , ZB-5HT Inferno	153, 134, 146
	G38 Phase G1 Plus A Tailing Inhibitor	ZB-1, ZB-1 <sup>PLUS</sup> , ZB-1HT Inferno	152, 132, 144
	G39 Polyethylene Glycol (Average MW 1,500)	ZB-WAX, ZB-WAX <sup>PLUS</sup> ™	160, 138
	G41 Phenylmethyldimethylsilicone (10 % Phenyl Substituted)	ZB-5, ZB-5 <sup>PLUS</sup> , ZB-5HT Inferno	153, 134, 146
	G42 35 % Phenyl 65 % Dimethylpolysiloxane	ZB-35, ZB-35HT Inferno	155, 148
	G43 6 % Cyanopropylphenyl 94 % Dimethylpolysiloxane	ZB-624, ZB-624 <sup>PLUS</sup> ™	157, 140
	G46 14 % Cyanopropylphenyl 86 % Methylpolysiloxane	ZB-1701, ZB-1701P	158, 159
	G47 Polyethylene glycol (average MW 8,000)	ZB-WAX <sup>PLUS</sup> , ZB-WAX	138, 160
	G48 Highly polar, partially cross-linked cyanopolysiloxane	ZB-FAME	112
	G51 50 % Phenyl - 50 % Dimethylpolysiloxane phase which is modified with an aromatic selector for optimized separation of polycyclic aromatic hydrocarbons (PAH)	ZB-PAH-CT, ZB-PAH-EU	110,106

Residual Solvents	USP <467> Procedure	USP Phase for Residual Solvents	Recommended Columns	Page
	Procedure A	G43 (6 % Cyanopropyl 94 % Dimethylpolysiloxane)	ZB-624, ZB-624 <sup>PLUS</sup>	157, 140
	Procedure B	G16 (Polyethylene Glycol)	ZB-WAX <sup>PLUS</sup>	138
	Procedure C	G43 or G16	ZB-624 <sup>PLUS</sup> or ZB-WAX <sup>PLUS</sup>	140, 138



Doing Headspace Testing?


Find the right headspace vial for your analysis and learn more about Verex™ Certified Vial Products with our interactive web tool.

[www.phenomenex.com/verex](http://www.phenomenex.com/verex)



# ASTM Method Selection Chart

Listed below are recommended Zebron columns for ASTM methods. Other columns may also be used for these analyses — please contact Phenomenex for your specific GC column needs.

ASTM	Method	Description	Recommended Columns	Page
	D 1946	Reformed gas	ZB-1	152
	D 2268	Analysis of n-heptane and iso-octane (high purity)	ZB-1	152
	D 2306-96	Xylene isomers	ZB-WAX, ZB-WAX <sub>PLUS</sub> <sup>™</sup>	160, 138
	D 2426	Butadiene and styrene in butadiene concentrates	ZB-1	152
	D 2504	Non-condensable gases in C1-C3 hydrocarbons	ZB-1 (thick phase)	152
	D 2580	Phenols in water	ZB-WAX <sub>PLUS</sub>	138
	D 2600	Aromatic traces in light saturated hydrocarbons	ZB-WAX	160
	D 2804	Purity of methyl ethyl ketone	ZB-WAX	160
	D 2887	SimDist analysis of petroleum fractions	ZB-1, ZB-1XT SimDist	152, 124
	D 2908	Volatile organics in water	ZB-WAX, ZB-WAX <sub>PLUS</sub>	160, 138
	D 2998	Polyhydric alcohols in alkyd resins	ZB-1	152
	D 2999	Monopentaerythritol in commercial pentaerythritol	ZB-1	152
	D 3009	Composition of turpentine	ZB-WAX <sub>PLUS</sub>	138
	D 3054	Purity and benzene content of cyclohexane	ZB-1	152
	D 3086	Organochlorine pesticides in water	ZB-CLPesticides-1 or -2, ZB-MultiResidue <sup>™</sup> -1 or -2	116, 118
	D 3168	Polymers in emulsion paints	ZB-1	152
	D 3271	Solvent analysis in paints	ZB-WAX <sub>PLUS</sub>	138
	D 3304	PCBs in environmental materials	ZB-MultiResidue-1 or -2	118
	D 3328	Comparison of waterborne petroleum oils	ZB-1	152
	D 3329	Purity of methyl isobutyl ketone	ZB-WAX <sub>PLUS</sub>	138
	D 3432	Toluene diisocyanates in urethane prepolymers	ZB-1	152
	D 3447	Purity of trichlorotrifluoroethane (CFC-113)	ZB-1, ZB-624	152, 157
	D 3452	Identification of rubber	ZB-1HT Inferno <sup>™</sup>	144
	D 3465	Purity of monomeric plasticizers	ZB-1	152
	D 3524	Diesel fuel in lubricating oil (SAE 30)	ZB-1HT Inferno	144
	D 3534	PCBs in water	ZB-5, ZB-5 <sub>PLUS</sub> <sup>™</sup>	153, 134
	D 3606	Benzene and toluene in gasoline	ZB-1	152
	D 3687	Volatile organic compounds	ZB-WAX, ZB-WAX <sub>PLUS</sub>	160, 138
	D 3710	Gasoline fractions	ZB-1XT SimDist	124
	D 3725	Fatty acids in drying oils	ZB-FFAP	161
	D 3760	Analysis of cumene	ZB-WAX, ZB-WAX <sub>PLUS</sub>	160, 138
	D 3797	Analysis of o-xylene	ZB-WAX, ZB-WAX <sub>PLUS</sub>	160, 138
	D 3798	Analysis of p-xylene impurities	ZB-WAX, ZB-WAX <sub>PLUS</sub>	160, 138
	D 3876	Methoxyl and hydroxypropyl substitution in cellulose ether products	ZB-1	152
	D 3962	Impurities in styrene	ZB-FFAP	161
	D 4059	PCBs in insulating liquids	ZB-5 <sub>PLUS</sub> , ZB-5HT Inferno	134, 146
	D 4275	Butylated hydroxy toluene in ethylene and ethylenevinylacetate polymers	ZB-1	152
	D 4367	Benzene in hydrocarbon solvent	ZB-1	152
	D 4420	Aromatics in gasoline	ZB-1	152
	D 4735	Thiophene impurities in benzene	ZB-FFAP	161
	D 4768	Phenol and cresol inhibitors in insulating oils	ZB-FFAP	161
	D 5060	Impurities in ethylbenzene	ZB-FFAP, ZB-WAX, ZB-WAX <sub>PLUS</sub>	161, 160, 138
	D 5134	Petroleum naphthas through n-nonane	ZB-1, ZB-DHA-PONA	152, 126
	D 5135-95	Analysis of styrene	ZB-WAX, ZB-WAX <sub>PLUS</sub>	160, 138
	D 5441	Analysis of Methyl Tert-Butyl Ether (MTBE)	ZB-DHA-PONA	126
	D 5501	Determination of denatured bioethanol	ZB-1, ZB-Bioethanol, ZB-DHA-PONA	152, 122, 126
	D 5580	Aromatics in finished gasoline	ZB-1	152
	D 6352	Extended SimDist	ZB-1HT Inferno, ZB-1XT SimDist	144, 124
	D 6584	Determination of glycerine in biodiesel	ZB-5HT Inferno	146
	D 6729-30, D 6733	Components in spark ignition fuels	ZB-DHA-PONA	126
D 7169	Crude Oil; Vacuum distillates	ZB-1XT SimDist	124	
E 0202	Analysis of glycols	ZB-WAX <sub>PLUS</sub> , ZB-1	138, 152	
E 1100	Analysis of denatured ethanol	ZB-WAX <sub>PLUS</sub> , ZB-Bioethanol	138, 122	



## Essentials

A collection of tried-and-true selectivities, Essentials phases are the smart starting point for the GC method developer.

ZB-1	ZB-624
ZB-5	ZB-1701
ZB-5ms	ZB-1701P
ZB-35	ZB-WAX
ZB-50	ZB-FFAP
	ZB-XLB



## PLUS

PLUS phases offer a suite of upgrades compared to their Essentials counterparts – from exceptional inertness to enhanced aqueous stability.

ZB-1 <sup>PLUS</sup> ™
ZB-5 <sup>PLUS</sup> ™
ZB-5MS <sup>PLUS</sup> ™
ZB-WAX <sup>PLUS</sup> ™
ZB-624 <sup>PLUS</sup> ™



## Inferno™

Resilient under even the most intense GC conditions, Inferno phases dare to defy high boilers, contaminants, and carry-overs.

ZB-1HT
ZB-5HT
ZB-35HT
ZB-XLB-HT



## Unlimited

Designed for the truly bold GC scientist, Unlimited phases unleash the power of selectivity for targeted performance that breaks from the mold.

ZB-PAH-EU
ZB-PAH-CT
ZB-Dioxin
ZB-FAME
ZB-SemiVolatiles
ZB-MultiResidue™ -1 & -2
ZB-CLPesticides -1 & -2
ZB-Drug-1
ZB-BAC-1 & -2
ZB-1XT SimDist
ZB-Bioethanol
ZB-DHA-PONA

# Meet Your GC Column Family



## Selected Zebron Polarities

Polarity	5	<b>ZB-1</b> ZB-DHA-PONA ZB-1PLUS™ ZB-1HT Inferno™ ZB-1XT SimDist	<b>For Non-Polar Analytes</b> <ul style="list-style-type: none"> <li>• Alkanes</li> <li>• Aromatics</li> <li>• Oils</li> <li>• Boiling Point Separations</li> </ul>
	8	<b>ZB-5</b> ZB-5ms ZB-5PLUS™ ZB-5MSPLUS™ ZB-5HT Inferno ZB-SemiVolatiles	
	9	<b>ZB-XLB</b> ZB-XLB-HT Inferno	
	11	<b>ZB-MultiResidue™-1</b>	
	13	<b>ZB-624</b> ZB-624PLUS™	<b>For Slightly Polar Analytes</b> <ul style="list-style-type: none"> <li>• Volatiles</li> <li>• Drugs</li> <li>• Pesticides</li> </ul>
	15	<b>ZB-MultiResidue-2</b>	
	18	<b>ZB-35</b> ZB-35HT Inferno	
	19	<b>ZB-1701</b> ZB-1701P	
	24	<b>ZB-50</b>	
	52	<b>ZB-WAXPLUS™</b>	<b>For Very Polar Analytes</b> <ul style="list-style-type: none"> <li>• Polar Volatiles</li> <li>• Alcohols</li> <li>• Phenols</li> <li>• Acids</li> </ul>
	57	<b>ZB-WAX</b>	
	58	<b>ZB-FFAP</b>	

## Meet Your GC Column Family Zebron Unlimited

### Food Testing

ZB-FAME .....	112
---------------	-----

### Environmental Testing

ZB-Dioxin.....	102
ZB-PAH-EU.....	106
ZB-PAH-CT.....	110
ZB-SemiVolatiles .....	114
ZB-CLPesticides-1 & -2 .....	116
ZB-MultiResidue™-1 & -2 .....	118

### Fuels

ZB-Dioxin.....	102
ZB-PAH-EU.....	106
ZB-PAH-CT.....	110
ZB-Bioethanol .....	122
ZB-1XT SimDist.....	124
ZB-DHA-PONA .....	126

### Forensics & Toxicology

ZB-Drug-1.....	128
ZB-BAC-1 & -2.....	130

### Zebron PLUS

ZB-1PLUS™.....	132
ZB-5PLUS™.....	134
ZB-5MSPLUS™.....	136
ZB-WAXPLUS™.....	138
ZB-624PLUS™.....	140

### Zebron Inferno™

ZB-1HT Inferno.....	144
ZB-5HT Inferno.....	146
ZB-35HT Inferno.....	148
ZB-XLB-HT Inferno .....	150

### Zebron Essentials

ZB-1 .....	152
ZB-5 .....	153
ZB-5ms .....	154
ZB-35 .....	155
ZB-50 .....	156
ZB-624 .....	157
ZB-1701 .....	158
ZB-1701P.....	159
ZB-WAX .....	160
ZB-FFAP .....	161
ZB-XLB .....	162

### Zebron Guard Columns

Guardian™ Integrated Guard Columns .....	163
Z-Guard™ Columns.....	164



## ZB-Dioxin

### Improve Lab Productivity by 50 %

- Fast PCB analysis
- Enhanced resolution of TCDD and TCDF
- Improved column lifetime with integrated guard column option
- MS certified, low bleed GC column

Upgrade to Zebron from traditional phases used for Dioxin analysis:

#### Agilent®

- DB®-5MSUI
- DB-Dioxin
- DB-225


#### Restek®

- Rtx®-Dioxin2

#### Supelco®

- SP®-2330

### Why Choose ZB-Dioxin?



**A Proprietary Phenyl phase** which provides improved resolution of critical dioxin isomers

Intermediate polarity and thin film (60 meter x 0.20 µm x 0.25 mm) to **reduce analysis times** and **increase productivity**

**Extensive cross-linking through ESC™** (Engineered Self Cross-Linking™) for low bleed and high temperature stability to minimize GC-MS maintenance and system downtime.

### ZB-Dioxin for Fast Dioxin and PCB Analysis

Zebron ZB-Dioxin columns are specifically tailored for the analysis of aromatic compounds like dioxins, furans and PCBs in Food and Environmental matrices. The current analysis of tetra-octa dioxin and furan is time consuming and requires two GC columns and two GC-HRMS instruments. By switching to the unique selectivity of ZB-Dioxin, you will gain enhanced resolution of 2,3,7,8-TCDD and 2,3,7,8-TCDF from its isomers in one run, with only a single ZB-Dioxin GC column. In addition, ZB-Dioxin serves as your single column solution for Dioxin and PCB analysis. Upgrade your existing GC column to a Zebron ZB-Dioxin GC column and get all the analytical benefits and productivity gains of a single column solution.

Learn more at:

[www.phenomenex.com/GCDioxin](http://www.phenomenex.com/GCDioxin)

#### Our Customer Says YES!

“

*The ZB-Dioxin achieves superior resolution for both 2,3,7,8-TCDD and 2,3,7,8-TCDF while not only maintaining chromatography for the hexes but actually improving it. This is all performed while reducing the overall runtime over traditional 5ms dioxin columns by as much as 25%. The ZB-Dioxin increases throughput by not only eliminating the need for a second column confirmation, but also by allowing additional samples to be analyzed in each 12-hour analytical sequence.*

”

**Andrew Patterson, Technical Director  
Eurofins Specialty Services, USA**



The opinions stated herein are solely those of the speaker and not necessarily those of any company or organization.







## ZB-Dioxin (cont'd)



### Increase Throughput and Improve Resolution of Dioxins and PCBs by using Zebron ZB-Dioxin!



Zebron ZB-Dioxin has a unique phase that allows for improved resolution of critical dioxin and PCB congeners and its consistent film thickness allows for optimal dioxin analysis on a single column.



Traditional	VS.	ZB-Dioxin Upgrade
		
<ul style="list-style-type: none"> <li>PCB and Dioxin required 2 different GC columns</li> </ul>		<ul style="list-style-type: none"> <li>Zebron ZB-Dioxin is a single column solution for Dioxin and PCB</li> </ul>

Traditional	VS.	ZB-Dioxin Upgrade
		
<ul style="list-style-type: none"> <li>Higher analysis cost: 2 GC-HRMS + 2 GC columns</li> </ul>		<ul style="list-style-type: none"> <li>Lower analysis cost: 1 GC-HRMS + 1 GC column</li> </ul>

Traditional	VS.	ZB-Dioxin Upgrade
		
<ul style="list-style-type: none"> <li>Long run time for Dioxin analysis</li> <li>First column (5% phenyl phase) ~60 minutes</li> <li>Second column (225 phase) ~30 minutes</li> </ul>		<ul style="list-style-type: none"> <li>Faster run time using one ZB-Dioxin ~40 minutes</li> </ul>

Traditional	VS.	ZB-Dioxin Upgrade
		
<ul style="list-style-type: none"> <li>Shorter column lifetime for difficult matrix like soil</li> </ul>		<ul style="list-style-type: none"> <li>Longer column lifetime with ZB-Dioxin Guardian™ option (Part No: <a href="#">7KG-G045-10-GGA</a>)</li> </ul>

Traditional	VS.	ZB-Dioxin Upgrade
		
<ul style="list-style-type: none"> <li>Lower throughput from customer perspective</li> </ul>		<ul style="list-style-type: none"> <li>HIGH throughput from customer perspective</li> </ul>

Traditional	VS.	ZB-Dioxin Upgrade
		
<ul style="list-style-type: none"> <li>Some GC Dioxin columns do not exceed 290 °C Temperature Limits</li> </ul>		<ul style="list-style-type: none"> <li>Low Bleed GC column temp with 320/340 °C Temperature Limits, this will provide higher sensitivity for later eluters and the ability to bake out major contaminants.</li> </ul>

## ZB-Dioxin (cont'd)


### TCDD on a Zebron ZB-Dioxin and a Popular Brand A

#### Zebron ZB-Dioxin GC Column


Part No. [7KG-G045-10](#)  
60 meter x 0.25 mm x 0.20 µm

#### Brand A Premium 5MS Phase

60 meter x 0.25 mm x 0.25 µm



High Resolution of 2,3,7,8-TCDD by using ZB-Dioxin which exceeds 25% valley EPA-1613 method requirement and provided extended lifetime



2,3,7,8-TCDD is not completely resolved which affects the column lifetime

Sample:	ZB-Dioxin	Brand A
1. 1,3,6,8-TCDD	25.65	23.20
2. 1,2,3,7-TCDD	30.58	30.33
3. 1,2,3,8-TCDD	30.75	30.55
4. 2,3,7,8-TCDD	31.07	30.78
5. 1,2,8,9-TCDD	32.41	32.13

#### Conditions for all separations:

Column 1: Zebron ZB-Dioxin  
 Column 1 Dimension: 60 meter x 0.25 mm x 0.20 µm  
 Column Part No.: [7KG-G045-10](#)  
 Column 2: Brand A Premium 5MS  
 Column 2 Dimension: 60 meter x 0.25 mm x 0.25 µm  
 Guard Column: 5 meter Z-Guard™ Kit  
 Guard Kit Part No.: [7AG-G000-00-GZK](#)  
 Injection: Pulse Splitless (2.0 min, 60 psi) @ 280 °C, 1 µL  
 Liner: Zebron PLUS 4 mm ID Single Taper Liner  
 Liner Part No.: [AG2-0A10-05](#)

Carrier Gas: Helium @ 1.25 mL/min (constant flow)  
 Oven Program: 160 °C for 2.4 min to 200 °C @ 25 °C/min to 220 °C @ 5 °C/min for 19 min to 288 °C @ 4 °C/min to 300 °C @ 5 °C/min for 7.6 min  
 Detector: HRMS  
 Transfer Line Temperature: 300 °C


### TCDF on a Zebron ZB-Dioxin and a Popular Brand A

#### Zebron ZB-Dioxin GC Column


Part No. [7KG-G045-10](#)  
60 meter x 0.25 mm x 0.20 µm

#### Brand A Premium 5MS Phase

60 meter x 0.25 mm x 0.25 µm



Complete resolution of 2,3,7,8-TCDF on a single column ZB-Dioxin—NO NEED FOR ADDITIONAL CONFIRMATION COLUMN




2,3,7,8-TCDF are not completely resolved and need an additional GC column to confirm isomers separation

Sample:	ZB-Dioxin	Brand A
1. 1,3,6,8-TCDF	23.20	23.16
2. 1,3,4,7-TCDF	29.57	29.23
3. 2,3,7,8-TCDF	29.98	29.53
4. 1,2,3,9-TCDF	31.14	30.76
5. 1,2,8,9-TCDF	32.79	32.33

Comparative separations may not be representative of all applications.

## ZB-Dioxin (cont'd)

### Fast GC-MS/MS Analysis of PCBs on a Single 40 Meter Zebron ZB-Dioxin GC Column



#### GC-MS/MS Conditions:

**Column:** Zebron ZB-Dioxin

**Dimension:** 40 meter x 0.18 mm x 0.14 µm

**Part No.:** [7PD-G045-47](#)

**Injection:** Splitless for 1.5 min @ 290 °C, 1 µL

**Recommended Liner:** Zebron PLUS Z-Liner™

(Compatible with Agilent® & Thermo® GC instrument)

**Part No.:** [AG2-0A13-05](#)

**Carrier Gas:** Helium @ 0.8 mL/min (constant flow)

**Oven Program:** 45 °C for 0 min to 175 °C @ 50 °C/min, to 220 °C @ 15 °C/min, to 250 °C @ 5 °C/min for 3 min, to 300 °C @ 50 °C/min for 10 min

**Detector:** GC-MS/MS

**Transfer Line Temperature:** 300 °C

**Mode:** Scan (100-450 m/z)

**Source Temperature:** 300 °C


**Quad Temperature:** 150 °C

**Solvent Delay:** 8.0 min

#### Analyte:

1. PCB-28
2. PCB-31
3. PCB-52
4. PCB-101
5. PCB-81
6. PCB-77
7. PCB-123
8. PCB-118
9. PCB-114
10. PCB-153
11. PCB-105
12. PCB-138
13. PCB-126
14. PCB-167
15. PCB-128
16. PCB-156
17. PCB-157
18. PCB-180
19. PCB-169
20. PCB-189

#### Column Profile



#### Engineered Self Cross-linking™ (ESC)

#### Phase Chemistry

- Proprietary

#### Recommended Applications

- Dioxin and PCB in Food, Environmental Samples
- POPs in Food



Zebron GC Columns MS Certification, see p. 437



Contact Phenomenex or your local Phenomenex distributor for additional GC products and applications.



Extend column lifetime. Add a Z-Guard™ to your next Zebron GC order.

#### Ordering Information

##### Zebron ZB-Dioxin GC Columns

ID (mm)	df (µm)	Temp. Limits °C	Part No.
<b>40-Meter</b>			
0.18	0.14	40 to 320/340	<a href="#">7PD-G045-47</a>
<b>60-Meter</b>			
0.25	0.20	40 to 320/340	<a href="#">7KG-G045-10</a>
<b>60-Meter with 5-Meter Guardian™ Integrated Guard</b>			
0.25	0.20	40 to 320/340	<a href="#">7KG-G045-10-GGA</a>

## ZB-PAH-EU

- Up to 70% faster PAH analysis
- Elevated temperature stability (340/360 °C)
- Great resolution of critical isomers, e.g. Benzo[b,j,k]fluoranthene

Upgrade to Zebron from traditional phases used for PAHs:

**Agilent®**

- DB®-EUPAH


**Restek®**

- Rxi®-PAH

### Priority PAH Analysis by GC

Zebron ZB-PAH-EU columns are designed to move conventional PAH testing to the exceptional, Zebron GC columns come to life through a coupling of innovative spirit and technical excellence. The Zebron ZB-PAH-EU and ZB-PAH-CT offer unparalleled performance through the power of targeted selectivity when analyzing Polycyclic Aromatic Hydrocarbons (PAHs).

#### Analysis of EU 15+1 PAHs



#### GC-MS conditions:

**Column:** Zebron ZB-PAH-EU  
**Dimensions:** 30 meter x 0.25 mm x 0.20 µm  
**Part No.:** 7HG-G043-10  
**Injection:** Split 5:1 @ 330 °C, 1 µL  
**Recommended Liner:** Zebron PLUS Single Taper Z-Liner™  
**Liner Part No.:** AG2-4B13-05 (for Shimadzu® 2010 GC)  
**Carrier Gas:** Helium @ 24 psi (constant pressure)  
**Oven Program:** 45 °C for 0.8 min to 200 °C @ 45 °C/min to 226 °C @ 3 °C/min for 0 min to 320 °C @ 10 °C/min for 20 min  
**Detector:** MSD, 50-500 m/z  
**Transfer Line Temperature:** 300 °C  
**Source Temperature:** 300 °C

#### Sample:

- |                          |                            |
|--------------------------|----------------------------|
| 1. Benzo[c]fluorene      | 9. Benzo[a]pyrene          |
| 2. Benz[aj]anthracene    | 10. Indeno[1,2,3-cd]pyrene |
| 3. Cyclopenta[c,d]pyrene | 11. Dibenzo[a,h]anthracene |
| 4. Chrysene              | 12. Benzo[g,h,i]perylene   |
| 5. 5-Methylchrysene      | 13. Dibenzo[a,l]pyrene     |
| 6. Benzo[b]fluoranthene  | 14. Dibenzo[a,e]pyrene     |
| 7. Benzo[k]fluoranthene  | 15. Dibenzo[a,i]pyrene     |
| 8. Benzo[j]fluoranthene  | 16. Dibenzo[a,h]pyrene     |




## ZB-PAH-EU (cont'd)

### Complete Resolution of EU 15+1 and EPA 610 PAHs

Zebron ZB-EU-PAH GC column demonstrates excellent resolution and accurate quantitation of European regulated EU 15+1 and EPA 610 PAHs.

#### Analysis of EU 15+1 and EPA 610 PAHs



App ID 25713

#### GC-MS conditions:

**Column:** Zebron ZB-PAH-EU  
**Dimensions:** 30 meter x 0.25 mm x 0.20  $\mu$ m  
**Part No.:** [7HG-G043-10](#)  
**Injection:** Split 5:1 @ 330 °C, 1  $\mu$ L  
**Recommended Liner:** Zebron PLUS Single Taper Z-Liner™  
**Liner Part No.:** [AG2-4B13-05](#) (for Shimadzu® 2010 GC)  
**Carrier Gas:** Helium @ 24 psi (constant pressure)  
**Oven Program:** 45 °C for 0.8 min to 200 °C @ 45 °C/min to 226 °C @ 3 °C/min for 0 min to 320 °C @ 10 °C/min for 20 min  
**Detector:** MSD, 50-500 m/z  
**Transfer Line Temperature:** 300 °C  
**Source Temperature:** 300 °C

#### Sample:

- |                   |                           |                            |
|-------------------|---------------------------|----------------------------|
| 1. Naphthalene    | 9. Benzo[c]fluorene       | 17. Benzo[a]pyrene         |
| 2. Acenaphthylene | 10. Benz[a]anthracene     | 18. Indeno[1,2,3-cd]pyrene |
| 3. Acenaphthene   | 11. Cyclopenta[c,d]pyrene | 19. Dibenzo[a,h]anthracene |
| 4. Fluorene       | 12. Chrysene              | 20. Benzo[g,h,i]perylene   |
| 5. Phenanthrene   | 13. 5-Methylchrysene      | 21. Dibenzo[a,i]pyrene     |
| 6. Anthracene     | 14. Benzo[b]fluoranthene  | 22. Dibenzo[a,e]pyrene     |
| 7. Fluoranthene   | 15. Benzo[k]fluoranthene  | 23. Dibenzo[a,i]pyrene     |
| 8. Pyrene         | 16. Benzo[j]fluoranthene  | 24. Dibenzo[a,h]pyrene     |


## ZB-PAH-EU (cont'd)

### Better Combination of Resolution and Speed

Zebron outperforms popular GC columns for the separation of EU 15+1 and EPA 610 PAHs.


#### Zebron ZB-PAH-EU

30 meter x 0.25 mm x 0.20 μm




#### Popular Brand A

30 meter x 0.25 mm x 0.25 μm



#### Popular Brand B

30 meter x 0.25 mm x 0.10 μm



#### GC-MS conditions for both applications:

Column: As Indicated  
 Part No.: **ZHG-G043-10** (for ZB-PAH-EU)  
 Dimension: As indicated  
 Injection: Split 5:1 @ 330 °C, 1 μL  
 Recommended Liner: Zebron PLUS Single Taper Z-Liner™  
 Liner Part No.: **AG2-4B13-05** (for Shimadzu® 2010 GC)  
 Carrier Gas: Helium @ 24 psi (constant pressure)  
 Oven Program: 45 °C for 0.8 min to 200 °C @ 45 °C/min to 226 °C @ 3 °C/min for 0 min to 320 °C @ 10 °C/min for 20 min  
 Detector: MSD, 50-500 m/z  
 Transfer Line Temperature: 300 °C  
 Source Temperature: 300 °C

#### Sample:


- |                   |                           |                            |
|-------------------|---------------------------|----------------------------|
| 1. Naphthalene    | 9. Benzo[c]fluorene       | 17. Benzo[a]pyrene         |
| 2. Acenaphthylene | 10. Benz[a]anthracene     | 18. Indeno[1,2,3-cd]pyrene |
| 3. Acenaphthene   | 11. Cyclopenta[c,d]pyrene | 19. Dibenzo[a,h]anthracene |
| 4. Fluorene       | 12. Chrysene              | 20. Benzo[g,h,i]perylene   |
| 5. Phenanthrene   | 13. 5-Methylchrysene      | 21. Dibenzo[a,h]pyrene     |
| 6. Anthracene     | 14. Benzo[b]fluoranthene  | 22. Dibenzo[a,e]pyrene     |
| 7. Fluoranthene   | 15. Benzo[k]fluoranthene  | 23. Dibenzo[a,i]pyrene     |
| 8. Pyrene         | 16. Benzo[j]fluoranthene  | 24. Dibenzo[a,h]pyrene     |

Comparative separations may not be representative of all applications.

## ZB-PAH-EU (cont'd)

### GC-MS Analysis of PAHs in Rubber and Plastic

Zebron ZB-PAH-EU GC column separates the 18 PAH isomers within 12 minutes, resolving all critical pairs while demonstrating consistent column inertness.



The high temperature limits of the ZB-PAH-EU (340/360 °C) is extremely helpful for easy bakeout of dirty matrices and allows for eluting heavy PAHs!

**Column:** Zebron ZB-PAH-EU  
**Dimensions:** 10 meter x 0.10 mm x 0.08 µm  
**Part No.:** [7CB-G043-59](#)  
**Injection:** Split 5:1 @ 320 °C, 1 µL  
**Recommended Liner:** Zebron PLUS Single Taper Z-Liner™  
**Liner Part No.:** [AG2-4B13-05](#) (for Shimadzu® 2010 GC System)  
**Carrier Gas:** Helium @ 0.68 mL/min (constant flow)  
**Oven Program:** 100 °C for 3.0 min to 200 °C @ 60 °C/min to 270 °C @ 22 °C/min to 300 °C @ 4.5 °C/min to 330 °C @ 80 °C/min for 0.5 min  
**Detector:** MSD, 50-500 m/z  
**Transfer Line Temperature:** 300 °C  
**Source Temperature:** 330 °C


**Sample:**

1. Naphthalene	10. Chrysene
2. Acenaphthylene	11. Benzo[b]fluoranthene
3. Acenaphthene	12. Benzo[k]fluoranthene
4. Fluorene	13. Benzo[j]fluoranthene
5. Phenanthrene	14. Benzo[a]pyrene
6. Anthracene	15. Benzo[e]pyrene
7. Fluoranthene	16. Indeno[1,2,3-cd]pyrene
8. Pyrene	17. Dibenzo[a,h]anthracene
9. Benz[a]anthracene	18. Benzo[g,h,i]perylene

#### Ordering Information

Zebron ZB-PAH-EU GC Columns			
ID (mm)	df (µm)	Temp. Limits °C	Part No.
<b>10-Meter</b>			
0.10	0.08	40 to 340/360	<a href="#">7CB-G043-59</a>
<b>20-Meter</b>			
0.18	0.14	40 to 340/360	<a href="#">7FD-G043-47</a>
<b>30-Meter</b>			
0.25	0.20	40 to 340/360	<a href="#">7HG-G043-10</a>
<b>60-Meter</b>			
0.25	0.20	40 to 340/360	<a href="#">7KG-G043-10</a>

#### Column Profile



#### Phase Chemistry

- Proprietary

#### Recommended Applications

- Analysis of 15+1 EU-regulated and EPA regulated PAHs in food testing, rubber, plastic, coal
- Sources include cigarette smoke, vehicle exhausts, asphalt roads, coal, coal tar, wildfires, agricultural burning, residential wood burning, municipal, industrial waste incineration



Contact Phenomenex or your local Phenomenex distributor for additional GC products and applications.

## ZB-PAH-CT

- Enhanced resolution for chrysene and triphenylene (PAH Interferences)
- Increased Benzo[b,k] fluoranthene separation
- Optimal performance and resolution of regulated PAH isomers
- Great resolution of critical isomers, e.g. Benzo[b,j,k]fluoranthene



Use the column finder to replace or upgrade an existing column:  
[www.phenomenex.com/FindGC](http://www.phenomenex.com/FindGC)


### Priority PAH Analysis by GC

Zebtron ZB-PAH-CT columns are manufactured to provide the most optimal performance for EU-regulated polycyclic aromatic hydrocarbons (PAHs) and EPA regulated PAHs.

### Excellent Resolution of Chrysene and Triphenylene

We designed the Zebtron ZB-PAH-CT GC column to achieve complete resolution of Chrysene from Triphenylene along with other EU 15+1 PAH compounds. Its unique selectivity helps eliminate false positives while resolving PAH isomers, providing easy, fast, and accurate quantification of PAHs in environmental and food samples.

### EU 15+1 PAH Analysis Using Zebtron ZB-PAH-CT



**Column:** Zebtron ZB-PAH-CT  
**Dimensions:** 40 meter x 0.18 mm x 0.14 µm  
**Part No.:** [7PD-G044-47](#)  
**Injection:** Split 30:1 @ 320 °C, 1 µL  
**Recommended Liner:** Zebtron PLUS Single Taper Z-Liner™  
**Liner Part No.:** [AG2-4B13-05](#) (for Shimadzu® 2010 GC)  
**Carrier Gas:** Helium @ 78 psi (constant pressure)  
**Oven Program:** 45 °C for 0.8 min to 200 °C @ 45 °C/min to 265 °C @ 3 °C/min for 5 min to 270 °C @ 1 °C/min to 320 °C @ 10 °C/min for 15 min  
**Detector:** MSD (Shimadzu GC-[MS-QP2010](#) Ultra)  
**Mode:** SIM  
**SIM Ions:** 216, 226, 228, 242, 252, 276, 278, 302 m/z  
**Transfer Line Temperature:** 300 °C  
**Source Temperature:** 300 °C

- Sample:**
- |                          |                             |
|--------------------------|-----------------------------|
| 1. Benzo[c]fluorene      | 10. Benzo[a]pyrene          |
| 2. Cyclopenta[c,d]pyrene | 11. Indeno[1,2,3-c,d]pyrene |
| 3. Benz[a]anthracene     | 12. Dibenzo[a,h]anthracene  |
| 4. Triphenylene          | 13. Benzo[g,h,i]perylene    |
| 5. Chrysene              | 14. Dibenzo[a,i]pyrene      |
| 6. 5-Methylchrysene      | 15. Dibenzo[a,e]pyrene      |
| 7. Benzo[b]fluoranthene  | 16. Dibenzo[a,j]pyrene      |
| 8. Benzo[j]fluoranthene  | 17. Dibenzo[a,h]pyrene      |
| 9. Benzo[k]fluoranthene  |                             |



## ZB-PAH-CT (cont'd)


### Comparison of ZB-PAH-CT vs. Popular GC PAH column

Shorter Run Time and Better Sensitivity

#### Analysis of EU 15+1 and Triphenylene PAHs


##### Zebron ZB-PAH-CT

30 meter x 0.25 mm x 0.20 μm



##### Popular Brand A

30 meter x 0.25 mm x 0.15 μm



#### GC-MS conditions for both applications:

Column: As Indicated  
Dimensions: As Indicated  
Part No.: [7HG-G044-10](#) (Zebron ZB-PAH-CT)

Injection: Split 15:1 @ 320 °C, 1 μL

Recommended Liner: Zebron PLUS Single Taper Z-Liner™

Liner Part No.: [AG2-4B13-05](#)

Carrier Gas: Helium @ 23.7 psi (constant pressure)

Oven Program: 45 °C for 0.8 min to 200 °C @ 45 °C/min to 266 °C @ 3 °C/min for 0 min to 320 °C @ 10 °C/min to 320 °C for 20 min

Detector: MSD (Shimadzu® GC-MS-QP2010 Ultra)

Mode: SIM

SIM Ions: 216, 226, 228, 242, 252, 276, 278, 302 m/z

Transfer Line Temperature: 300 °C

Source Temperature: 300 °C

- Sample:
1. Benzo[c]fluorene
  2. Cyclopenta[c,d]pyrene
  3. Benz[a]anthracene
  4. Triphenylene
  5. Chrysene
  6. 5-Methylchrysene
  7. Benzo[b]fluoranthene
  8. Benzo[j]fluoranthene
  9. Benzo[k]fluoranthene
  10. Benzo[a]pyrene
  11. Indeno[1,2,3-c,d]pyrene
  12. Dibenzo[a,h]anthracene
  13. Benzo[g,h,i]perylene
  14. Dibenzo[a,i]pyrene
  15. Dibenzo[a,e]pyrene
  16. Dibenzo[a,i]pyrene
  17. Dibenzo[a,h]pyrene

#### Ordering Information

##### Zebron ZB-PAH-CT GC Columns

ID (mm)	df (μm)	Temp. Limits °C	Part No.
<b>10-Meter</b>			
0.10	0.10	40 to 320/340	<a href="#">7CB-G044-02</a>
<b>20-Meter</b>			
0.18	0.14	40 to 320/340	<a href="#">7FD-G044-47</a>
<b>30-Meter</b>			
0.25	0.20	40 to 320/340	<a href="#">7HG-G044-10</a>
<b>40-Meter</b>			
0.18	0.14	40 to 320/340	<a href="#">7PD-G044-47</a>

Comparative separations may not be representative of all applications.

**Column Profile**

**Phase Chemistry**

- Proprietary

**Recommended Applications**

- Analysis of 15+1 EU-regulated and EPA regulated PAHs in food testing, rubber, plastic, coal
- Sources include cigarette smoke, vehicle exhausts, asphalt roads, coal, coal tar, wildfires, agricultural burning, residential wood burning, municipal, industrial waste incineration.



Contact Phenomenex or your local Phenomenex distributor for additional GC products and applications.



Extend column lifetime. Add a Z-Guard™ to your next Zebron GC order.

## ZB-FAME

- Reduce traditional run times up to 75%
- Improve separation of cis/trans FAME isomers
- Suitable with AOAC, AOCS, and IOC methods

Upgrade to Zebron from any high cyanopropyl phase:


**Agilent®**

- CP-Sil 88
- HP-88
- DB®-23

**Supelco®**

- SP®-2380
- [SP-2560](#)

### Baseline Separation of Common Isomers



**Column:** Zebron ZB-FAME  
**Dimensions:** 30 meter x 0.25 mm x 0.20 µm  
**Part No.:** [ZHG-G033-10](#)  
**Injection:** Split 50:1 @ 240 °C, 1 µL  
**Recommended Liner:** Zebron PLUS Single Taper with Wool, 4 mm ID  
**Liner Part No.:** [AG2-0A11-05](#) (for Agilent® systems)  
**Carrier Gas:** Helium @ 1.2 mL/min (constant flow)  
**Oven Program:** 100 °C for 2 min to 140 °C @ 10 °C/min to 190 °C @ 3 °C/min to 260 °C @ 30 °C/min for 2 min  
**Detector:** FID @ 260 °C  
**Sample:** 37 FAME standard

### Easy Liner Selection




Our GC liner finder tool makes liner selection a breeze. You can even search by application, injection type, GC system, or your current liner part number.

[www.phenomenex.com/FindLiner](http://www.phenomenex.com/FindLiner)

## The Fast FAME GC Column

Traditionally, cis/trans FAME separations require the use of long (100 meters or more) columns and can run up to 60 minutes, resulting in a bottleneck to higher productivity. Zebron ZB-FAME provides targeted selectivity that allows for reduced column length – run times as short as 11 minutes without compromising your results!

### 37 FAMES In A Short 11 Minute Run



**Column:** Zebron ZB-FAME

**Dimensions:** 20 meter x 0.18 mm x 0.15 µm

**Part No.:** [7FD-G033-05](#)

**Injection:** Split 100:1 @ 250 °C, 1 µL

**Recommended Liner:** Zebron PLUS Single Taper Z-Liner™

**Liner Part No.:** [AG2-0A13-05](#) (for Agilent® systems)

**Carrier Gas:** Helium @ 1.0 mL/min (constant flow)

**Oven Program:** 80 °C for 1.5 min to 160 °C @ 40 °C/min to 185 °C @ 5 °C/min to 260 °C @ 30 °C/min

**Detector:** FID @ 260 °C

**Sample:** 37 FAME standard




### Ordering Information


Zebron ZB-FAME GC Columns			
ID (mm)	df (µm)	Temp. Limits °C	Part No.
<b>20-Meter</b>			
0.18	0.15	-20 to 280	<a href="#">7FD-G033-05</a>
<b>30-Meter</b>			
0.25	0.20	-20 to 280	<a href="#">7HG-G033-10</a>
<b>30-Meter with 5-Meter Guardian™ Integrated Guard</b>			
0.25	0.20	-20 to 280	<a href="#">7HG-G033-10-GGA</a>
<b>60-Meter</b>			
0.25	0.20	-20 to 280	<a href="#">7KG-G033-10</a>
<b>100-Meter</b>			
0.25	0.20	-20 to 280	<a href="#">7MG-G033-10</a>

Note: If you need a 5 in. cage, contact technical support at [www.phenomenex.com/chat](http://www.phenomenex.com/chat) or simply reach out to your Technical consultant. Conditions may apply. Agilent 6850, some SRI and process GC systems use only 5 in. cages.

### Column Profile



### Phase Chemistry



High Cyanopropyl

### Recommended Applications

- Fatty Acid Methyl Ester (FAMES)
- cis/trans FAME isomers
- Omega 3, Omega 6 FAMES

## ZB-SemiVolatiles

### Maximize Inertness

- Specifically designed to overcome obstacles for sensitive semi-volatiles methods
- **Enviro-Inert™ Technology** provides a rugged 5% phenyl-arylene phase – reduce activity without compromising selectivity
- Rugged QC test includes EPA 8270 tuning standard to ensure column is ready to pass suitability requirements
- Popular for EPA Methods 525, 610, 625, 8100, and 8270D

Upgrade to Zebron from any 5%-phenyl or 5% phenyl-arylene / 95% dimethylpolysiloxane phase:

#### Agilent®

- DB®-5ms
- DB-5ms Ultra Inert
- DB-5.625
- DB-UI 8270D

- HP-5ms
- HP-5ms Ultra Inert
- VP-5ms
- CP-5il 8 CB MS

#### Restek®

- Rxi®-5Sil MS
- Rxi-5ms


#### Supelco®

- SLB®-5ms


### 135 Compounds in Under 16 Minutes

ZB-SemiVolatiles provides improved productivity with shorter run times for EPA 8270D, while maintaining resolution of key critical pairs.


### Semivolatile Organic Compounds




### Improved Peak Shapes




1,4-Dioxane-D8 and 1,4-Dioxane




N-Nitrosodimethylamine and Pyridine




Indeno[1,2,3-cd]pyrene and Dibenz[a,h]anthracene, both share mass 276



Pentachlorophenol



2,4-Dinitrophenol



Benzo[b]fluoranthene and Benzo[k]fluoranthene

- Column:** Zebron ZB-SemiVolatiles
- Dimensions:** 30 meter x 0.25 mm x 0.25 μm
- Part No.:** 7HG-G027-11
- Injection:** Split 10:1 @ 280 °C, 1 μL
- Liner:** AGO-8499 (Single Taper with Wool)
- Septum:** AGO-4697 (PhenoRed™-400)
- Inlet Seal:** AGO-8620 (Easy Seals™ Inlet Base Seal)
- Carrier Gas:** Helium @ 1.4 mL/min (constant flow)
- Oven Program:** 40 °C for 0.5 min to 260 °C @ 40 °C/min to 295 °C @ 6 °C/min to 325 °C @ 25 °C/min for 2 min
- Detector:** MSD @ 340 °C; 45 – 450 amu
- Sample:** Analytes are 25 ppm in Dichloromethane  
135 compounds in EPA Method 8270D





## ZB-SemiVolatiles

### We QC Test For the Compounds You Analyze

We take the guesswork out of meeting method requirements by aggressively testing ZB-SemiVolatiles with two different test mixes. We incorporated troublesome analytes from your samples and compounds in the EPA 8270D tuning standard into our QC test, so you can be sure your column is ready to meet suitability requirements for the method.

### Meet Requirements Out-of-the-Box

Test Probe	Criteria	EPA Requirement	Our Requirement
<b>Pyridine</b> Very active amine that exposes even the smallest amount of column activity. This ensures that our Enviro-Inert™ deactivated column performs at the highest possible level for difficult basic compounds.	Peak Response	Not Specified	≥ 0.6
<b>Pentachlorophenol</b> Disappears and tails on active columns; it is important to measure relative response and peak skew criteria.	Peak Skew Peak Response	≤ 2.0 Not Specified	≤ 2.0 ≥ 0.3
<b>Benzidine</b> Active amine that tails when column activity is present, complicating peak quantification.	Peak Skew	≤ 2.0	≤ 2.0
<b>DDT</b> Breaks down in an active system to DDE and DDD. With our QC test, you are assured that your column will meet the EPA requirements upon installation.	Breakdown	< 20 %	< 20 %
<b>Injection</b> To ensure trace-level sensitivity, QC is performed with a 20 ppm mix using a 100:1 split injection – effectively 250 times less than the EPA maximum allowed.	Sensitivity	50 ng or less on column	0.2 ng on column

### Stands Up to Tough Samples for Increased Lifetime

“ I have found the Phenomenex ZB-SemiVolatiles columns to be superior in quality and durability than any other columns we have previously used. The columns not only last longer, but the reproducibility of column is extraordinary. The column holds calibrations particularly well, even after multiple injections of samples with far less than desirable matrices. All of this equates to less downtime and maintenance and more productivity for TestAmerica. ”

**Ryan McKernan, GC-MS Semi-Volatile Analyst**  
TestAmerica Laboratories, Inc. Buffalo

The opinions stated herein are solely those of the speaker and not necessarily those of any company or organization.

**Column Profile**

**Engineered Self Cross-linking™ (ESC)**

**Phase Chemistry**

5 % Phenyl-Arylene

C[Si](C)(C)c1ccc(cc1)[Si](C)(C)Oc2ccc(cc2)[Si](C)(C)Oc3ccc(cc3)C

95 % Dimethylpolysiloxane

**Recommended Applications**

- Semivolatiles (SVOCs)
- EPA Methods (525, 610, 625, 8100, 8270D)
- PAHs
- PBDEs

Contact Phenomenex or your local Phenomenex distributor for additional GC products and applications.

### Ordering Information

Zebtron ZB-SemiVolatiles GC Columns			
ID (mm)	df (µm)	Temp. Limits °C	Part No.
<b>15-Meter</b>			
0.25	0.25	-60 to 325/350	<a href="#">7EG-G027-11</a>
<b>20-Meter</b>			
0.18	0.18	-60 to 325/350	<a href="#">7FD-G027-08</a>
0.18	0.36	-60 to 325/350	<a href="#">7FD-G027-53</a>
<b>30-Meter</b>			
0.25	0.25	-60 to 325/350	<a href="#">7HG-G027-11</a>
0.25	0.50	-60 to 325/350	<a href="#">7HG-G027-17</a>
0.32	0.25	-60 to 325/350	<a href="#">7HM-G027-11</a>
<b>30-Meter with 5-Meter Guardian™ Integrated Guard</b>			
0.25	0.25	-60 to 325/350	<a href="#">7HG-G027-11-GGA</a>
0.25	0.50	-60 to 325/350	<a href="#">7HG-G027-17-GGA</a>
<b>30-Meter with 10-Meter Guardian Integrated Guard</b>			
0.25	0.25	-60 to 325/350	<a href="#">7HG-G027-11-GGC</a>
0.25	0.50	-60 to 325/350	<a href="#">7HG-G027-17-GGC</a>
<b>60-Meter</b>			
0.25	0.25	-60 to 325/350	<a href="#">7KG-G027-11</a>

Note: If you need a 5 in. cage, contact technical support at [www.phenomenex.com/chat](http://www.phenomenex.com/chat) or simply reach out to your Technical consultant. Conditions may apply. Agilent 6850, some SRI and process GC systems use only 5 in. cages.

## ZB-CLPesticides-1 and -2

### 7 EPA Methods, One Column Set

- **Guaranteed alternative to Restek Rtx-CLPesticides**
- **Optimized, versatile selectivity for chlorinated pesticides and herbicides**
- **Well-suited for dual-column configurations using GC-ECD**
- **Run EPA Methods 8081 and 8081 extended, 8082, 8151, 504, 505, 508, and 552 without changing columns – save time**

Upgrade to Zebron from these similar\* phases:

#### Restek®

- Rtx®-CLPesticides
- Rtx-CLPesticides2
- Stx®-CLPesticides
- Stx-CLPesticides2

\*not exact equivalent, selectivity may differ

### Direct Replacement for Restek Rtx-CLPesticides Phases


You asked for optimized performance for pesticides by GC-ECD detectors, without time-consuming method development. We've delivered a direct replacement\*\*! ZB-CLPesticides-1 and -2 provide

guaranteed drop-in performance compared to your current Rtx-CLPesticides column set, without the hassle.

### EPA 8081


#### Zebron ZB-CLPesticides-1

30 m x 0.32 mm x 0.32 µm




#### Zebron ZB-CLPesticides-2

30 m x 0.32 mm x 0.25 µm




#### Restek Rtx-CLPesticides

30 m x 0.32 mm x 0.32 µm



#### Restek Rtx-CLPesticides2

30 m x 0.32 mm x 0.25 µm



#### Conditions for all columns:

**Columns:** As listed  
**Dimensions:** As listed  
**Part No.:** [ZHM-G028-51](#) (ZB-CLPesticides-1)  
[ZHM-G029-11](#) (ZB-CLPesticides-2)  
**Injection:** Splitless (hold 0.3 min) @ 250 °C, 1 µL  
**Carrier Gas:** Helium @ 3.9 mL/min (constant flow)  
**Oven Program:** 120 °C to 200 °C @ 45 °C/min to 230 °C @ 15 °C/min to 330 °C @ 30 °C/min for 2 min

**Detector:** ECD @ 330 °C  
**Y-Connector:** [AGO-4717](#) (Fused Quartz)  
**Guard Column:** [ZAM-G000-00-GZ0](#) (5 m Z-Guard™)  
**Liner:** [AGO-8499](#) (Single Taper with Wool at Bottom)  
**Septum:** [AGO-4696](#) (PhenoRed™-400)  
**Inlet Seal:** [AGO-8620](#) (Gold-Plated Easy Seals™)  
**Sample:** Analytes are 250 ng/mL in hexane.

See page 117 for compound list.

\*\*Direct replacement: this category indicates an alternative column which will likely give a similar selectivity. Conditions for each method were the same for all columns tested. Comparative separations may not be representative of all applications.

## ZB-CLPesticides-1 and -2

### Five-Point Calibration Curve at 5, 15, 25, 100, and 250 ng/mL

Peak No.	Analyte	ZB-CLPesticides-1 % RSD*	ZB-CLPesticides-2 % RSD*	US EPA Specifications
1	2,4,5,6-TCMX (Surr)	3.8	3.0	< 20
2	α-BHC	8.3	3.8	< 20
3	γ-BHC	5.9	5.6	< 20
4	β-BHC	6.9	6.9	< 20
5	δ-BHC	4.9	5.7	< 20
6	Heptachlor	8.0	6.5	< 20
7	Aldrin	4.2	2.3	< 20
8	Heptachlor epoxide	3.8	2.3	< 20
9	trans-Chlordane	4.1	3.8	< 20
10	cis-Chlordane	4.0	3.3	< 20
11	4,4'-DDE	4.8	2.9	< 20
12	Endosulfan I	6.0	2.5	< 20
13	Dieldrin	7.7	4.9	< 20
14	Endrin	9.4	6.6	< 20
15	4,4'-DDD	9.2	3.6	< 20
16	Endosulfan II	6.6	4.1	< 20
17	4,4'-DDT	11.6	6.9	< 20
18	Endrin aldehyde	8.3	7.3	< 20
19	Endosulfan sulfate	8.0	7.1	< 20
20	Methoxychlor	6.7	6.1	< 20
21	Endrin ketone	6.5	7.2	< 20
22	Decachlorobiphenyl (Surr)	6.7	6.6	< 20
Average		6.6%	4.9%	< 20


\*Calculated using response factors as per EPA guidelines

## ZB-CLPesticides GC Column Kits


### Ordering Information

0.32 mm ID Kit (includes 1 of each below) Part No.: <a href="#">KG0-9286</a>		
Description	Dimension	Part No.
ZB-CLPesticides-1	30 meter x 0.32 mm x 0.32 μm	<a href="#">7HM-G028-51</a>
ZB-CLPesticides-2	30 meter x 0.32 mm x 0.25 μm	<a href="#">7HM-G029-11</a>
Z-Guard Column	5 meter x 0.32 mm	<a href="#">7AM-G000-00-G20</a>
Y-Connector	Fused Quartz	<a href="#">AGO-4717</a>
Polyimide Resin	0.5 mL, rated to 350 °C	<a href="#">AGO-5722</a>

0.53 mm ID Kit (includes 1 of each below) Part No.: <a href="#">KG0-9290</a>		
Description	Dimension	Part No.
ZB-CLPesticides-1	30 meter x 0.53 mm x 0.50 μm	<a href="#">7HK-G028-17</a>
ZB-CLPesticides-2	30 meter x 0.53 mm x 0.42 μm	<a href="#">7HK-G029-16</a>
Z-Guard Column	5 meter x 0.53 mm	<a href="#">7AK-G000-00-G20</a>
Y-Connector	Fused Quartz	<a href="#">AGO-4717</a>
Polyimide Resin	0.5 mL, rated to 350 °C	<a href="#">AGO-5722</a>

 Contact Phenomenex or your local Phenomenex distributor for additional GC products and applications.

**Column Profile**




**CLP-1**

**POLARITY\***

**TEMP LIMITS (°C)**  
40 to 320/340

\*Similar polarity to ZB-35.



**CLP-2**

**POLARITY\*\***

**TEMP LIMITS (°C)**  
40 to 320/340

\*\*Similar polarity to ZB-MultiResidue-2

**Phase Chemistry**

- Proprietary

**Recommended Applications**

- Dual-Column Chlorinated Pesticide Methods
- EPA Methods (8081 and 8081 extended, 8082, 8151, 504, 505, 508, 552)

## ZB-CLPesticides GC Columns

### Ordering Information

ZB-CLPesticides-1 GC Columns			
ID (mm)	df (μm)	Temp. Limits °C	Part No.
<b>30-Meter</b>			
0.25	0.25	40 to 320/340	<a href="#">7HG-G028-11</a>
0.32	0.32	40 to 320/340	<a href="#">7HM-G028-51</a>
0.32	0.50	40 to 320/340	<a href="#">7HM-G028-17</a>
0.53	0.50	40 to 320/340	<a href="#">7HK-G028-17</a>
<b>ZB-CLPesticides-2 GC Columns</b>			
ID (mm)	df (μm)	Temp. Limits °C	Part No.
<b>30-Meter</b>			
0.25	0.20	40 to 320/340	<a href="#">7HG-G029-10</a>
0.32	0.25	40 to 320/340	<a href="#">7HM-G029-11</a>
0.32	0.50	40 to 320/340	<a href="#">7HM-G029-17</a>
0.53	0.42	40 to 320/340	<a href="#">7HK-G029-16</a>

## ZB-MultiResidue™ -1 and -2

### Optimized Selectivity for Pesticides

- Specially designed for the separation of all types of pesticides, herbicides, and insecticides
- Baseline resolution and confirmation of all 20 chlorinated pesticides regulated under EPA Method 8081 in  $\leq 10$  min
- Decreased breakdown of sensitive pesticides such as DDT
- Robust performance for high temperature bakeouts
- Low bleed performance for pesticide confirmation by MS

Upgrade to Zebron from these similar\* phases:

#### Agilent®


- DB®-CLP1
- DB-CLP2

#### Restek®

- Rtx®-CLPesticides
- Rtx-CLPesticides2
- Stx®-CLPesticides
- Stx-CLPesticides2

\*not exact equivalent, selectivity may differ

### Improved Multi-Residue Pesticide Screening by GC-MS




**Column:** Zebron MultiResidue-1  
**Dimensions:** 30 meter x 0.25 mm x 0.25  $\mu$ m  
**Part No.:** ZHG-G016-11  
**Injection:** Splitless @ 260 °C, 1  $\mu$ L  
**Carrier Gas:** Helium @ 0.9 mL/min (constant flow)  
**Oven Program:** 80 °C for 0.5 min to 150 °C @ 10 °C/min to 240 °C @ 4 °C/min to 320 °C @ 15 °C/min for 3 min  
**Detector:** MSD @ 320 °C; 45-400 amu  
**Sample:** Analytes were 1 ppm in Dichloromethane




## ZB-MultiResidue™ -1 and -2 (cont'd)

### Resolve Common Pesticide Isomers


ZB-MultiResidue optimized selectivities improve resolution of complex pesticide, herbicide, and insecticide isomers. Our extremely stable siloxane-based polymer contains absolutely no nitrogen or halogenated functionality, which can be unfriendly to NPD and ECD detectors. Engineered Self-Crosslinking™ (ESC) bonding incorporates ladders into the phase backbone for low bleed and unmatched spectral integrity – even for trace-level samples.



#### Acquisition Spectra




#### NIST Reference Spectra




**Column:** Zebron ZB-MultiResidue-1  
**Dimensions:** 30 meter x 0.25 mm x 0.25 µm  
**Part No.:** 7HG-G016-11  
**Injection:** Splitless (hold 0.5 min) @ 260 °C, 1 µL  
**Carrier Gas:** Helium @ 0.8 mL/min (constant flow)  
**Oven Program:** 100 °C for 0.5 min to 200 °C @ 25 °C/min to 320 °C @ 15 °C/min for 2 min  
**Detector:** MSD @ 320 °C, 45-400 amu  
**Sample:**


1. o,p-DDE




2. p,p-DDE




3. o,p-DDD




4. o,p-DDT





5. p,p-DDD



6. p,p-DDT



#### Column Profile



#### Engineered Self Cross-linking™ (ESC)

#### Phase Chemistry




- Proprietary

#### Recommended Applications

- Haloacetic Acids (HAAs)
- Herbicides / Insecticides
- Multi-Pesticide Screening
- Nitrogen Containing Pesticides
- Organochlorine Pesticides
- Organophosphorous Pesticides
- PCBs / Aroclors

## ZB-MultiResidue™ -1 and -2 (cont'd)

### Great Results for Organophosphate Pesticides



ZEBRON UNLIMITED GC COLUMNS | GC COLUMNS

**Conditions for both column**

**Columns:** Zebron MultiResidue-1  
Zebron MultiResidue-2

**Dimensions:** 30 meter x 0.32 mm x 0.50 µm  
30 meter x 0.32 mm x 0.25 µm

**Part No.:** [7HM-G016-17](#)  
[7HM-G017-11](#)

**Injection:** On-Column @ 103 °C, 1 µL

**Carrier Gas:** Helium @ 2.8 mL/min (constant flow)

**Oven Program:** 100 °C for 0.5 min to 180 °C @ 20 °C/min to 240 °C @ 6 °C/min to 320 °C @ 15 °C/min for 2 min

**Detector:** FID @ 340 °C

Note: Columns connected using a 5 m Z-Guard™ Column and a 'Y' splitter.

**Sample:** Analytes are 2 ppm in Dichloromethane.

- |                                    |                         |                              |
|------------------------------------|-------------------------|------------------------------|
| 1. Dichlorvos                      | 17. Fonofos             | 33. Chlorfenvinphos          |
| 2. Mevinphos                       | 18. Phosphamidon Isomer | 34. Crotoxyphos              |
| 3. Trichlorfon                     | 19. Disulfoton          | 35. Stirofos                 |
| 4. TEPP (Tetraethyl Pyrophosphate) | 20. Dichlofenthion      | 36. Tokuthion                |
| 5. Demeton Isomer                  | 21. Phosphamidon        | 37. Merphos Oxide (Tribusof) |
| 6. Thionazin                       | 22. Chlorpyrifos Methyl | 38. Ethion                   |
| 7. Ethoprop                        | 23. Ronnel              | 39. Fensulfiothion           |
| 8. Sulfotep                        | 24. Aspon               | 40. Contaminant              |
| 9. Naled                           | 25. Methyl Parathion    | 41. Carbophenothion          |
| 10. Dicrotophos                    | 26. Malathion           | 42. Famfur                   |
| 11. Phorate                        | 27. Fenitrothion        | 43. EPN                      |
| 12. Monocrotophos                  | 28. Chlorpyrifos        | 44. Phosmet                  |
| 13. Demeton                        | 29. Fenthion            | 45. Leptophos                |
| 14. Terbufos                       | 30. Trichloronate       | 46. Azinphos Methyl          |
| 15. Diazinon                       | 31. Parathion           | 47. Azinphos Ethyl           |
| 16. Dimethoate                     | 32. Merphos             | 48. Coumaphos                |

## ZB-MultiResidue™ -1 and -2 (cont'd)

### Ordering Information

#### Zebron ZB-MultiResidue -1 GC Columns

ID(mm)	df(µm)	Temp. Limits °C	Part No.
<b>20-Meter</b>			
0.18	0.18	-60 to 320/340	<a href="#">7FD-G016-08</a>
<b>30-Meter</b>			
0.25	0.25	-60 to 320/340	<a href="#">7HG-G016-11</a>
0.32	0.25	-60 to 320/340	<a href="#">7HM-G016-11</a>
0.32	0.50	-60 to 320/340	<a href="#">7HM-G016-17</a>
0.53	0.50	-60 to 320/340	<a href="#">7HK-G016-17</a>

### Ordering Information

#### Zebron ZB-MultiResidue -2 GC Columns

ID(mm)	df(µm)	Temp. Limits °C	Part No.
<b>30-Meter</b>			
0.25	0.20	-60 to 320/340	<a href="#">7HG-G017-10</a>
0.32	0.25	-60 to 320/340	<a href="#">7HM-G017-11</a>

Note: If you need a 5 in. cage, contact technical support at [www.phenomenex.com/chat](http://www.phenomenex.com/chat) or simply reach out to your Technical consultant. Conditions may apply. Agilent 6850, some SRI and process GC systems use only 5 in. cages.



## ZB-MultiResidue Column Kits

### Ordering Information

<b>0.25 mm ID</b> (kit consists of products below)			
	<b>Part No.:</b>	<a href="#">KG0-8237</a>	
<b>Description</b>	<b>Dimension</b>	<b>Part No.</b>	
ZB-MultiResidue-1 Column	30 meter x 0.25 mm x 0.25 µm df	<a href="#">7HG-G016-11</a>	
ZB-MultiResidue-2 Column	30 meter x 0.25 mm x 0.20 µm df	<a href="#">7HG-G017-10</a>	
Z-Guard™	5 meter x 0.25 mm	<a href="#">7AG-G000-00-GZ0</a>	
Universal Capillary Column Y-connector, Fused Quartz	—	<a href="#">AGO-4717</a>	
Polyimide Resin	0.5 mL, rated to 350 °C	<a href="#">AGO-5722</a>	
<b>0.32 mm ID</b> (kit consists of products below)			
	<b>Part No.:</b>	<a href="#">KG0-8238</a>	
<b>Description</b>	<b>Dimension</b>	<b>Part No.</b>	
ZB-MultiResidue-1 Column	30 meter x 0.32 mm x 0.50 µm df	<a href="#">7HM-G016-17</a>	
ZB-MultiResidue-2 Column	30 meter x 0.32 mm x 0.25 µm df	<a href="#">7HM-G017-11</a>	
Z-Guard	5 meter x 0.32 mm	<a href="#">7AM-G000-00-GZ0</a>	
Universal Capillary Column Y-connector, Fused Quartz	—	<a href="#">AGO-4717</a>	
Polyimide Resin	0.5 mL, rated to 350 °C	<a href="#">AGO-5722</a>	
<b>0.53 mm ID</b> (kit consists of products below)			
	<b>Part No.:</b>	<a href="#">KG0-8239</a>	
<b>Description</b>	<b>Dimension</b>	<b>Part No.</b>	
ZB-MultiResidue-1 Column	30 meter x 0.53 mm x 0.50 µm df	<a href="#">7HK-G016-17</a>	
ZB-MultiResidue-2 Column	30 meter x 0.53 mm x 0.50 µm df	*	
Z-Guard	5 meter x 0.53 mm	<a href="#">7AK-G000-00-GZ0</a>	
Universal Capillary Column Y-connector, Fused Quartz	—	<a href="#">AGO-4717</a>	
Polyimide Resin	0.5 mL, rated to 350 °C	<a href="#">AGO-5722</a>	

\* Only available with kit order



Contact Phenomenex or your local Phenomenex distributor for additional GC products and applications.



Extend column lifetime. Add a Z-Guard™ to your next Zebron GC order.

## ZB-Bioethanol

### Quicker Bioethanol Testing

- Specially designed for fast and accurate bioethanol testing
- Provides accurate and reproducible results for Certificate of Analysis (COA)
- Resolve methanol and ethanol from all other denaturant peaks
- Great resolution of fusel alcohols
- Allows for quick bake out in between runs to eliminate contaminants



Upgrade to Zebron from traditional phases used for bioethanol:

Agilent®	Restek®	SGE®	Supelco®
• DB®-1	• Rtx®-1	• BP1	• SPB®-1
• HP-1	• Rxi®-1ms		• SE-30
• CP-Sil 5 CB			



### Fast, Accurate Analysis

Determination of Denatured Bioethanol: ASTM Method D5501



**Column:** Zebron ZB-1

**Dimensions:** 100 meter x 0.25 mm x 0.50 µm

**Part No.:** ZMG-G001-17

**Injection:** Split 50:1 @ 300 °C, 1 µL

**Carrier Gas:** Helium @ 35 cm/sec (constant flow)

**Oven Program:** 45 °C for 7 min to 255 °C @ 30 °C/min for 6 min

**Detector:** FID @ 300 °C

**Instrument:** Shimadzu® GC-2010 with Flame Ionization

**Sample:** 1. Methanol  
2. Ethanol  
3. Denaturant

**Column:** Zebron ZB-Bioethanol

**Dimensions:** 15 meter x 0.25 mm x 1.00 µm

**Part No.:** ZEG-G020-22

**Injection:** Split 50:1 @ 300 °C, 1 µL

**Carrier Gas:** Hydrogen @ 25 cm/sec (constant flow)

**Oven Program:** 55 °C for 1.7 min to 260 °C @ 40 °C/min (hold 2.67 min)

**Detector:** FID @ 300 °C

**Instrument:** Shimadzu® GC-2010 with Flame Ionization Detection and AOC-20i Automatic Liquid

**Sample:** 1. Methanol  
2. Ethanol  
3. Denaturant






2009 R&D 100  
Award Recipient

## ZB-Bioethanol

### Resolve Fusel Alcohols



**Column:** Zebron ZB-Bioethanol  
**Dimensions:** 30 meter x 0.25 mm x 1.00 µm  
**Part No.:** [7HG-G020-22](#)  
**Injection:** Split 100:1 @ 240 °C, 0.1 µL  
**Carrier Gas:** Helium @ 1.2 mL/min (constant flow)  
**Oven Program:** 40 °C for 5 min to 300 °C @ 25 °C/min  
**Detector:** MSD @ 230 °C; 30-450 amu

- Sample:**
- |                        |                            |
|------------------------|----------------------------|
| 1. Methanol            | 17. Methylcyclopentane     |
| 2. Ethanol             | 18. 2,4-Dimethylpentane    |
| 3. Acrolein            | 19. Benzene                |
| 4. Acetone             | 20. Cyclohexane            |
| 5. 2-Methylbutane      | 21. 2-Methylhexane         |
| 6. Isopropyl alcohol   | 22. 2,3-Dimethylpentane    |
| 7. Pentane             | 23. 3-Methylhexane         |
| 8. t-Butanol           | 24. 2,2,4-Trimethylpentane |
| 9. Allyl alcohol       | 25. Heptane                |
| 10. n-Propanol         | 26. Acetal                 |
| 11. 2,3-Dimethylbutane | 27. Toluene                |
| 12. 2-Methylpentane    | 28. Xylene                 |
| 13. 3-Methylpentane    | 29. Trimethylbenzene       |
| 14. 2-Butanol          |                            |
| 15. Ethyl acetate      |                            |
| 16. Hexane             |                            |


#### Ordering Information

##### Zebron ZB-Bioethanol GC Columns

ID(mm)	df(µm)	Temp. Limits °C	Part No.
<b>15-Meter</b>			
0.25	1.00	-60 to 340/360	<a href="#">7EG-G020-22</a>
<b>30-Meter</b>			
0.25	1.00	-60 to 340/360	<a href="#">7HG-G020-22</a>

Note: If you need a 5 in. cage, contact technical support at [www.phenomenex.com/chat](http://www.phenomenex.com/chat) or simply reach out to your Technical consultant. Conditions may apply. Agilent 6850, some SRI and process GC systems use only 5 in. cages.

#### Column Profile



Engineered Self Cross-linking™ (ESC)

#### Phase Chemistry

- Proprietary

#### Recommended Applications

- Alcohols
- Ethanol Testing
- Fusel Alcohols



Contact Phenomenex or your local Phenomenex distributor for additional GC products and applications.



For Bioethanol fermentation monitoring, use Rezex-ROA HPLC columns, see p. 325

## ZB-1XT SimDist

### High Efficiency Metal Column Performance

- Glass Infusion™ technology for higher efficiency and greater column-to-column reproducibility
- Individual QC testing for every column
- Up to 70% higher efficiency than other columns
- Increased accuracy for high temperature simulated distillation

Upgrade to Zebron from any 100% dimethylpolysiloxane phase:

**Agilent®**

- DB®-1
- DB-HT SimDis
- DB-PS1
- DB-PS2887
- CP-SimDist
- CP-SimDist UltiMetal


**Restek®**

- Rtx®-1
- Rxi®-1HT
- MXT®-1HT SimDist

**SGE®**

- BP1
- BPX1-SimD


### Glass Infusion Technology for Improved Performance



## ZB-1XT SimDist

### Improve Results for Simulated Distillation

Hydrocarbons C7–C100+: ASTM Method D7169



**Column:** Zebron ZB-1XT SimDist  
**Dimensions:** 5 meter x 0.53 mm x 0.15 μm  
**Part No.:** [7AK-G026-05](#)  
**Injection:** On-Column @ 33 °C, 1 μL  
**Carrier Gas:** Helium @ 7 mL/min (constant flow)  
**Oven Program:** 30 °C to 450 °C @ 10 °C/min for 10 min  
**Detector:** FID @ 450 °C  
**Sample:** C7 to C44 hydrocarbons and POLYWAX® 655 in CS<sub>2</sub>  
 Note: Chromatogram is baseline subtracted.


#### Ordering Information

##### Zebron ZB-1XT SimDist GC Columns

ID(mm)	df(μm)	Temp. Limits °C	Part No.
<b>5-Meter</b>			
0.53	0.09	-60 to 450	<a href="#">7AK-G026-55</a>
0.53	0.15	-60 to 450	<a href="#">7AK-G026-05</a>
0.53	0.88	-60 to 450	<a href="#">7AK-G026-49</a>
<b>5-Meter with 2-Meter Guardian™ Integrated Guard</b>			
0.53	0.09	-60 to 450	<a href="#">7AK-G026-55-GGT</a>
0.53	0.15	-60 to 450	<a href="#">7AK-G026-05-GGT</a>
<b>10-Meter</b>			
0.53	0.15	-60 to 450	<a href="#">7CK-G026-05</a>
0.53	0.88	-60 to 450	<a href="#">7CK-G026-49</a>
0.53	2.65	-60 to 400	<a href="#">7CK-G026-35</a>
<b>15-Meter</b>			
0.53	0.25	-60 to 450	<a href="#">7EK-G026-11</a>

Note: If you need a 5 in. cage, contact technical support at [www.phenomenex.com/chat](http://www.phenomenex.com/chat) or simply reach out to your Technical consultant. Conditions may apply. Agilent 6850, some SRI and process GC systems use only 5 in. cages.


#### Column Profile



\*Thicker film (2.65 μm) is rated to 400 °C.

#### Engineered Self Cross-linking™ (ESC)

##### Phase Chemistry



100 % Dimethylpolysiloxane

#### Recommended Applications

- ASTM Methods (D2887, D3710, D6352, D7169)
- Crude Oil
- Gasoline Fractions
- Petroleum Distillates
- Petroleum Fractions
- Simulated Distillation
- Vacuum Distillates



**ZB-1XT SimDist Test Mix**  
 Part No.: [AG0-8645](#)



Contact Phenomenex or your local Phenomenex distributor for additional GC products and applications.



Guard Column Connections  
 SiTite™ Mini-Unions for 0.8mm ID columns (P/N: [AG0-8825](#)) and Replacement Ferrules (P/N: [AG0-8824](#))

## ZB-DHA-PONA

- Excellent peak shape for polar and nonpolar compounds
- Temperature stability and flexibility
- Highly efficient dimension and consistent film thickness delivers excellent separation of paraffins, iso-paraffins, olefins, naphthenes, aromatics and polar compounds
- Extensive ESC™ provides intact stationary phase and MS certified low bleed
- Well-suited for true boiling point separation
- Excellent resolving power for critical pairs in complex petrochemical samples

### The Choice for PIONA, DHA and PONA

Zebron ZB-DHA-PONA is the choice for the analysis of Detailed Hydrocarbon Analysis (DHA) within the fuel industry. The Engineered Self Cross-linking™(ESC) stationary phase provides low bleed and exceptional column life for separation of DHA critical pairs with symmetric peaks. In addition, Zebron ZB-DHA-PONA GC columns provide excellent response and peak symmetry for polar oxygenates.

The Zebron ZB-DHA-PONA GC column provides excellent separation of DHA critical pairs with symmetric peaks suitable for ASTM Method D6730, see table below.

Upgrade to Zebron from any 100% dimethylpolysiloxane phase:

#### Agilent®


- HP-PONA
- DB®-PETRO
- CP-Sil PONA CB

#### Restek®

- Rtx®-DHA

#### Supelco®

- Petrocol®-DH




#### Easy ZB-DHA-PONA Selection for Your ASTM Method

Method	Description	Recommended Column	Recommended Dimensions	Part Number
ASTM D5134	Standard Test Method for Detailed Analysis of Petroleum Naphtha's through n-Nonane by Capillary Gas Chromatography	ZB-DHA-PONA	50 m x 0.20 mm x 0.5 µm	<a href="#">7JE-G042-17</a>
ASTM D5441	Standard Test Method for Analysis of Methyl Tert-Butyl Ether (MTBE) by GC	ZB-DHA-PONA	100 m x 0.25 mm x 0.5 µm	<a href="#">7MG-G042-17</a>
		ZB-DHA-PONA	50 m x 0.20 mm x 0.5 µm	<a href="#">7JE-G042-17</a>
		ZB-DHA-PONA	100 m x 0.25 mm x 0.5 µm	<a href="#">7MG-G042-17</a>
ASTM D5501	Standard Test Method for Determination of Ethanol and Methanol Content in Fuels Containing Greater than 20% Ethanol by Gas Chromatography	ZB-DHA-PONA	150 m x 0.25 mm x 1 µm	<a href="#">7QG-G042-22</a>
		ZB-DHA-PONA	100 m x 0.25 mm x 0.5 µm	<a href="#">7MG-G042-17</a>
ASTM D6729	Standard Test Method for Determination of Individual Components in Spark Ignition Engine Fuels by 100 Meter Capillary High Resolution Gas Chromatography	ZB-DHA-PONA	150 m x 0.25 mm x 1 µm	<a href="#">7QG-G042-22</a>
ASTM D6730	Standard Test Method for Determination of Individual Components in Spark Ignition Engine Fuels by 100-Meter Capillary (with Pre-column) High-Resolution Gas Chromatography	ZB-DHA-PONA	100 m x 0.25 mm x 0.5 µm	<a href="#">7MG-G042-17</a>
		ZB-DHA-PONA	50 m x 0.20 mm x 0.5 µm	<a href="#">7JE-G042-17</a>
		ZB-DHA-PONA	100 m x 0.25 mm x 0.5 µm	<a href="#">7MG-G042-17</a>
		ZB-DHA-PONA-TUNE	150 m x 0.25 mm x 1 µm	<a href="#">7QG-G042-22</a>
ASTM D6733	Standard Test Method for Determination of Individual Components in Spark Ignition Engine Fuels by 50-Meter Capillary High Resolution Gas Chromatography	ZB-DHA-PONA	5 m x 0.25 mm x 1 µm	<a href="#">7AG-G042-22</a>
ASTM D6733	Standard Test Method for Determination of Individual Components in Spark Ignition Engine Fuels by 50-Meter Capillary High Resolution Gas Chromatography	ZB-DHA-PONA	50 m x 0.20 mm x 0.5 µm	<a href="#">7JE-G042-17</a>

## ZB-DHA-PONA

Analysis of ASTM D6730 Components by GC-FID on ZB-DHA-PONA & ZB-DHA-PONA-TUNE GC Column



**Column 1 (Tuning):** Zebron ZB-DHA-PONA-TUNE  
**Phase:** 5% Phenyl 95% Dimethylpolysiloxane  
**Dimensions:** 5 meter x 0.25 mm x 1.00 µm  
**Part No.:** [7AG-G042-22](#)  
**Column 2:** Zebron ZB-DHA-PONA  
**Phase:** 100% Dimethylpolysiloxane  
**Dimensions:** 100 meter x 0.25 mm x 0.50 µm  
**Part No.:** [7MG-G042-17](#)

**Recommended Column Union:** [AG0-4716](#)  
**Injection:** Split 150:1 @ 200 °C, 0.2 µL  
**Recommended Liner:** Zebron PLUS Straight Z-Liner™  
**Part No.:** [AG2-0A03-05](#)  
**Carrier Gas:** Hydrogen @ 2 mL/min (constant flow)  
**Oven Program:** 30 °C for 8.5 min, to 48 °C @ 22 °C/min for 27 min, to 141 °C @ 3 °C/min for 1 min, to 275 °C @ 1 °C/min for 2 min  
**Detection:** Flame Ionization (FID) @ 275 °C

**Sample:**

1. Ethanol	16. C8 (n-octane)
2. C5 (n-pentane)	17. Ethylbenzene
3. 2-Methylbutane	18. 2,3-Dimethylheptane
4. Tert-Butanol	19. p-Xylene
5. 2,3-Dimethylbutane	20. C9 (n-nonane)
6. Methyl tert-butyl ether (MTBE)	21. 5-Methylnonane
7. C6 (n-hexane)	22. 1-Methyl-2-ethylbenzene
8. 1-Methylcyclopentene	23. C10 (n-decane)
9. Benzene	24. C11 (undecane)
10. Cyclohexane	25. 1,2,3,5-Tetramethylbenzene
11. 3-Ethylpentane	26. Naphthalene
12. trans-1,2-Dimethylcyclopentane	27. C12 (dodecane)
13. C7 (n-heptane)	28. 1-Methylnaphthalene
14. 2,3,3-Trimethylpentane	29. C13 (Tridecane)
15. Toluene	


### Ordering Information

#### Zebron ZB-DHA-PONA GC Columns

ID (mm)	df (µm)	Temp. Limits °C	Part No.
<b>5-Meter (PONA-TUNE)</b>			
0.25	1.00	-60 to 340/360	<a href="#">7AG-G042-22</a>
<b>50-Meter</b>			
0.20	0.50	-60 to 340/360	<a href="#">7GE-G042-17</a>
<b>100-Meter</b>			
0.25	0.50	-60 to 360/370	<a href="#">7MG-G042-17</a>
<b>150-Meter</b>			
0.25	1.00	-60 to 340/360	<a href="#">7QG-G042-22</a>

Note: If you need a 5 in. cage, contact technical support at [www.phenomenex.com/chat](http://www.phenomenex.com/chat) or simply reach out to your Technical consultant. Conditions may apply. Agilent 6850, some SRI and process GC systems use only 5 in. cages.


### Column Profile



\*Thicker films (≥ 1.0 µm) are rated to 340/360 °C.

### Engineered Self Cross-linking™ (ESC)

### Phase Chemistry



100% Dimethylpolysiloxane

### Recommended Applications

- DHA
- PONA
- PIONA
- PIANO
- ASTM D5134, D5441, D5501, D6729, D6730 and D6733

**i** Zebron GC Columns MS Certification, see p. 437

**i** Contact Phenomenex or your local Phenomenex distributor for additional GC products and applications.

**△** Extend column lifetime. Add a Z-Guard™ to your next Zebron GC order.



## ZB-Drug-1

### Faster Drugs of Abuse Testing

- Optimized phase for the separation of drugs of abuse
- Provides fast analysis with great peak shape
- Improves resolution of target analytes from matrix interferences
- Specially deactivated to improve quantitation for drug compounds

Upgrade to Zebron from traditional phases used for drugs of abuse:

#### Agilent®

- DB®-1ms
- DB-5ms
- DB-35

#### Restek®


- Rxi®-1ms
- Rtx®-5
- Rtx-5ms
- Rtx-35ms

#### Supelco®

- SPB®-1

### Optimized Selectivity for Multiple Drug Classes


#### Traditional 5% Phenyl Phase



**Column:** As listed  
**Dimensions:** 10 meter x 0.18 mm x 0.18 µm  
**Injection:** Split 10:1 @ 240 °C, 1 µL  
**Carrier Gas:** Helium @ 1.2 mL/min (constant flow)  
**Oven Program:** 140 °C to 320 °C @ 20 °C for 1 min  
**Detector:** MSD @ 230 °C

#### Zebron ZB-Drug-1


Unique Selectivity For Drugs



**Column:** As listed  
**Dimensions:** 10 meter x 0.18 mm x 0.18 µm  
**Part No.:** [ZCD-G023-08](#)  
**Injection:** Split 10:1 @ 280 °C, 1 µL  
**Carrier Gas:** Helium @ 55 cm/sec (constant flow)  
**Oven Program:** 180 °C to 340 °C @ 20 °C/min  
**Detector:** MSD @ 230 °C

**Sample:** Analytes @ 10 ppm BSTFA/Pyridine  
 1. Codeine-d6  
 2. Codeine  
 3. Morphine-d3  
 4. Morphine

### Common Drug Screen by GC-MS



**Column:** Zebron ZB-Drug-1  
**Dimensions:** 10 meter x 0.18 mm x 0.18 µm  
**Part No.:** [ZCD-G023-08](#)  
**Injection:** Split 10:1 @ 260 °C, 1 µL  
**Carrier Gas:** Helium @ 1 mL/min (constant flow)  
**Oven Program:** 50 °C to 150 °C @ 15 °C/min to 240 °C @ 7 °C/min to 320 °C @ 25 °C/min for 2 min  
**Detector:** MSD @ 320 °C; 45-450 amu

**Sample:** Analytes are 25 ppm in Methanol

- |                   |                      |                    |
|-------------------|----------------------|--------------------|
| 1. Acetophenone   | 15. Meprobamate      | 29. Trimipramine   |
| 2. Nicotine       | 16. Diphenhydramine  | 30. Chlorcyclizine |
| 3. Benzocaine     | 17. Lidocaine        | 31. Cocaine        |
| 4. Ibuprofen      | 18. Hexobarbital     | 32. Desipramine    |
| 5. Allobarbitol   | 19. Doxylamine       | 33. Codeine        |
| 6. Aprobarbital   | 20. Glutethimide     | 34. Morphine       |
| 7. Butalbital     | 21. Caffeine         | 35. Diazepam       |
| 8. Amobarbital    | 22. Chlorpheniramine | 36. Hydrocodone    |
| 9. Phenacetin     | 23. Methapyrilene    | 37. 6-MAM          |
| 10. Pentobarbital | 24. Phenobarbital    | 38. Oxycodone      |
| 11. Acetaminophen | 25. Procaine         | 39. Heroin         |
| 12. Benzphetamine | 26. Methadone        | 40. Fentanyl       |
| 13. Secobarbital  | 27. Brompheniramine  | 41. Ibogaine       |
| 14. Phencyclidine | 28. Propoxyphene     | 42. Triazolam      |
|                   |                      | 43. LSD            |




2009 R&D 100 Award Recipient

## ZB-Drug-1

### Faster Run Times and Improved Resolution

#### Zebron ZB-Drug-1




**Column:** Zebron ZB-Drug-1  
**Dimensions:** 10 meter x 0.18 mm x 0.18 μm  
**Part No.:** [7CD-G023-08](#)  
**Injection:** Split 10:1 @ 280 °C, 1 μL  
**Carrier Gas:** Helium @ 0.7 mL/min (constant flow)  
**Oven Program:** 200 °C to 210 °C @ 20 °C/min at 320 °C @ 30 °C/min for 1 min  
**Detector:** MSD @ 320 °C

#### Traditional Mid-Polar Phase

**Sample:**

1. Diazepam
2. Nordiazepam-d5
3. Nordiazepam
4. Oxazepam-d5
5. Oxazepam
6. Temazepam-d5
7. Temazepam
8. Clonazepam
9. Lorazepam
10. Alprazolam
11. α-Hydroxyalprazolam-d5
12. α-Hydroxyalprazolam



**Dimensions:** 10 meter x 0.18 mm x 0.18 μm  
**Injection:** Split 10:1 @ 250 °C, 1 μL  
**Carrier Gas:** Helium @ 0.6 mL/min (constant flow)  
**Oven Program:** 180 °C to 340 °C @ 20 °C/min for 2 min  
**Detector:** MSD @ 320 °C


#### Ordering Information

##### Zebron ZB-Drug-1 GC Columns

ID(mm)	df(μm)	Temp. Limits °C	Part No.
<b>10-Meter</b>			
0.18	0.18	40 to 320/340	<a href="#">7CD-G023-08</a>
<b>15-Meter</b>			
0.25	0.25	40 to 320/340	<a href="#">7EG-G023-11</a>
<b>15-Meter with 5-Meter Guardian™ Integrated Guard</b>			
0.25	0.25	40 to 320/340	<a href="#">7EG-G023-11-GGA</a>
<b>30-Meter</b>			
0.25	0.25	40 to 320/340	<a href="#">7HG-G023-11</a>

Note: If you need a 5 in. cage, contact technical support at [www.phenomenex.com/chat](http://www.phenomenex.com/chat) or simply reach out to your Technical consultant. Conditions may apply. Agilent 6850, some SRI and process GC systems use only 5 in. cages.

#### Column Profile



\*Similar polarity to ZB-MultiResidue™-2.


#### Engineered Self Cross-linking™ (ESC)

#### Phase Chemistry

- Proprietary

#### Recommended Applications

- Drug Screening
- 6-MAM
- Amphetamines
- Barbiturates
- Benzodiazepines
- PCP
- THC



**ZB-Drug-1 Test Mix**  
**Part No.:** [AGO-8431](#)



Contact Phenomenex or your local Phenomenex distributor for additional GC products and applications.



Extend column lifetime. Add a Z-Guard™ to your next Zebron GC order.

## ZB-BAC-1 and -2

### Optimized Pair for Blood Alcohol Testing

- Enhanced accuracy for post mortem samples
- Fast run time with baseline resolution of key components in just 2 minutes
- Enhanced resolution of ethanol and acetone peaks
- Achieve confirmation with two elution order changes when running columns in parallel
- Allows for the use of t-butanol or n-propanol as an internal standard

Upgrade to Zebron from these similar\* phases:


- |                 |                |
|-----------------|----------------|
| <b>Agilent®</b> | <b>Restek®</b> |
| • DB®-ALC1      | • Rtx®-BAC1    |
| • DB-ALC2       | • Rtx-BAC2     |

\*not exact equivalent, selectivity may differ

### Faster, More Sensitive Blood Alcohol Analysis


#### Zebron ZB-BAC-1

30 meter x 0.53 mm x 3.00 µm



#### Zebron ZB-BAC-2

30 meter x 0.53 mm x 2.00 µm



#### Conditions for both columns:


**Column:** As listed  
**Dimensions:** As listed  
**Part No.:** 7HK-G021-36 (ZB-BAC-1)  
 7HK-G022-32 (ZB-BAC-2)  
**Injection:** Split 0.8:1 @ 150 °C, 1 mL  
**Carrier Gas:** Helium @ 80 cm/sec (constant flow)  
**Oven Program:** 40 °C (Isothermal)  
**Detector:** FID @ 250 °C

**Sample:** Analytes 0.025% and internal standards 0.100% in water

1. Methanol
2. Acetaldehyde
3. Ethanol
4. Isopropanol
5. Acetone
6. t-Butanol (IS)
7. n-Propanol (IS)
8. 2-Butanol (IS)


#### Restek Rtx-BAC1

30 meter x 0.53 mm x 3.00 µm



#### Restek Rtx-BAC2

30 meter x 0.53 mm x 2.00 µm



#### Conditions for both columns:

**Column:** As listed  
**Dimensions:** As listed  
**Injection:** Split 5:1 @ 150 °C, 1 mL  
**Carrier Gas:** Helium @ 80 cm/sec (constant flow)  
**Oven Program:** 40 °C (Isothermal)  
**Detector:** FID @ 220 °C

**Sample:** Analytes and internal standards 0.100% in water


- |                 |                    |
|-----------------|--------------------|
| 1. Methanol     | 5. Acetone         |
| 2. Acetaldehyde | 6. t-Butanol (IS)  |
| 3. Ethanol      | 7. n-Propanol (IS) |
| 4. Isopropanol  | 8. 2-Butanol (IS)  |

Comparative separations may not be representative of all applications.


## ZB-BAC-1 and -2

### Run on Helium or Hydrogen

#### BAC-1



#### BAC-2



#### Conditions for both columns:

- Injection:** Split 5:1 @ 150 °C, 1 mL
- Carrier Gas:** Hydrogen @ 80 cm/sec (constant flow)
- Oven Program:** 40 °C (Isothermal)
- Detector:** FID @ 250 °C
- Sample:** Analytes are 0.100% in water
  1. Methanol
  2. Acetaldehyde
  3. Ethanol
  4. Isopropanol
  5. Acetone
  6. n-Propanol


#### Ordering Information

Zebron ZB-BAC-1 GC Columns			
ID(mm)	df( $\mu$ m)	Temp. Limits °C	Part No.
<b>30-Meter</b>			
0.32	1.80	-20 to 260/280	<a href="#">7HM-G021-31</a>
0.53	3.00	-20 to 260/280	<a href="#">7HK-G021-36</a>


Zebron ZB-BAC-2 GC Columns			
ID(mm)	df( $\mu$ m)	Temp. Limits °C	Part No.
<b>30-Meter</b>			
0.32	1.20	-20 to 260/280	<a href="#">7HM-G022-25</a>
0.53	2.00	-20 to 260/280	<a href="#">7HK-G022-32</a>

Note: If you need a 5 in. cage, contact technical support at [www.phenomenex.com/chat](http://www.phenomenex.com/chat) or simply reach out to your Technical consultant. Conditions may apply. Agilent 6850, some SRI and process GC systems use only 5 in. cages.

#### Column Profile



\*Similar polarity to ZB-35.



\*\*Similar polarity to ZB-624.

#### Engineered Self Cross-linking™ (ESC)

#### Phase Chemistry

- Proprietary

#### Recommended Applications

- Abused Inhalant Anesthetics
- Blood Alcohol Analysis



Contact Phenomenex or your local Phenomenex distributor for additional GC products and applications.



Extend column lifetime. Add a Z-Guard™ to your next Zebron GC order.

## ZB-1PLUS™

### MS Certified “1” Phase


- Very low bleed (MS Certified) phase especially suited to high sensitivity GC-MS
- Extremely inert for active compounds such as drugs, pesticides, or acids and bases
- Improved signal-to-noise ratio for better sensitivity and mass spectral integrity
- Identical selectivity to 100% dimethylpolysiloxane phases

Upgrade to Zebron from any 100% dimethylpolysiloxane phase:

Agilent®	Restek®	SGE®	Supelco®
• DB®-1	• Rtx®-1	• BP1	• SPB®-1
• DB-1ms	• Rtx-1ms	• SolGel-1ms™	• SE-30
• DB-1ms Ultra Inert	• Rxi®-1ms		• MET-1
• HP-1			• MDN-1
• HP-1ms			• Equity®-1
• HP-1ms Ultra Inert			
• VF-1ms			
• CP-Sil 5 CB			
• Ultra 1			

### Lower Overall Column Activity


Activity is a key measure of column quality. ZB-1PLUS columns are aggressively tested to ensure full deactivation. Our QC test below demonstrates low tailing on ZB-1PLUS for even the most active compounds, like 2-ethylhexanoic acid.



**Column:** Zebron ZB-1PLUS  
**Dimensions:** 30 meter x 0.25 mm x 0.25 µm  
**Part No.:** ZHG-G031-11  
**Injection:** Split 100:1 @ 250 °C, 1.0 µL  
**Carrier Gas:** Hydrogen @ 1.18 mL/min (constant flow)  
**Oven Program:** 140 °C (Isothermal)  
**Detector:** FID @ 325 °C  
**Sample:** 1. Decane  
 2. 2-Ethylhexanoic Acid  
 3. 4-Chlorophenol  
 4. Naphthalene  
 5. Tridecane  
 6. 1-Undecanol  
 7. Dicyclohexylamine  
 8. Pentadecane

### Lower Column Bleed

We tested the ZB-1PLUS column bleed profile against other “MS” columns on the market – ZB-1PLUS shows the lowest bleed, even at temperatures up to 360 °C.



Test conditions were stopped at 350 °C for all competitor columns so as not to cause damage to the stationary phase by exceeding their maximum temperature limit.

#### Conditions for all columns:


**Columns:** As listed  
**Dimensions:** 30 meter x 0.25 mm x 0.25 µm  
**Injection:** Null Injection @ 250 °C  
**Carrier Gas:** Hydrogen @ 100 mL/min (constant flow)  
**Oven Program:** 240 °C for 9 min to 280 °C for 6.3 min to 320 °C for 6.4 min to 340 °C for 5.8 min to 350 °C for 5.5 min to 360 °C  
**Detector:** FID @ 320 °C



## ZB-1PLUS™

### Well-Suited for Food & Flavors

#### Cold Pressed Orange Oil by GC-MS



**Column:** Zebron ZB-1PLUS  
**Dimensions:** 10 meter x 0.10 mm x 0.10 μm  
**Part No.:** 7CB-G031-02  
**Injection:** Split 120:1 @ 160 °C, 0.2 μL  
**Carrier Gas:** Helium @ 0.3 mL/min (constant flow)  
**Oven Program:** 60 °C to 130 °C @ 10 °C/min to 280 °C @ 30 °C/min for 3 min  
**Detector:** MSD  
**Sample:** Sample was 10% in dichloromethane


1. α-Pinene	9. cis-Limonene oxide	17. α-Cubebene
2. β-Phellandrine	10. trans-Limonene oxide	18. β-Cubebene
3. β-Myrcene	11. Citronellal	19. Dodecanal
4. Octanal	12. α-Terpineol	20. Valencene
5. 3-Carene	13. Decanal	21. Cadinene
6. Limonene	14. Carvone	22. Nootkatone
7. Nonanal	15. Neral	
8. Linalool	16. Geranial	

#### Ordering Information

Zebron ZB-1PLUS GC Columns			
ID(mm)	df(μm)	Temp. Limits °C	Part No.
<b>5-Meter</b>			
0.18	0.18	-60 to 360/370	<a href="#">7AD-G031-08</a>
<b>10-Meter</b>			
0.10	0.10	-60 to 360/370	<a href="#">7CB-G031-02</a>
0.18	0.18	-60 to 360/370	<a href="#">7CD-G031-08</a>
<b>12-Meter</b>			
0.20	0.33	-60 to 360/370	<a href="#">7DE-G031-14</a>
<b>15-Meter</b>			
0.25	0.25	-60 to 360/370	<a href="#">7EG-G031-11</a>
0.32	0.25	-60 to 360/370	<a href="#">7EM-G031-11</a>
<b>20-Meter</b>			
0.18	0.18	-60 to 360/370	<a href="#">7FD-G031-08</a>
<b>30-Meter</b>			
0.25	0.10	-60 to 360/370	<a href="#">7HG-G031-02</a>
0.25	0.25	-60 to 360/370	<a href="#">7HG-G031-11</a>
0.32	0.25	-60 to 360/370	<a href="#">7HM-G031-11</a>
<b>30-Meter with 5-Meter Guardian Integrated Guard</b>			
0.25	0.25	-60 to 360/370	<a href="#">7HG-G031-11-GGA</a>
<b>30-Meter with 10-Meter Guardian Integrated Guard</b>			
0.25	0.25	-60 to 360/370	<a href="#">7HG-G031-11-GGC</a>
<b>60-Meter</b>			
0.25	0.25	-60 to 360/370	<a href="#">7KG-G031-11</a>
0.25	1.00	-60 to 360/370	<a href="#">7KG-G031-22</a>


Note: If you need a 5 in. cage, contact technical support at [www.phenomenex.com/chat](http://www.phenomenex.com/chat) or simply reach out to your Technical consultant. Conditions may apply. Agilent 6850, some SRI and process GC systems use only 5 in. cages.

#### Column Profile



#### Engineered Self Cross-linking™ (ESC)

#### Phase Chemistry



100 % Dimethylpolysiloxane

#### Recommended Applications

- Acids
- Amines
- Diesel Fuel
- Drugs
- EPA Methods (1668)
- Essential Oils
- Flavors & Fragrances
- Oxygenates and GROs
- PCBs
- Pesticides
- Solvent Impurities
- Sulfur Compounds (Light)



Contact Phenomenex or your local Phenomenex distributor for additional GC products and applications.



Extend column lifetime. Add a Z-Guard™ to your next Zebron GC order.

## ZB-5PLUS™

### Inert 5% Phenyl Selectivity

- Highly inert—improved peak shape for acidic/basic compounds, drugs of abuse, and pesticides
- Very low bleed (MS certified) levels provide maximum sensitivity
- Intense QC specifications ensure column-to-column performance
- ESC™ bonding results in phase stability and high temperature limits
- Traditional bonding chemistry provides the same selectivity as the ZB-5 columns

Upgrade to Zebron from any

5% phenyl / 95% dimethylpolysiloxane phase:

#### Agilent®

- DB®-5
- HP-5
- HP-5ms
- HP-5msi

#### Restek®

- Rtx®-5
- Rtx-5MS
- Rtx-5Amine
- Rxi®-5ms

#### SGE®

- BP5
- BPX5


#### Supelco®

- MDN-5S
- SPB®-5
- Equity®-5

#### OV®

- OV-5


### Phenols



**Column:** Zebron ZB-5PLUS  
**Dimensions:** 30 meter x 0.25 mm x 0.25 µm  
**Part No.:** 7HG-G032-11  
**Injection:** Split 5:1 @ 240 °C, 1 µL  
**Carrier Gas:** Helium @ 1.2 mL/min (constant flow)  
**Oven Program:** 60 °C to 140 °C @ 5 °C/min to 280 °C @ 10 °C/min  
**Detector:** MSD @ 230 °C, 45-450 amu  
**Sample:**

1. Phenol	16. 2,3-Dimethylphenol
2. 2-Chlorophenol	17. 3,4-Dimethylphenol
3. 2-Methylphenol	18. 2,6-Dichlorophenol
4. 4-Methylphenol	19. 4-Chloro-3-methylphenol
5. 3-Methylphenol	20. 2,4,6-Trichlorophenol
6. 2,6-Dimethylphenol	21. 2,4,5-Trichlorophenol
7. 2-Nitrophenol	22. 2,4-Dinitrophenol
8. 2-Ethylphenol	23. 4-Nitrophenol
9. 2,4-Dimethylphenol	24. 2,3,4,6-Tetrachlorophenol
10. 3,5-Dimethylphenol	25. 4,6-Dinitro-2-methylphenol
11. 2,5-Dimethylphenol	26. Pentachlorophenol
12. 4-Ethylphenol	27. Dinoseb
13. 3-Ethylphenol	
14. 2,4-Dichlorophenol	
15. Benzoic Acid	


### Underivatized Antihistamines by GC-FID



**Column:** Zebron ZB-5PLUS  
**Dimensions:** 30 meter x 0.25 mm x 1.0 µm  
**Part No.:** 7HG-G032-22  
**Injection:** Split 50:1 @ 305 °C, 1 µL  
**Carrier Gas:** Helium @ 1.3 mL/min (constant flow)  
**Oven Program:** 40 °C for 1 min to 240 °C @ 25 °C/min for 2 min to 305 °C @ 25 °C/min for 8 min  
**Detector:** FID @ 320 °C  
**Sample:**

1. Phenylpropanolamine	7. Phenyltoloxamine
2. Ephedrine	8. Methapyrilene
3. Pseudoephedrine	9. Chlorpheniramine
4. Pheniramine	10. Brompheniramine
5. Diphenhydramine	11. Triprolidine
6. Doxylamine	

### Endocrine Disruptors by GC-MS




**Column:** Zebron ZB-5PLUS  
**Dimensions:** 30 meter x 0.25 mm x 0.25 µm  
**Part No.:** 7HG-G032-11  
**Injection:** Split 40:1 @ 250 °C, 1 µL  
**Carrier Gas:** Helium @ 1.2 mL/min (constant flow)  
**Oven Program:** 100 °C to 180 °C @ 5 °C/min to 320 °C @ 15 °C/min  
**Detector:** MSD @ 180 °C, 45-450 amu  
**Sample:** Analytes are 50 ppm in acetone

1. Dimethyl phthalate	9. Terbutylazine	17. 4,4'-DDD
2. Diethyl phthalate	10. Secbumeton	18. Di-n-hexyl phthalate
3. Atraton	11. Simetryn	19. 4,4'-DDT
4. Simazine	12. Ametryn	20. Dicyclohexyl phthalate
5. Prometon	13. Prometryn	21. bis(2-Ethylhexyl)phthalate
6. Atrazine	14. Terbutryn	22. Di-n-octyl phthalate
7. Propazine	15. Dibutyl phthalate	23. Ethinyl estradiol
8. Dipropyl phthalate	16. 4,4'-DDE	

## ZB-5PLUS™

### Good Results for Drugs

Drug Screening by GC-MS



**Column:** Zebron ZB-5PLUS  
**Dimensions:** 30 meter x 0.25 mm x 0.25 µm  
**Part No.:** 7HG-G032-11  
**Injection:** Split 15:1 @ 240 °C, 1 µL  
**Carrier Gas:** Helium @ 1.1 mL/min (constant flow)  
**Oven Program:** 140 °C to 240 °C @ 10 °C/min for 5 min to 320 °C @ 25 °C/min for 2.25 min  
**Detector:** MSD @ 230 °C, 45-450 amu

**Sample:** Analytes (underivatized) are 25 ppm in dichloromethane

- |                      |                          |                               |
|----------------------|--------------------------|-------------------------------|
| 1. Nicotine          | 14. Caffeine             | 27. Morphine                  |
| 2. Methylecgonine    | 15. Benzphetamine        | 28. Diazepam                  |
| 3. Ibuprofen         | 16. Hexobarbital         | 29. Hydrocodone               |
| 4. Allobarbitol      | 17. Dimenhydrinate       | 30. 6-Monoacetylmorphine      |
| 5. Aprobital         | 18. Doxylamine           | 31. Oxymorphone               |
| 6. Butobarbital      | 19. Phenobarbital        | 32. Diacetylmorphine (Heroin) |
| 7. Acetaminophen     | 20. 8-Chlorotheophylline |                               |
| 8. Phenacetin        | 21. Methapyrilene        |                               |
| 9. Amobarbital       | 22. Chlorpheniramine     |                               |
| 10. Pentobarbital    | 23. Brompheniramine      |                               |
| 11. Secobarbital     | 24. Cocaine              |                               |
| 12. Meprobamate      | 25. Chlorcyclizine       |                               |
| 13. Methyl benzilate | 26. Codeine              |                               |


### Ordering Information

#### Zebron ZB-5PLUS GC Columns

ID(mm)	df(µm)	Temp. Limits °C	Part No.
<b>15-Meter</b>			
0.25	0.25	-60 to 360/370	<a href="#">7EG-G032-11</a>
<b>20-Meter</b>			
0.18	0.18	-60 to 360/370	<a href="#">7FD-G032-08</a>
<b>20-Meter with 5-Meter Guardian™ Integrated Guard</b>			
0.18	0.18	-60 to 360/370	<a href="#">7FD-G032-08-GGA</a>
<b>30-Meter</b>			
0.25	0.25	-60 to 360/370	<a href="#">7HG-G032-11</a>
0.25	0.50	-60 to 360/370	<a href="#">7HG-G032-17</a>
0.25	1.00	-60 to 360/370	<a href="#">7HG-G032-22</a>
0.32	0.25	-60 to 360/370	<a href="#">7HM-G032-11</a>
<b>30-Meter with 5-Meter Guardian Integrated Guard</b>			
0.25	0.10	-60 to 360/370	<a href="#">7HG-G032-02-GGA</a>
0.25	0.25	-60 to 360/370	<a href="#">7HG-G032-11-GGA</a>
<b>60-Meter</b>			
0.25	0.25	-60 to 360/370	<a href="#">7KG-G032-11</a>


Note: If you need a 5 in. cage, contact technical support at [www.phenomenex.com/chat](http://www.phenomenex.com/chat) or simply reach out to your Technical consultant. Conditions may apply. Agilent 6850, some SRI and process GC systems use only 5 in. cages.

### Column Profile



Engineered Self Cross-linking™ (ESC)

### Phase Chemistry



5 % Phenyl 95 % Dimethylpolysiloxane

### Recommended Applications

- Barbiturates
- Benzodiazepines
- Drugs of Abuse
- EPA Methods
- FAMES
- Nitrosamines
- Pesticides
- Phenols
- THC Metabolites



For high temperature analysis, consider using a ZB-5HT, see p. 146



Contact Phenomenex or your local Phenomenex distributor for additional GC products and applications.



Extend column lifetime. Add a Z-Guard™ to your next Zebron GC order.

## ZB-5MS<sup>PLUS</sup>™

### The Next Generation of Inertness

- The next generation of inertness for specialty chemical, forensic, toxicology, and food testing applications
- Specialized deactivation for versatile 5% phenyl-arylene selectivity with improved sensitivity
- Low bleed (MS Certified) and well-suited to high sensitivity GC-MS and GC-MS/MS work

Upgrade to Zebron from any 5% phenyl or 5% phenyl-arylene / 95% dimethylpolysiloxane phase:

#### Agilent®

- DB®-5ms
- DB-5ms Ultra Inert
- HP-5ms
- HP-5ms Ultra Inert
- VF-5ms

#### Restek®

- Rxi®-5Sil MS


#### Supelco®

- SLB®-5ms

### Engineered for High Performance

Active sites on a GC column's surface can result in analyte adsorption and degradation, negatively affecting peak shape and response. To reduce potential surface activity, Zebron ZB-5MS<sup>PLUS</sup> is designed with a rigorous fused silica deactivation process that improves inertness for troublesome compounds. Instantly achieve higher responses for active compounds compared to your current 5ms phase column, without changing your selectivity.

#### Melamine in Dog Food by GC-MS



- Extraction Protocol:**
1. Combine 0.5 g of homogenized dog food with 10 mL of DEA/Water/ Acetonitrile (1:4:5) in a 15 mL centrifuge tube
  2. Sonicate for 30 min
  3. Centrifuge at 5000 rpm for 10 min
  4. Transfer 100 µL of supernatant to an autosampler vial and evaporate to dryness using nitrogen gas
  5. Reconstitute with 100 µL of Acetonitrile/Pyridine (1:1) and then derivatize using 100 µL BSTFA with 1% TCMS at 70 °C for 45 min

**Column:** Zebron ZB-5MS<sup>PLUS</sup>

**Dimensions:** 30 meter x 0.25 mm x 0.25 µm

**Part No.:** 7HG-G030-11


**Injection:** Splitless @ 280 °C, 1 µL

**Carrier Gas:** Helium @ 1 mL/min (constant flow)

**Oven Program:** 75 °C for 1 min to 320 °C @ 15 °C/min hold for 4 min

**Detector:** MSD @ 320 °C


- Sample:**
1. Cyanuric acid
  2. Ammelide
  3. Ammeline
  4. Melamine
  5. Benzoguanamine



## ZB-5MSPLUS™

### Versatile Performance for Drugs and Chemicals

#### Cannabinoids in Brownies by GC-MS




#### Extraction Protocol:

1. Combine 1 g of chocolate brownie with 10 mL of water in a 50 mL centrifuge tube
2. Shake using a mechanical shaker until dissolved
3. Add roQ™ QuEChERS EN15662 extraction salt packet (KSO-8909) and 10 mL of acetonitrile
4. Shake tube for 3 min using mechanical shaker
5. Centrifuge at 2700 rpm for 5 min
6. Transfer 1 mL of supernatant to an autosampler vial for GC-MS analysis

**Column:** Zebron ZB-5MSPLUS  
**Dimensions:** 30 meter x 0.25 mm x 0.25 µm  
**Part No.:** [7HG-G030-11](#)  
**Injection:** Splitless @ 250 °C, 1 µL  
**Carrier Gas:** Helium @ 1.5 mL/min (constant flow)  
**Oven Program:** 100 °C for 1 min to 320 °C @ 50 °C/min, hold for 2 min  
**Detector:** MSD @ 320 °C  
**Sample:** 1. Cannabidiol  
 2. Δ-9-Tetrahydrocannabinol  
 3. Cannabinol

#### Ethanolamines by GC-MS



**Column:** Zebron ZB-5MSPLUS  
**Dimensions:** 30 meter x 0.25 mm x 1.00 µm  
**Part No.:** [7HG-G030-22](#)  
**Injection:** Split 200:1 @ 250 °C, 1 µL  
**Carrier Gas:** Helium @ 1.4 mL/min (constant flow)  
**Oven Program:** 30 °C to 300 °C @ 40 °C/min  
**Detector:** MSD @ 320 °C  
**Sample:** 1. Monoethanolamine  
 2. Diethanolamine  
 3. Triethylene glycol monomethyl ether (IS)  
 4. Triethanolamine

#### Ordering Information

##### Zebron ZB-5MSPLUS GC Columns

ID(mm)	df(µm)	Temp. Limits °C	Part No.
<b>1.5-Meter</b>			
0.25	0.25	-60 to 325/350	<a href="#">7XG-G030-11</a>
<b>15-Meter</b>			
0.25	0.25	-60 to 325/350	<a href="#">7EG-G030-11</a>
<b>20-Meter</b>			
0.18	0.18	-60 to 325/350	<a href="#">7FD-G030-08</a>
0.18	0.36	-60 to 325/350	<a href="#">7FD-G030-53</a>
<b>30-Meter</b>			
0.25	0.25	-60 to 325/350	<a href="#">7HG-G030-11</a>
0.25	0.50	-60 to 325/350	<a href="#">7HG-G030-17</a>
0.25	1.00	-60 to 325/350	<a href="#">7HG-G030-22</a>
0.32	0.25	-60 to 325/350	<a href="#">7HM-G030-11</a>
0.32	0.50	-60 to 325/350	<a href="#">7HM-G030-17</a>
0.32	1.00	-60 to 325/350	<a href="#">7HM-G030-22</a>
0.32	1.50	-60 to 325/350	<a href="#">7HM-G030-28</a>
0.53	1.00	-60 to 325/350	<a href="#">7HK-G030-22</a>


#### Ordering Information

##### Zebron ZB-5MSPLUS GC Columns (cont'd)

ID(mm)	df(µm)	Temp. Limits °C	Part No.
<b>30-Meter with 5-Meter Guardian™ Integrated Guard</b>			
0.25	0.25	-60 to 325/350	<a href="#">7HG-G030-11-GGA</a>
0.25	0.50	-60 to 325/350	<a href="#">7HG-G030-17-GGA</a>
<b>30-Meter with 10-Meter Guardian Integrated Guard</b>			
0.25	0.25	-60 to 325/350	<a href="#">7HG-G030-11-GGC</a>
<b>60-Meter</b>			
0.25	0.25	-60 to 325/350	<a href="#">7KG-G030-11</a>
0.25	1.00	-60 to 325/350	<a href="#">7KG-G030-22</a>
0.32	1.00	-60 to 325/350	<a href="#">7KM-G030-22</a>

Note: If you need a 5 in. cage, contact technical support at [www.phenomenex.com/chat](http://www.phenomenex.com/chat) or simply reach out to your Technical consultant. Conditions may apply. Agilent 6850, some SRI and process GC systems use only 5 in. cages.


#### Column Profile



#### Engineered Self Cross-linking™ (ESC)

#### Phase Chemistry

5 % Phenyl-Arylene



95 % Dimethylpolysiloxane

#### Recommended Applications

- Acids
- Alkaloids
- Amines
- Drugs
- Essential Oils
- Flavors
- Halo-hydrocarbons
- Pesticides
- Phenols
- Residual Solvents
- Solvent Impurities



Contact Phenomenex or your local Phenomenex distributor for additional GC products and applications.



Extend column lifetime. Add a Z-Guard™ to your next Zebron GC order.



## ZB-WAXPLUS™

### Enhanced Aqueous Stability

- 100% aqueous stable, excellent for aqueous samples
- Extremely inert for acidic compounds
- Enhanced selectivity for low boiling solvents
- High retention of alcohols and chlorinated solvents
- Increased efficiency at 20 °C

Upgrade to Zebron from any polyethylene glycol phase:

#### Agilent®

- DB®-WAX
- CAM
- HP-20M
- Carbowax 20M
- CP-Wax 52 CB

#### Restek®

- Stabilwax®

#### SGE®


- BP20

#### Supelco®

- SUPELCO WAX® 10


### Water Reproducibility of ZB-WAXPLUS

Historically, polyethylene glycol (PEG) phases have been unstable with aqueous samples such as beverages or glycols, resulting in poor reproducibility and decreased lifetime. ZB-WAXPLUS bonding procedure results in exceptional stability to repeated injections of aqueous matrices.




**Column:** Zebron ZB-WAXPLUS  
**Dimensions:** 30 meter x 0.25 mm x 0.25 µm  
**Part No.:** [7HG-G013-11](#)  
**Injection:** Split 30:1 @ 140 °C, 0.2 µL  
**Carrier Gas:** Helium @ 1.4 mL/min (constant flow)  
**Oven Program:** 35 °C for 5 min to 85 °C @ 10 °C/min to 200 °C @ 25 °C/min for 1 min  
**Detector:** FID @ 200 °C  
**Sample:** 1. Acetaldehyde  
 2. Ethyl Acetate  
 3. Methanol  
 4. Ethanol  
 5. Propanol  
 6. Isobutanol  
 7. 2-Methylbutanol  
 8. 3-Methylbutanol

### Improve Resolution for ASTM D5135-95: Styrene Impurities



#### Restek Stabilwax



Conditions same for both columns:

**Dimensions:** 30 meter x 0.25 mm x 0.25 µm  
**Part No.:** [7HG-G013-11](#)  
**Injection:** Split 100:1 @ 250 °C, 1 µL  
**Carrier Gas:** Hydrogen @ 1.0 mL/min (constant flow)  
**Oven Program:** 35 °C for 2.5 min to 85 °C @ 10 °C/min and hold until last peak elutes  
**Detector:** FID @ 225 °C


**Sample:** 1. Methyl Formate  
 2. Acetone  
 3. Ethyl Acetate  
 4. Methyl Ethyl Ketone  
 5. Methanol  
 6. 2-Methyl-2-propanol  
 7. Methylene Chloride  
 8. Benzene  
 9. Ethanol  
 10. 2-Butanol  
 11. Toluene  
 12. n-Propanol  
 13. Ethyl Benzene  
 14. p-Xylene  
 15. m-Xylene  
 16. 1-Butanol  
 17. o-Xylene

Comparative separations may not be representative of all applications.

## ZB-WAXPLUS™

### A Food Testing Must-Have

#### Food Industry FAMES




**Column:** Zebron ZB-WAXPLUS  
**Dimensions:** 30 meter x 0.25 mm x 0.25 µm  
**Part No.:** [7HG-G013-11](#)  
**Injection:** Split 5:1 @ 220 °C, 1 µL  
**Carrier Gas:** Helium @ 3 mL/min (constant flow)  
**Oven Program:** 60 °C for 2 min to 150 °C @ 13 °C/min to 240 °C @ 2 °C/min  
**Detector:** FID @ 250 °C  
**Sample:** See the full compound list at [www.phenomenex.com/GC](http://www.phenomenex.com/GC)

#### Ordering Information

Zebron ZB-WAXPLUS GC Columns			
ID(mm)	df(µm)	Temp. Limits °C	Part No.
<b>10-Meter</b>			
0.10	0.10	20 to 250/260	<a href="#">7CB-G013-02</a>
<b>15-Meter</b>			
0.25	0.25	20 to 250/260	<a href="#">7EG-G013-11</a>
0.53	1.00	20 to 230/240	<a href="#">7EK-G013-22</a>
<b>20-Meter</b>			
0.18	0.18	20 to 250/260	<a href="#">7FD-G013-08</a>
<b>30-Meter</b>			
0.25	0.25	20 to 250/260	<a href="#">7HG-G013-11</a>
0.25	0.50	20 to 250/260	<a href="#">7HG-G013-17</a>
0.32	0.25	20 to 250/260	<a href="#">7HM-G013-11</a>
0.32	0.50	20 to 250/260	<a href="#">7HM-G013-17</a>
0.32	1.00	20 to 230/240	<a href="#">7HM-G013-22</a>
0.53	0.25	20 to 250/260	<a href="#">7HK-G013-11</a>
0.53	1.00	20 to 230/240	<a href="#">7HK-G013-22</a>
<b>60-Meter</b>			
0.25	0.15	20 to 250/260	<a href="#">7KG-G013-05</a>
0.25	0.25	20 to 250/260	<a href="#">7KG-G013-11</a>
0.25	0.50	20 to 250/260	<a href="#">7KG-G013-17</a>
0.32	0.25	20 to 250/260	<a href="#">7KM-G013-11</a>
0.32	0.50	20 to 250/260	<a href="#">7KM-G013-17</a>
0.53	1.00	20 to 230/240	<a href="#">7KK-G013-22</a>


Note: If you need a 5 in. cage, contact technical support at [www.phenomenex.com/chat](http://www.phenomenex.com/chat) or simply reach out to your Technical consultant. Conditions may apply. Agilent 6850, some SRI and process GC systems use only 5 in. cages.

#### Column Profile



\*Thicker films (≥ 1.0 µm) are rated to 230/240 °C.

#### Phase Chemistry



100 % Polyethylene Glycol

#### Recommended Applications

- Alcohols
- Aldehydes
- Aromatics
- Essential Oils
- Flavors & Fragrances
- Free Fatty Acids
- Glycols
- OVIs
- Pharmaceuticals
- Solvents / Residual Solvents
- Styrene
- Xylene Isomers



**ZB-WAXPLUS Test Mix**  
**Part No.:** [AG0-7869](#)



Contact Phenomenex or your local Phenomenex distributor for additional GC products and applications.



Extend column lifetime.  
 Add a Z-Guard™ to your next Zebron GC order.






## ZB-624PLUS™

### We QC Test for the Compounds You Analyze

We added challenging and troublesome analytes to our QC test to make sure each ZB-624PLUS column has superior deactivation.


Test Probe	The PLUS Advantage	Property
2,4-Dimethylphenol 2,4-Dimethylaniline	We screen challenging analytes, like acids and bases, to mimic your most challenging compounds.	Inertness



#### Conditions for all separations:

- Column:** Zebron ZB-624PLUS
- Dimensions:** 30 meter x 0.32 mm x 1.80 μm
- Part No.:** [7HM-G040-31](#)
- Injection:** Split 50:1 @ 250°C, 1 μL
- Recommended Liner:** Zebron PLUS Straight Z-Liner™
- Liner Part No.:** [AG2-0A03-05](#) (for Agilent® & Thermo Scientific® systems)
- Carrier Gas:** Hydrogen @ 6 psi (constant pressure)
- Oven Program:** 85°C for 10.5 min
- Detector:** FID @ 305°C
- Sample:**
  1. Methane
  2. Dodecane
  3. 2,4-Dimethylphenol
  4. 2,4-Dimethylaniline
  5. Tridecane
  6. 1-Methylnaphthalene
  7. 1-Undecanol
  8. Pentadecane

#### Column Profile



#### Engineered Self Cross-linking™ (ESC)

#### Phase Chemistry

- Proprietary

#### Recommended Applications

- Cannabis
- Terpenes
- Residual Solvents
- Volatile Amines
- EPA Method 8260
- EPA Method 524
- EPA Method 624
- Food
- Flavors and Fragrances
- Solvent Purity
- Alcohols



**ZB-624PLUS Test Mix**  
Part No.: [AG0-9203](#)




Contact Phenomenex or your local Phenomenex distributor for additional GC products and applications.

## ZB-624PLUS™ (cont'd)

### Improved Peak Shape of Volatile Amines

Volatile amines are challenging analytes for GC analysis. They can adsorb to even the smallest imperfections in fused silica. ZB-624PLUS undergoes a superior deactivation process which minimizes active compound adsorption leading to gains in peak response and shape.

#### Comparison of Volatile Amines on Various 624 Columns




**Volatile Amines on a Zebron ZB-624PLUS - 500 ppm**  
**Volatile Amines on a Restek® Rxi®-624Sil MS - 500 ppm**  
**Volatile Amines on a Agilent® DB®-624UI Ultra Inert - 500 ppm**

#### Conditions for all separations:

- Column:** Zebron ZB-624PLUS  
Restek Rxi-624Sil MS  
Agilent DB-624UI Ultra Inert
- Dimensions:** 30 meter x 0.32 mm x 1.80 µm
- Injection:** Split 20:1 @ 200 °C, 1 µL
- Recommended Liner:** Zebron PLUS Straight Z-Liner™
- Liner Part No.:** AG2-0A03-05 (for Agilent® & Thermo Scientific® systems)
- Carrier Gas:** Helium @ 1.8 mL/min (constant flow)
- Oven Program:** 50 °C for 1 min, to 200 °C @ 20 °C/min for 5 min
- Detector:** FID @ 250 °C
- Sample:** 1. Isopropylamine  
2. Diethylamine  
3. Triethylamine

### Volatile Organic Compounds in EPA Methods 8260, 524, and 624

Our high efficiency dimension and superior deactivation can stand real world samples. In addition, MS certification provides extreme low bleed to your GC-MS analysis.



- Column:** Zebron ZB-624PLUS
- Dimensions:** 30 meter x 0.25 mm x 1.40 µm
- Part No.:** ZHG-G040-27
- Injection:** Split 50:1 @ 230 °C, 1 µL
- Recommended Liner:** Zebron PLUS Straight Z-Liner™
- Liner Part No.:** AG2-0A03-05 (for Agilent® & Thermo Scientific® systems)
- Carrier Gas:** Helium @ 0.7 mL/min (constant flow)
- Oven Program:** 40 °C for 2 min, to 210 °C @ 17 °C/min for 3 min
- Detection:** Mass Spec transfer line @ 250 °C

Comparative separations may not be representative of all applications.




## ZB-624PLUS™ (cont'd)


### Exceeding USP <467> System Suitability

USP <467> method requires resolution of 1.5 for critical pairs. Zebron ZB-624PLUS took the challenge and succeeded with even greater resolution!


#### Class 1 Residual Solvents




#### Class 2B Residual Solvents



#### Class 2A Residual Solvents



#### Class 1 and 2A Residual Solvents



#### Same conditions for all separations:

- Column: Zebron ZB-624PLUS
- Dimensions: 30 meter x 0.32 mm x 1.80 µm
- Part No.: [7HM-G040-31](#)
- Injection: Split 5:1 @ 140 °C, 1 µL
- Recommended Liner: Zebron PLUS Straight Z-Liner™
- Liner Part No.: [AG2-0A03-05](#) (for Agilent® & Thermo Scientific® systems)
- Carrier Gas: Helium @ 2.2 mL/min (constant flow)
- Oven Program: 40 °C for 20 min to 240 °C @ 10 °C/min
- Detector: FID @ 250 °C

#### Ordering Information

##### Zebron ZB-624PLUS GC Columns

ID(mm)	df(µm)	Temp. Limits °C	Part No.
<b>20-Meter</b>			
0.18	1.00	-20 to 300/320	<a href="#">7FD-G040-22</a>
<b>30-Meter</b>			
0.25	1.40	-20 to 300/320	<a href="#">7HG-G040-27</a>
0.32	1.80	-20 to 300/320	<a href="#">7HM-G040-31</a>
0.53	3.00	-20 to 300/320	<a href="#">7HK-G040-36</a>
<b>60-Meter</b>			
0.25	1.40	-20 to 300/320	<a href="#">7KG-G040-27</a>
0.32	1.80	-20 to 300/320	<a href="#">7KM-G040-31</a>
0.53	3.00	-20 to 300/320	<a href="#">7KK-G040-36</a>
<b>75-Meter</b>			
0.53	3.00	-20 to 300/320	<a href="#">7LK-G040-36</a>

Note: If you need a 5 in. cage, contact technical support at [www.phenomenex.com/chat](http://www.phenomenex.com/chat) or simply reach out to your Technical consultant. Conditions may apply. Agilent 6850, some SRI and process GC systems use only 5 in. cages. 0.18mm, 0.25mm, and 0.32mm IDs are MS certified.

## ZB-1HT Inferno™

### Robust Results up to 430 °C

- First non-metal columns stable to 430 °C
- Provides true boiling point separation for hydrocarbon distillation methods
- Longer lifetime with rugged high temperature, polyimide coated, fused silica tubing
- Low activity, provides good peak shape for acidic and basic samples
- Provides robust column performance for high temperature bake outs

Upgrade to Zebron from any 100% dimethylpolysiloxane phase:

**Agilent®**

- DB®-1
- DB-1ht
- HP-1
- CP-Sil 5 CB
- CP-SimDist

**Restek®**

- Rtx®-1
- Rxi®-1HT

**SGE®**


- BP1

**Supelco®**

- SPB®-1
- Petrocol® 2887

### Rugged, High-Temperature Performance

#### Great Separation of High Boiling Hydrocarbons (ASTM D6352)




**Column:** Zebron ZB-1HT Inferno  
**Dimensions:** 5 meter x 0.53 mm x 0.10 µm  
**Part No.:** [7AK-G014-02](#)  
**Injection:** On-Column @ 43 °C, 0.1 µL  
**Carrier Gas:** Helium @ 4.4 mL/min (constant flow)  
**Oven Program:** 40 °C for 0.5 min to 430 °C @ 20 °C/min for 10 min  
**Detector:** FID @ 430 °C

Sample:	1. C10	16. C23	31. C38	46. C66
	2. C11	17. C24	32. C39	47. C68
	3. C12	18. C25	33. C40	48. C70
	4. C13	19. C26	34. C42	49. C72
	5. C14	20. C27	35. C44	50. C74
	6. C15	21. C28	36. C46	51. C76
	7. C16	22. C29	37. C48	52. C78
	8. C17	23. C30	38. C50	53. C80
	9. Pristane	24. C31	39. C52	54. C82
	10. C18	25. C32	40. C54	55. C84
	11. Phytane	26. C33	41. C56	56. C86
	12. C19	27. C34	42. C58	57. C88
	13. C20	28. C35	43. C60	58. C90
	14. C21	29. C36	44. C62	
	15. C22	30. C37	45. C64	

Note: Sample was a combination of PolyWax® 655 and retention time markers C8-C40 in CS<sub>2</sub>/Chloroform

#### Bake Off Contaminants from Dirty Matrices



**Column:** Zebron ZB-1HT Inferno  
**Dimensions:** 30 meter x 0.25 mm x 0.10 µm  
**Part No.:** [7HG-G014-02](#)  
**Injection:** On-Column @ 153 °C, 1 µL  
**Carrier Gas:** Helium @ 1 mL/min (constant flow)  
**Oven Program:** 150 °C to 400 °C @ 14 °C/min for 6 min  
**Detector:** FID @ 400 °C  
**Sample:** Sample was 1% in dichloromethane Mobil® 1 10W-30 Fully Synthetic Motor Oil




2007 R&D 100 Award Recipient

## ZB-1HT Inferno™

### Run Versatile Samples

Hydrocarbons from Water by GC-FID  
DIN EN ISO 9377-2 (DEV H53)



**Column:** Zebron ZB-1HT Inferno  
**Dimensions:** 15 meter x 0.32 mm x 0.25 µm  
**Part No.:** [ZEM-G014-11](#)  
**Injection:** Splitless @ 300 °C, 20 µL  
**Carrier Gas:** Helium @ 2.0 mL/min (constant flow)  
**Oven Program:** 50 °C for 2 min to 320 °C @ 30 °C/min for 5 min  
**Detector:** FID @ 330 °C  
**Sample:** 1. Decane (C10)  
 2. Tetracontane (C40)


#### Ordering Information

##### Zebron ZB-1HT Inferno GC Columns

ID(mm)	df(µm)	Temp. Limits °C	Part No.
<b>5-Meter</b>			
0.53	0.10	-60 to 400	<a href="#">7AK-G014-02</a>
<b>10-Meter</b>			
0.32	0.25	-60 to 400/430	<a href="#">7CM-G014-11</a>
<b>15-Meter</b>			
0.25	0.10	-60 to 400/430	<a href="#">7EG-G014-02</a>
0.25	0.25	-60 to 400/430	<a href="#">7EG-G014-11</a>
0.32	0.10	-60 to 400/430	<a href="#">7EM-G014-02</a>
0.32	0.25	-60 to 400/430	<a href="#">7EM-G014-11</a>
0.53	0.15	-60 to 400	<a href="#">7EK-G014-05</a>
<b>20-Meter</b>			
0.18	0.18	-60 to 400/430	<a href="#">7FD-G014-08</a>
<b>30-Meter</b>			
0.25	0.10	-60 to 400/430	<a href="#">7HG-G014-02</a>
0.25	0.25	-60 to 400/430	<a href="#">7HG-G014-11</a>
0.32	0.10	-60 to 400/430	<a href="#">7HM-G014-02</a>
0.32	0.25	-60 to 400/430	<a href="#">7HM-G014-11</a>

Note: If you need a 5 in. cage, contact technical support at [www.phenomenex.com/chat](http://www.phenomenex.com/chat) or simply reach out to your Technical consultant. Conditions may apply. Agilent 6850, some SRI and process GC systems use only 5 in. cages.


#### Column Profile



\*0.53 mm ID columns are rated to 400 °C.

#### Engineered Self Cross-linking™ (ESC)

#### Phase Chemistry



100 % Dimethylpolysiloxane

#### Recommended Applications

- Diesel Fuel
- High Boiling Petroleum Products
- High Molecular Weight Waxes
- Hydrocarbons
- Motor Oils
- Polymers/Plastics
- Simulated Distillation



**ZB-1HT Test Mix**  
Part No.: [AGO-5155](#)



Contact Phenomenex or your local Phenomenex distributor for additional GC products and applications.



Extend column lifetime. Add a Z-Guard™ to your next Zebron GC order.

## ZB-5HT Inferno™


### Robust Results Up to 430 °C

- First non-metal columns stable to 430 °C
- Robust column for high temperature bake outs and analysis, such as biodiesel, long-chain hydrocarbons, polymers, and high molecular weight compounds
- Provides true boiling point separation for hydrocarbon distillation methods
- Longer lifetime with rugged high temperature, polyimide coated, fused silica tubing
- Low activity, provides good peak shape for acidic and basic samples

### Zebron Inferno Columns Win in the Lifetime Test

#### How does the lifetime test work?

All columns were held at 400 °C for 2 hours and then the oven was lowered to 120 °C for pentadecane analysis. The VF-5ht column broke just after 40 hours at 400 °C. The ZB-5HT had the same retention for pentadecane at 100 hours as the DB-5ht column at 40 hours — over 2X the lifetime!



Note that the VF-5ht column died around 40 hours at 400 °C whereas the Zebron ZB-5HT Inferno column maintained great retention of Pentadecane over 100 hours.

#### Conditions for all columns:


- Dimensions:** 30 meter x 0.25 mm x 0.10 μm
- Injection:** 1.0 μL of test mix [AG0-7578](#)
- Carrier Gas:** Helium @ 1.9 mL/min (constant flow)
- Oven Program:** 120 °C (Isothermal)
- Detector:** FID @ 400 °C
- Sample:** Pentadecane

Comparative separations may not be representative of all applications.

Upgrade to Zebron from any 5 % phenyl / 95 % dimethylpolysiloxane phase:

Agilent®	Restek®	SGE®
<ul style="list-style-type: none"> <li>• DB®-5ht</li> <li>• VF-5ht</li> </ul>	<ul style="list-style-type: none"> <li>• Rxi®-5HT</li> <li>• Stx®-5HT</li> <li>• XT®-5HT</li> <li>• Rtx®-5HT</li> </ul>	<ul style="list-style-type: none"> <li>• HT-5</li> </ul>

### Paraffin Wax




**Column:** Zebron ZB-5HT Inferno  
**Dimensions:** 15 meter x 0.32 mm x 0.10 μm  
**Part No.:** [7EM-G015-02](#)  
**Injection:** On Column @ 43 °C, 0.1 μL  
**Carrier Gas:** Helium @ 1.9 mL/min (constant flow)  
**Oven Program:** 40 °C for 2 min to 430 °C @ 20 °C for 10 min  
**Detector:** FID @ 430 °C  
**Sample:** Paraffin Wax

## ZB-5HT Inferno™

### Well-Suited for Fuels Analysis

#### Free Total Glycerin in B100 Biodiesel by GC-FID



**Column:** Zebron ZB-5HT Inferno  
**Dimensions:** 15 meter x 0.32 mm x 0.10 µm  
 + 2 meter x 0.53 mm Z-Guard™  
**Part No.:** [ZEM-G015-02](#)  
**Injection:** On-Column @ 53 °C, 1 µL  
**Carrier Gas:** Helium @ 3.0 mL/min (constant flow)  
**Oven Program:** 50 °C for 1 min to 180 °C @ 15 °C/min to 230 °C @ 7 °C/min to 380 °C @ 30 °C/min for 10 min  
**Detector:** FID @ 380 °C  
**Note:** A 2 m x 0.53 mm Guard Column was connected to the analytical column per ASTM method requirement  
**Sample:**

1. Glycerol	6. Tricarpin (ISTD2)
2. Butanetriol (ISTD1)	7. Diglycerides
3. Esters	8. 1,3-Diolein
4. Monoglycerides	9. Triglycerides
5. 1-Monooleoyl-rac-glycerol	10. Triolein


#### Ordering Information

##### Zebron ZB-5HT Inferno GC Columns

ID(mm)	df(µm)	Temp. Limits °C	Part No.
<b>10-Meter</b>			
0.32	0.10	-60 to 400/430	<a href="#">ZCM-G015-02</a>
<b>15-Meter</b>			
0.25	0.10	-60 to 400/430	<a href="#">ZEG-G015-02</a>
0.25	0.25	-60 to 400/430	<a href="#">ZEG-G015-11</a>
0.32	0.10	-60 to 400/430	<a href="#">ZEM-G015-02</a>
0.32	0.25	-60 to 400/430	<a href="#">ZEM-G015-11</a>
0.53	0.15	-60 to 400	<a href="#">ZEK-G015-05</a>
<b>15-Meter with 2-Meter Spliced Guard (0.53 mm ID)</b>			
0.32	0.10	-60 to 400/430	<a href="#">ZEM-G015-02-GST</a>
<b>20-Meter</b>			
0.18	0.18	-60 to 400/430	<a href="#">ZFD-G015-08</a>
<b>30-Meter</b>			
0.25	0.10	-60 to 400/430	<a href="#">ZHG-G015-02</a>
0.25	0.25	-60 to 400/430	<a href="#">ZHG-G015-11</a>
0.32	0.10	-60 to 400/430	<a href="#">ZHM-G015-02</a>
0.32	0.25	-60 to 400/430	<a href="#">ZHM-G015-11</a>
0.53	0.15	-60 to 400	<a href="#">ZHK-G015-05</a>
<b>60-Meter</b>			
0.25	0.25	-60 to 400/430	<a href="#">ZKG-G015-11</a>

Note: If you need a 5 in. cage, contact technical support at [www.phenomenex.com/chat](http://www.phenomenex.com/chat) or simply reach out to your Technical consultant. Conditions may apply. Agilent 6850, some SRI and process GC systems use only 5 in. cages.


#### Column Profile



\*0.53 mm ID columns are rated to 400 °C.

#### Engineered Self Cross-linking™ (ESC)


##### Phase Chemistry



5 % Phenyl 95 % Dimethylpolysiloxane

#### Recommended Applications

- Diesel Fuels
- High Boiling Petroleum Products
- High Molecular Weight Waxes
- Hydrocarbons
- Motor Oils
- Polymers/Plastics
- Simulated Distillation
- Surfactants
- Triglycerides



**ZB-5HT Test Mix**  
 Part No.: [AGO-5155](#)



Contact Phenomenex or your local Phenomenex distributor for additional GC products and applications.



## ZB-35HT Inferno™

### High Temperature Stability for Mid-Polarity

- First non-metal, 35% phenyl columns stable to 400 °C
- Longer lifetime with rugged high temperature, polyimide coated, fused silica tubing
- Robust column for high temperature analysis
- Great for high molecular weight compounds
- Eliminate carry-over with high temperature bake outs
- Low activity, provides good peak shape for acidic and basic samples

Upgrade to Zebron from any

35% phenyl / 65% dimethylpolysiloxane phase:

Agilent®	Restek®	SGE®	Supelco®	OV®
• DB®-35	• Rtx®-35	• BPX35	• MDN-35	• OV-11
• HP-35ms	• Rtx-35ms	• BPX608	• SPB®-35	
• HP-35			• SPB-608	

### Lower Bleed Than Other Columns!

Conditions for all columns:


**Dimensions:** 30 meter x 0.25 mm x 0.25 µm

**Injection:** Split 20:1 @ 200 °C, 1 µL

**Carrier Gas:** Helium @ 1.7 mL/min (constant flow)

**Oven Program:** 100 °C to 320 °C @ 30 °C/min for 5 min to 340 °C @ 20 °C/min for 5 min to 360 °C @ 20 °C/min for 5 min to 380 °C @ 20 °C/min for 5 min to 400 °C @ 20 °C/min for 5 min to 100 °C @ 30 °C/min for 8 min


**Detector:** FID @ 405 °C



### See What You've Been Missing

#### A) ZB-35HT Inferno


30 meter x 0.25 mm x 0.25 µm



VS.

#### B) Restek Rtx-35

30 meter x 0.25 mm x 1.00 µm



**Column:** As listed

**Dimensions:** As listed

**Part No.:** 7HG-G025-11 (ZB-35HT Inferno)

**Injection:** A) Split 50:1 @ 350 °C, 1 µL

B) Split 50:1 @ 300 °C, 1 µL

**Carrier Gas:** Helium @ 2.1 mL/min (constant flow)

**Oven Program:** A) 140 °C to 400 °C @ 10 °C/min

B) 140 °C to 300 °C @ 10 °C/min

**Detector:** A) FID @ 400 °C

B) FID @ 320 °C


**Sample:** 1. Hexadecylamine

**Note:** Chromatogram is courtesy of Northeastern Chemical Company.

## ZB-35HT Inferno™

### Well-Suited for Environmental Contaminants

#### PAHs and PCBs in a Single Run




**Column:** Zebron ZB-35HT Inferno  
**Dimensions:** 30 meter x 0.25 mm x 0.25 µm  
**Part No.:** [7HG-G025-11](#)  
**Injection:** Splitless @ 265 °C, 2 µL  
**Carrier Gas:** Helium @ 1 mL/min (constant flow)  
**Oven Program:** 85 °C for 3 min to 320 °C @ 7 °C /min for 8 min  
**Detector:** MSD @ 280 °C  
**Sample:** Compounds are 5 ppm  
 1. Naphthalene      9. PCB 101      17. Chrysene  
 2. Acenaphthylene    10. Fluoranthene    18. PCB 194  
 3. Acenaphthene      11. Pyrene      19. Benzo[b]fluoranthene  
 4. Fluorene          12. PCB 118      20. Benzo[k]fluoranthene  
 5. Phenanthrene      13. PCB 153      21. Benzo[a]pyrene  
 6. Anthracene        14. PCB 138      22. Indeno[1,2,3-cd]pyrene  
 7. PCB 28            15. PCB 180      23. Dibenz[a,h]anthracene  
 8. PCB 52            16. Benz[a]anthracene    24. Benzo[g,h,i]perylene

#### Ordering Information

Zebron ZB-35HT Inferno GC Columns			
ID(mm)	df(µm)	Temp. Limits °C	Part No.
<b>15-Meter</b>			
0.25	0.10	40 to 400	<a href="#">7EG-G025-02</a>
0.25	0.25	40 to 400	<a href="#">7EG-G025-11</a>
0.32	0.25	40 to 400	<a href="#">7EM-G025-11</a>
<b>20-Meter</b>			
0.18	0.18	40 to 400	<a href="#">7FD-G025-08</a>
<b>30-Meter</b>			
0.25	0.10	40 to 400	<a href="#">7HG-G025-02</a>
0.25	0.25	40 to 400	<a href="#">7HG-G025-11</a>
0.32	0.25	40 to 400	<a href="#">7HM-G025-11</a>


Note: If you need a 5 in. cage, contact technical support at [www.phenomenex.com/chat](http://www.phenomenex.com/chat) or simply reach out to your Technical consultant. Conditions may apply. Agilent 6850, some SRI and process GC systems use only 5 in. cages.

#### Column Profile



#### Engineered Self Cross-linking™ (ESC)

#### Phase Chemistry



35 % Phenyl      65 % Dimethylpolysiloxane

#### Recommended Applications

- Amines
- Chemicals
- Drugs
- EPA Methods (508, 608, 8081, 8141, 8151)
- PCBs / Aroclors
- Pesticides
- Pharmaceuticals
- Steroids



Contact Phenomenex or your local Phenomenex distributor for additional GC products and applications.


Extend column lifetime. Add a Z-Guard™ to your next Zebron GC order.

## ZB-XLB-HT Inferno™

### High Temp Stability, Low Bleed

- Rugged, non-metal si-arylene GC column stable to 400 °C
- Robust column for high temperature bake outs and analysis, such as high molecular weight compounds
- Provides unique selectivity for conformational analyses
- Longer lifetime with high temperature, polyimide coated, fused silica tubing
- Low activity, provides good peak shape for acidic and basic samples
- Good tool for general screening to identify unknown samples

#### Aroclor 1242: DIN Method 51527



**Column:** Zebron ZB-XLB-HT Inferno  
**Dimensions:** 30 meter x 0.25 mm x 0.25 µm  
**Part No.:** 7HG-G024-11  
**Injection:** Split 2:1 @ 250 °C, 1 µL, pressure pulse @ 40 psi for first 0.25 min  
**Carrier Gas:** Helium @ 1.5 mL/min (constant flow)  
**Oven Program:** 50 °C for 0.5 min to 210 °C @ 40 °C/min for 3 min to 230 °C @ 30 °C/min for 5 min to 250 °C @ 30 °C/min for 5 min to 320 °C @ 40 °C/min for 2 min  
**Detector:** ECD @ 350 °C  
**Sample:** Total concentration of aroclors was 90 ppm in isooctane

1. TCMX	12. BZ# 138
2. BZ# 31	13. BZ# 126
3. BZ# 28	14. BZ# 167
4. BZ# 52	15. BZ# 156
5. BZ# 101	16. BZ# 180
6. BZ# 77	17. BZ# 157
7. BZ# 123	18. BZ# 170
8. BZ# 118	19. BZ# 169
9. BZ# 153	20. BZ# 189
10. BZ# 114	21. DCB
11. BZ# 105	

Upgrade to Zebron from these similar\* phases:

**Agilent®**

- DB®-XLB
- VF-XMS

**Restek®**


- DB®-XLB

**Supelco®**

- MDN-12

\*not exact equivalent, selectivity may differ

#### Melamine and Cyanuric Acid by GC-MS




**Column:** Zebron ZB-XLB-HT Inferno  
**Dimensions:** 15 meter x 0.25 mm x 0.25 µm  
**Part No.:** 7EG-G024-11  
**Injection:** On-Column @ 103 °C, 1 µL  
**Carrier Gas:** Helium @ 1.4 mL/min (constant flow)  
**Oven Program:** 100 °C for 0.5 min to 320 °C @ 25 °C/min  
**Detector:** MSD @ 325 °C  
**Sample:** Analytes are 200 ng / 100 µL in BSTFA / Pyridine (1:1)

1. Cyanuric Acid 13C3 (IS)
2. Cyanuric Acid
3. Melamine 13C3 15N3 (IS)
4. Melamine

## ZB-XLB-HT Inferno™

### Good Results for Difficult Samples

#### Explosives by GC-MS



- Column:** Zebron ZB-XLB-HT Inferno  
**Dimensions:** 15 meter x 0.25 mm x 0.25 µm  
**Part No.:** [7EG-G024-11](#)  
**Injection:** On-Column @ 73 °C, 0.5 µL  
**Carrier Gas:** Helium @ 1.4 mL/min (constant flow)  
**Oven Program:** 70 °C for 1 min to 140 °C @ 25 °C/min for 4 min to 280 °C @ 25 °C/min  
**Detector:** MSD @ 300 °C, 40-400 amu  
**Sample:** Analytes are 10 ppm in dichloromethane
- |                                 |                                       |
|---------------------------------|---------------------------------------|
| 1. Nitrobenzene                 | 9. 2,4,6-Trinitrotoluene (2,4,6-TNT)  |
| 2. 2-Nitrotoluene               | 10. PETN                              |
| 3. 3-Nitrotoluene               | 11. 1,3,5-Trinitrobenzene (1,3,5-TNB) |
| 4. 4-Nitrotoluene               | 12. RDX                               |
| 5. Nitroglycerin                | 13. 4-Amino-2,6-dinitrotoluene        |
| 6. 2,6-Dinitrotoluene (2,6-DNT) | 14. 3,5-Nitroaniline                  |
| 7. 1,3-Dinitrobenzene (1,3-DNB) | 15. 2-Amino-4,6-dinitrotoluene        |
| 8. 2,4-Dinitrotoluene           | 16. Tetryl                            |

#### Ordering Information

Zebron ZB-XLB-HT Inferno GC Columns			
ID(mm)	df(µm)	Temp. Limits °C	Part No.
<b>15-Meter</b>			
0.25	0.10	30 to 400	<a href="#">7EG-G024-02</a>
0.25	0.25	30 to 400	<a href="#">7EG-G024-11</a>
<b>20-Meter</b>			
0.18	0.18	30 to 400	<a href="#">7FD-G024-08</a>
<b>30-Meter</b>			
0.25	0.10	30 to 400	<a href="#">7HG-G024-02</a>
0.25	0.25	30 to 400	<a href="#">7HG-G024-11</a>
0.32	0.25	30 to 400	<a href="#">7HM-G024-11</a>
0.32	0.10	30 to 400	<a href="#">7HM-G024-02</a>
<b>60-Meter</b>			
0.25	0.25	30 to 400	<a href="#">7KG-G024-11</a>

Note: If you need a 5 in. cage, contact technical support at [www.phenomenex.com/chat](http://www.phenomenex.com/chat) or simply reach out to your Technical consultant. Conditions may apply. Agilent 6850, some SRI and process GC systems use only 5 in. cages.

**Column Profile**

**Engineered Self Cross-linking™ (ESC)**

**Phase Chemistry**

- Proprietary

**Recommended Applications**

- Herbicides / Insecticides
- PCBs
- Pesticides
- Unknown Samples

**ZB-XLB-HT Test Mix**

Part No.: [AGO-7578](#)

Contact Phenomenex or your local Phenomenex distributor for additional GC products and applications.

Extend column lifetime. Add a Z-Guard™ to your next Zebron GC order.

## ZB-1

### Low Polarity for Versatile Applications

- Low polarity phase suited for true boiling point compounds
- Low bleed (MS Certified), low activity, and high efficiency
- Excellent resolving power of critical pairs in complex petrochemical samples
- Used for “fingerprinting” and routine quality control analyses

Upgrade to Zebron from any 100% dimethylpolysiloxane phase:

#### Agilent®

- DB®-1
- DB-2887
- DB-1 EVDX
- HP-1
- HP-101
- HP-PONA
- Ultra 1
- CP-Sil 5 CB

#### Restek®

- Rtx®-1
- Rtx-1PONA
- Rtx-1 F&F

#### SGE®

- BP1
- BP1-PONA
- BPX1-SimD

#### Supelco®

- SPB®-1
- SPB-1 TG
- SE-30
- MET-1
- SPB-1 Sulfur
- SPB-HAP


### Ordering Information

#### Zebron ZB-1 GC Columns

ID(mm)	df(µm)	Temp. Limits °C	Part No.
<b>10-Meter</b>			
0.53	2.65	-60 to 340/360	<a href="#">7CK-G001-35</a>
<b>15-Meter</b>			
0.25	0.10	-60 to 360/370	<a href="#">7EG-G001-02</a>
0.25	0.25	-60 to 360/370	<a href="#">7EG-G001-11</a>
0.25	1.00	-60 to 340/360	<a href="#">7EG-G001-22</a>
0.32	0.25	-60 to 360/370	<a href="#">7EM-G001-11</a>
0.32	1.00	-60 to 340/360	<a href="#">7EM-G001-22</a>
0.53	0.15	-60 to 360/370	<a href="#">7EK-G001-05</a>
0.53	0.50	-60 to 360/370	<a href="#">7EK-G001-17</a>
0.53	1.50	-60 to 340/360	<a href="#">7EK-G001-28</a>
<b>30-Meter</b>			
0.25	0.10	-60 to 360/370	<a href="#">7HG-G001-02</a>
0.25	0.25	-60 to 360/370	<a href="#">7HG-G001-11</a>
0.25	0.50	-60 to 360/370	<a href="#">7HG-G001-17</a>
0.25	1.00	-60 to 340/360	<a href="#">7HG-G001-22</a>
0.32	0.25	-60 to 360/370	<a href="#">7HM-G001-11</a>
0.32	0.50	-60 to 360/370	<a href="#">7HM-G001-17</a>
0.32	1.00	-60 to 340/360	<a href="#">7HM-G001-22</a>
0.32	3.00	-60 to 340/360	<a href="#">7HM-G001-36</a>
0.32	5.00	-60 to 340/360	<a href="#">7HM-G001-39</a>
0.53	0.50	-60 to 360/370	<a href="#">7HK-G001-17</a>
0.53	1.50	-60 to 340/360	<a href="#">7HK-G001-28</a>
0.53	3.00	-60 to 340/360	<a href="#">7HK-G001-36</a>
0.53	5.00	-60 to 340/360	<a href="#">7HK-G001-39</a>
<b>50-Meter</b>			
0.25	0.50	-60 to 360/370	<a href="#">7JG-G001-17</a>
<b>60-Meter</b>			
0.25	0.25	-60 to 360/370	<a href="#">7KG-G001-11</a>
0.25	1.00	-60 to 340/360	<a href="#">7KG-G001-22</a>
0.32	0.25	-60 to 360/370	<a href="#">7KM-G001-11</a>
0.32	1.00	-60 to 340/360	<a href="#">7KM-G001-22</a>
0.32	3.00	-60 to 340/360	<a href="#">7KM-G001-36</a>
0.53	1.50	-60 to 340/360	<a href="#">7KK-G001-28</a>
<b>100-Meter</b>			
0.25	0.50	-60 to 360/370	<a href="#">7MG-G001-17</a>

Note: If you need a 5 in. cage, contact technical support at [www.phenomenex.com/chat](http://www.phenomenex.com/chat) or simply reach out to your Technical consultant. Conditions may apply. Agilent 6850, some SRI and process GC systems use only 5 in. cages.


### Column Profile



\*Thicker films (≥ 1.0 µm) are rated to 340/360 °C.

### Engineered Self Cross-linking™ (ESC)

### Phase Chemistry



100% Dimethylpolysiloxane

### Recommended Applications

- Ethanol
- Hydrocarbons
- Mercaptans
- MTBE
- Natural Gas Odorants
- Oxygenates and GROs
- Solvent Impurities
- Sulfur Compounds (Light)



**ZB-1 Test Mix**  
Part No.: [AGO-5155](#)



Zebron GC Columns MS Certification, see p. 437



Contact Phenomenex or your local Phenomenex distributor for additional GC products and applications.



Extend column lifetime. Add a Z-Guard™ to your next Zebron GC order.



## ZB-5

### Low Polarity for a Wide Application Range

- Rugged, versatile low polarity column for general lab purpose
- Resilient to dirty samples—long column life
- Low bleed (MS Certified) especially suited to high sensitivity work using GC-MS
- Extremely inert for active compounds such as drugs or pesticides
- Great column for unknown samples

Upgrade to Zebron from any 5% phenyl / 95% dimethylpolysiloxane phase:


Agilent®	Restek®	SGE®	Supelco®	OV®
<ul style="list-style-type: none"> <li>• DB®-5</li> <li>• HP-5</li> <li>• HP-PAS-5</li> <li>• CP-Sil 8 CB</li> <li>• Ultra 2</li> </ul>	<ul style="list-style-type: none"> <li>• Rtx®-5</li> </ul>	<ul style="list-style-type: none"> <li>• BP5</li> <li>• BPX5</li> </ul>	<ul style="list-style-type: none"> <li>• MDN-5</li> <li>• SPB®-5</li> <li>• PTE-5</li> <li>• SE-54</li> <li>• PTA-5</li> <li>• Equity®-5</li> <li>• Sac-5</li> </ul>	<ul style="list-style-type: none"> <li>• OV-5</li> </ul>

### Ordering Information

Zebron ZB-5 GC Columns			
ID(mm)	df(µm)	Temp. Limits °C	Part No.
<b>15-Meter</b>			
0.25	0.10	-60 to 360/370	<a href="#">7EG-G002-02</a>
0.25	0.25	-60 to 360/370	<a href="#">7EG-G002-11</a>
0.25	0.50	-60 to 360/370	<a href="#">7EG-G002-17</a>
0.25	1.00	-60 to 340/360	<a href="#">7EG-G002-22</a>
0.32	0.25	-60 to 360/370	<a href="#">7EM-G002-11</a>
0.32	1.00	-60 to 340/360	<a href="#">7EM-G002-22</a>
0.53	0.50	-60 to 360/370	<a href="#">7EK-G002-17</a>
0.53	1.50	-60 to 340/360	<a href="#">7EK-G002-28</a>
0.53	3.00	-60 to 340/360	<a href="#">7EK-G002-36</a>
<b>20-Meter</b>			
0.18	0.18	-60 to 360/370	<a href="#">7FD-G002-08</a>
<b>30-Meter</b>			
0.25	0.10	-60 to 360/370	<a href="#">7HG-G002-02</a>
0.25	0.25	-60 to 360/370	<a href="#">7HG-G002-11</a>
0.25	0.50	-60 to 360/370	<a href="#">7HG-G002-17</a>
0.25	1.00	-60 to 340/360	<a href="#">7HG-G002-22</a>
0.32	0.25	-60 to 360/370	<a href="#">7HM-G002-11</a>
0.32	0.50	-60 to 360/370	<a href="#">7HM-G002-17</a>
0.32	1.00	-60 to 340/360	<a href="#">7HM-G002-22</a>
0.53	0.50	-60 to 360/370	<a href="#">7HK-G002-17</a>
0.53	1.50	-60 to 340/360	<a href="#">7HK-G002-28</a>
0.53	3.00	-60 to 340/360	<a href="#">7HK-G002-36</a>
0.53	5.00	-60 to 340/360	<a href="#">7HK-G002-39</a>
<b>60-Meter</b>			
0.25	0.10	-60 to 360/370	<a href="#">7KG-G002-02</a>
0.25	0.25	-60 to 360/370	<a href="#">7KG-G002-11</a>
0.25	1.00	-60 to 340/360	<a href="#">7KG-G002-22</a>
0.32	0.25	-60 to 360/370	<a href="#">7KM-G002-11</a>
0.32	1.00	-60 to 340/360	<a href="#">7KM-G002-22</a>
0.53	1.50	-60 to 340/360	<a href="#">7KK-G002-28</a>

Note: If you need a 5 in. cage, contact technical support at [www.phenomenex.com/chat](http://www.phenomenex.com/chat) or simply reach out to your Technical consultant. Conditions may apply. Agilent 6850, some SRI and process GC systems use only 5 in. cages.


### Column Profile



\*Thicker films (≥ 1.0 µm) are rated to 340/360 °C.

### Engineered Self Cross-linking™ (ESC)

#### Phase Chemistry



5 % Phenyl 95 % Dimethylpolysiloxane

#### Recommended Applications

- Alkaloids
- Dioxins
- Drugs
- Essential Oils
- Flavors
- FAMES
- Halo-Hydrocarbons
- Herbicides
- PCBs / Aroclors
- Pesticides
- Phenols
- Residual Solvents



For ultra low bleed, consider using a ZB-5ms, see p. 154  
For high temperature analysis, consider using a ZB-5HT, see p. 146



Contact Phenomenex or your local Phenomenex distributor for additional GC products and applications.



Extend column lifetime. Add a Z-Guard™ to your next Zebron GC order.

## ZB-5ms

### Robust Results, Versatile Performance

- Popular rugged column for general purpose use
- Fully conditioned within 35 minutes
- High response for acids and bases
- Enhanced resolution of polyaromatic hydrocarbons (PAHs) and other multi-ring aromatic compounds

Upgrade to Zebron from any 5% phenyl-arylene / 95% dimethylpolysiloxane phase:

#### Agilent®

- DB®-5ms
- DB-5.625
- DB-5ms EVDX
- CP-Sil 8 CB MS
- VF-5ms

#### Restek®


- Rtx®-5Sil MS
- Rxi®-5Sil MS

#### Supelco®

- SLB®-5ms

### Long Lifetime

Consistent response after more than 700 samples at pH 2!




#### Ordering Information

##### Zebron ZB-5ms GC Columns

ID (mm)	df (µm)	Temp. Limits °C	Part No.
<b>10-Meter</b>			
0.10	0.10	-60 to 325/350	<a href="#">7CB-G010-02</a>
0.18	0.18	-60 to 325/350	<a href="#">7CD-G010-08</a>
<b>12-Meter</b>			
0.20	0.33	-60 to 325/350	<a href="#">7DE-G010-14</a>
<b>15-Meter</b>			
0.25	0.25	-60 to 325/350	<a href="#">7EG-G010-11</a>
<b>20-Meter</b>			
0.18	0.18	-60 to 325/350	<a href="#">7FD-G010-08</a>
0.18	0.32	-60 to 325/350	<a href="#">7FD-G010-51</a>
<b>25-Meter</b>			
0.20	0.33	-60 to 325/350	<a href="#">7GE-G010-14</a>
<b>30-Meter</b>			
0.25	0.25	-60 to 325/350	<a href="#">7HG-G010-11</a>
0.25	0.50	-60 to 325/350	<a href="#">7HG-G010-17</a>
0.25	1.00	-60 to 325/350	<a href="#">7HG-G010-22</a>
0.32	0.25	-60 to 325/350	<a href="#">7HM-G010-11</a>
0.32	0.50	-60 to 325/350	<a href="#">7HM-G010-17</a>
0.32	1.00	-60 to 325/350	<a href="#">7HM-G010-22</a>
<b>60-Meter</b>			
0.25	0.10	-60 to 325/350	<a href="#">7KG-G010-02</a>
0.25	0.25	-60 to 325/350	<a href="#">7KG-G010-11</a>
0.32	0.25	-60 to 325/350	<a href="#">7KM-G010-11</a>

Note: If you need a 5 in. cage, contact technical support at [www.phenomenex.com/chat](http://www.phenomenex.com/chat) or simply reach out to your Technical consultant. Conditions may apply. Agilent 6850, some SRI and process GC systems use only 5 in. cages.


#### Column Profile



#### Engineered Self Cross-linking™ (ESC)

#### Phase Chemistry


5 % Phenyl-Arylene



95 % Dimethylpolysiloxane

#### Recommended Applications

- Acids
- Alkaloids
- Amines
- Dioxins
- Drugs
- Essential Oils
- Flavors
- FAMES
- Halo-hydrocarbons
- Herbicides
- PCBs/Aroclors
- Pesticides
- Phenols
- Residual Solvents
- Solvent Impurities



Contact Phenomenex or your local Phenomenex distributor for additional GC products and applications.

Extend column lifetime. Add a Z-Guard™ to your next Zebron GC order.

## ZB-35


### Intermediate Polarity for GC-MS

- Intermediate polarity column with temperature limits up to 360 °C allows high molecular weight analysis
- Excellent inertness to minimize analyte adsorption, improve efficiency, and reproducibility
- More rugged (longer column life) than other polar phases
- Excellent for trace analysis with bleed-sensitive detectors (MS, FID, ECD, NPD)

Upgrade to Zebron from any 35% phenyl / 65% dimethylpolysiloxane phase:

Agilent®	Restek®	SGE®	Supelco®	OV®
• DB®-35	• Rtx®-35	• BPX35	• MDN-35	• OV-11
• DB-35ms	• Rtx-35ms	• BPX608	• SPB®-35	
• HP-35			• SPB-608	
• HP-35ms				

### Common Drug Screen by GC-FID



**Sample:** All analytes are 25 ppm except nicotine at 100 ppm

1. Nicotine	16. Caffeine
2. Ibuprofen	17. Chlorpheniramine
3. Allobarbitol	18. Methapyrilene
4. Acetaminophen	19. Phenobarbital
5. Aprobital	20. Procaine
6. Butalbitol	21. Brompheniramine
7. Amobarbital	22. Chlorcyclizine
8. Pentobarbital	23. Cocaine
9. Phenacetin	24. Benactyzine
10. Secobarbital	25. Codeine
11. Benzphetamine	26. Diazepam
12. Meprobamate	27. Morphine
13. Dimenhydrinate	28. Hydrocodone
14. Hexobarbital	29. Oxymorphone
15. Doxylamine	30. Heroin


**Column:** Zebron ZB-35  
**Dimensions:** 30 meter x 0.25 mm x 0.25 µm  
**Part No.:** [7HG-G003-11](#)  
**Injection:** Split 10:1 @ 225 °C, 1.5 µL  
**Carrier Gas:** Helium @ 1.4 mL/min (constant flow)  
**Oven Program:** 120 °C to 180 °C @ 25 °C/min to 200 °C @ 6 °C/min to 300 °C @ 20 °C/min for 3 min  
**Detector:** FID @ 300 °C

### Ordering Information

Zebron ZB-35 GC Columns			
ID(mm)	df(µm)	Temp. Limits °C	Part No.
<b>10-Meter</b>			
0.10	0.10	40 to 340/360	<a href="#">7CB-G003-02</a>
<b>15-Meter</b>			
0.25	0.25	40 to 340/360	<a href="#">7EG-G003-11</a>
0.25	0.50	40 to 340/360	<a href="#">7EG-G003-17</a>
0.53	1.00	40 to 340/360	<a href="#">7EK-G003-22</a>
<b>30-Meter</b>			
0.25	0.25	40 to 340/360	<a href="#">7HG-G003-11</a>
0.25	0.50	40 to 340/360	<a href="#">7HG-G003-17</a>
0.32	0.25	40 to 340/360	<a href="#">7HM-G003-11</a>
0.32	0.50	40 to 340/360	<a href="#">7HM-G003-17</a>
0.53	0.50	40 to 340/360	<a href="#">7HK-G003-17</a>
0.53	1.00	40 to 340/360	<a href="#">7HK-G003-22</a>
<b>60-Meter</b>			
0.25	0.25	40 to 340/360	<a href="#">7KG-G003-11</a>
0.32	0.25	40 to 340/360	<a href="#">7KM-G003-11</a>


Note: If you need a 5 in. cage, contact technical support at [www.phenomenex.com/chat](http://www.phenomenex.com/chat) or simply reach out to your Technical consultant. Conditions may apply. Agilent 6850, some SRI and process GC systems use only 5 in. cages.

### Column Profile



Engineered Self Cross-linking™ (ESC)

### Phase Chemistry



35 % Phenyl 65 % Dimethylpolysiloxane

### Recommended Applications

- Amines
- Drugs
- EPA Methods (508, 608, 8081, 8141, 8151)
- PCBs / Aroclors
- Pesticides
- Pharmaceuticals



**ZB-35 Test Mix**  
 Part No.: [AGO-5156](#)



Contact Phenomenex or your local Phenomenex distributor for additional GC products and applications.



Extend column lifetime. Add a Z-Guard™ to your next Zebron GC order.

## ZB-50


### Robust Results, Rugged Performance

- High polarity column with temperature limits up to 340 °C allows high temperature bake out to remove contaminants
- Excellent inertness to minimize analyte adsorption, improve efficiency, and reproducibility
- More rugged (longer column life) than other polar phases
- Excellent for trace analysis with bleed-sensitive detectors
- Great for drug screening and environmental compounds

Upgrade to Zebron from any 50% phenyl / 50% dimethylpolysiloxane phase:

Agilent®	Restek®	SGE®	Supelco®
<ul style="list-style-type: none"> <li>• DB®-17</li> <li>• DB-17ht</li> <li>• DB-17ms</li> <li>• DB-17 EVDX</li> </ul>	<ul style="list-style-type: none"> <li>• Rtx®-50</li> </ul>	<ul style="list-style-type: none"> <li>• BPX50</li> </ul>	<ul style="list-style-type: none"> <li>• SP®-2250</li> <li>• SPB®-17</li> <li>• SPB-50</li> </ul>

### Antihistamine by GC-FID



**Column:** Zebron ZB-50  
**Dimensions:** 30 meter x 0.32 mm x 0.50 µm  
**Part No.:** [7HM-G004-17](#)  
**Injection:** Split 20:1 @ 250 °C, 1 µL  
**Carrier Gas:** Helium @ 40 cm/sec (constant flow)  
**Oven Program:** 190 °C to 260 °C @ 25 °C/min for 10 min  
**Detector:** FID @ 270 °C  
**Sample:** Brompheniramine


### Ordering Information

#### Zebron ZB-50 GC Columns

ID (mm)	df (µm)	Temp. Limits °C	Part No.
<b>10-Meter</b>			
0.10	0.10	40 to 320/340	<a href="#">7CB-G004-02</a>
0.53	2.00	40 to 320/340	<a href="#">7CK-G004-32</a>
<b>15-Meter</b>			
0.25	0.15	40 to 320/340	<a href="#">7EG-G004-05</a>
0.25	0.25	40 to 320/340	<a href="#">7EG-G004-11</a>
0.53	1.00	40 to 320/340	<a href="#">7EK-G004-22</a>
<b>30-Meter</b>			
0.25	0.25	40 to 320/340	<a href="#">7HG-G004-11</a>
0.25	0.50	40 to 320/340	<a href="#">7HG-G004-17</a>
0.32	0.25	40 to 320/340	<a href="#">7HM-G004-11</a>
0.32	0.50	40 to 320/340	<a href="#">7HM-G004-17</a>
0.53	1.00	40 to 320/340	<a href="#">7HK-G004-22</a>
<b>60-Meter</b>			
0.25	0.25	40 to 320/340	<a href="#">7KG-G004-11</a>
0.25	0.50	40 to 320/340	<a href="#">7KG-G004-17</a>


Note: If you need a 5 in. cage, contact technical support at [www.phenomenex.com/chat](http://www.phenomenex.com/chat) or simply reach out to your Technical consultant. Conditions may apply. Agilent 6850, some SRI and process GC systems use only 5 in. cages.

### Column Profile



### Engineered Self Cross-linking™ (ESC)

### Phase Chemistry



50 % Phenyl      50 % Dimethylpolysiloxane

### Recommended Applications

- Antidepressants
- Cholesterols
- Drugs of Abuse
- EPA Methods (508, 608, 8081, 8141, 8151)
- Glycols
- Herbicides
- Pesticides
- Steroids
- Triglycerides



Contact Phenomenex or your local Phenomenex distributor for additional GC products and applications.



Extend column lifetime. Add a Z-Guard™ to your next Zebron GC order.

## ZB-624

### Robust Results for VOCs and Residual Solvents


- Widely used phase to separate volatile organic flavor and fragrance additives and residual solvents in industrial or pharmaceutical products (OVIs)
- Popular choice for residual solvent testing
- Excellent for US EPA Methods 501.3, 502.2, 503.1, 524.2, 601, 602, 624, 8010, 8015, 8020, 8021, 8240, 8260
- Specifically designed for the separation of volatile organic compounds (VOCs)
- Increased temperature limit speeds run times and re-equilibration

Upgrade to Zebron from any 6% cyanopropylphenyl / 94% dimethylpolysiloxane phase:

Agilent®	Restek®	SGE®	Supelco®	OV®
• DB®-624	• Rtx®-624	• BPX624	• SPB®-624	• OV-624
• DB-1301	• Rtx-1301		• SPB-1301	
• DB-VRX	• Rtx-VMS			
• HP-VOC				
• CP-1301				
• CP-Select 624 CB				

### Good Performance for Pharmaceuticals

#### USP <467> Residual Solvents Procedure A — Class 1



**Column:** Zebron ZB-624  
**Dimensions:** 30 meter x 0.32 mm x 1.80 µm  
**Part No.:** [7HM-G005-31](#)  
**Injection:** Split 5:1 @ 140 °C, 1 mL  
**Carrier Gas:** Helium @ 35 cm/sec (constant flow)  
**Oven Program:** 40 °C for 20 min to 240 °C @ 10 °C/min for 20 min  
**Detector:** FID @ 250 °C


**Sample:** 1. 1,1-Dichloroethene  
 2. 1,1,1-Trichloroethane  
 3. Carbon tetrachloride  
 4. Benzene  
 5. 1,2-Dichloroethane

#### Ordering Information

Zebron ZB-624 GC Columns			
ID (mm)	df (µm)	Temp. Limits °C	Part No.
<b>20-Meter</b>			
0.18	1.00	-20 to 260	<a href="#">7FD-G005-22</a>
<b>30-Meter</b>			
0.25	1.40	-20 to 260	<a href="#">7HG-G005-27</a>
0.32	1.80	-20 to 260	<a href="#">7HM-G005-31</a>
0.53	3.00	-20 to 260	<a href="#">7HK-G005-36</a>
0.53	1.00	-20 to 260	<a href="#">7HK-G005-22</a>
<b>60-Meter</b>			
0.25	1.40	-20 to 260	<a href="#">7KG-G005-27</a>
0.32	1.80	-20 to 260	<a href="#">7KM-G005-31</a>
0.53	3.00	-20 to 260	<a href="#">7KK-G005-36</a>
<b>75-Meter</b>			
0.53	3.00	-20 to 260	<a href="#">7LK-G005-36</a>
<b>105-Meter</b>			
0.53	3.00	-20 to 260	<a href="#">7NK-G005-36</a>


Note: If you need a 5 in. cage, contact technical support at [www.phenomenex.com/chat](http://www.phenomenex.com/chat) or simply reach out to your Technical consultant. Conditions may apply. Agilent 6850, some SRI and process GC systems use only 5 in. cages.

#### Column Profile



#### Phase Chemistry

6% Cyanopropylphenyl



94% Dimethylpolysiloxane

#### Recommended Applications

- Pharmaceuticals
- Residual Solvents
- Volatile Organic Compounds (VOCs)
- EPA Methods (501.3, 502.2, 503.1, 524.2, 601, 602, 624, 8010, 8015, 8020, 8021, 8240, 8260)



**ZB-624 Test Mix**  
 Part No.: [AGO-5159](#)



Contact Phenomenex or your local Phenomenex distributor for additional GC products and applications.



Extend column lifetime. Add a Z-Guard™ to your next Zebron GC order.



## ZB-1701

### Alternate Selectivity for Mid-Polarity Analyses


- Fast run and re-equilibration times for enhanced sample throughput and productivity
- Provides alternate selectivity to phenyl phases with similar polarity

Upgrade to Zebron from any 14% cyanopropylphenyl / 86% dimethylpolysiloxane phase:

Agilent®	Restek®	SGE®	Supelco®	OV®
• DB®-1701	• Rtx®-1701	• BP10	• SPB®-1701	• OV-1701
• CP-Sil 19 CB	• Rtx-VMS		• Equity®-1701	

### Good Peak Shape for Active Analytes

#### Aromatic Amines by GC-MS



**Column:** Zebron ZB-1701  
**Dimensions:** 30 meter x 0.25 mm x 0.25 µm  
**Part No.:** ZHG-G006-11  
**Injection:** Split 15:1 @ 220 °C, 1 µL  
**Carrier Gas:** Helium @ 1.0 mL/min (constant flow)  
**Oven Program:** 60 °C for 1 min to 110 °C @ 30 °C/min to 135 °C @ 9 °C/min to 260 °C @ 30 °C/min for 2 min  
**Detector:** MSD @ 180 °C  
**Sample:** Analytes are at 1.58 mg/mL each


1. Piperidine	8. o-Toluidine
2. 2-Methylpiperidine	9. N,N-Dimethylaniline
3. Aniline	10. β-Phenylethylamine
4. Benzylamine	11. N-Ethylaniline
5. α-Phenylethylamine	12. 2,4-Dimethylaniline
6. N-Methylaniline	13. N,N-Diethylaniline
7. m-Toluidine	14. Dibenzylamine

#### Ordering Information

Zebron ZB-1701 GC Columns			
ID (mm)	df (µm)	Temp. Limits °C	Part No.
<b>15-Meter</b>			
0.25	0.25	-20 to 280/300	ZEG-G006-11
0.32	0.25	-20 to 280/300	ZEM-G006-11
<b>30-Meter</b>			
0.25	0.25	-20 to 280/300	ZHG-G006-11
0.25	1.00	-20 to 260/280	ZHG-G006-22
0.32	0.25	-20 to 280/300	ZHM-G006-11
0.32	1.00	-20 to 260/280	ZHM-G006-22
0.53	1.00	-20 to 260/280	ZHK-G006-22
<b>60-Meter</b>			
0.25	0.25	-20 to 280/300	ZKG-G006-11
0.32	0.25	-20 to 280/300	ZKM-G006-11

Note: If you need a 5 in. cage, contact technical support at [www.phenomenex.com/chat](http://www.phenomenex.com/chat) or simply reach out to your Technical consultant. Conditions may apply. Agilent 6850, some SRI and process GC systems use only 5 in. cages.


#### Column Profile



\*Thicker films (≥ 1.0 µm) are rated to 260/280 °C.

#### Phase Chemistry

14% Cyanopropylphenyl



86% Dimethylpolysiloxane

#### Recommended Applications

- Alcohols
- Amines
- Aromatic Hydrocarbons
- Drugs
- Esters
- PAHs
- PCBs
- Pharmaceutical Intermediates
- Phenols
- Solvents
- Steroids
- TMS Sugars
- Tranquilizers



**ZB-1701 Test Mix**  
**Part No.:** AGO-5156



For enhanced response to Endrin and DDT, consider using ZB-1701P, See p. 159  
 Contact Phenomenex or your local Phenomenex distributor for additional GC products and applications.



Extend column lifetime. Add a Z-Guard™ to your next Zebron GC order.

## ZB-1701P


### Enhanced Response for DDT and Endrin

- Specially tested to ensure response of DDT, Endrin, Endrin Aldehyde, and Endrin Ketone
- Fast run and re-equilibration times for enhanced sample throughput and productivity
- Guaranteed column for pesticide analysis

Upgrade to Zebron from any 14% cyanopropylphenyl / 86% dimethylpolysiloxane phase:

Agilent®	Restek®	SGE®	Supelco®	OV®
• DB®-1701	• Rtx®-1701	• BP10	• SPB®-1701	• OV-1701
• DB-1701P	• Rtx-VMS		• Equity®-1701	
• CP-Sil 19 CB				

### Chlorinated Pesticides by GC-ECD: EPA Method 8081



**Column:** Zebron ZB-1701P  
**Dimensions:** 30 meter x 0.25 mm x 0.25 µm  
**Part No.:** [7HG-G012-11](#)  
**Injection:** Splitless @ 240 °C, 1 µL  
**Carrier Gas:** Helium @ 1.6 mL/min (constant flow)  
**Oven Program:** 100 °C to 200 °C @ 25 °C/min to 240 °C @ 6 °C/min to 265 °C @ 20 °C/min for 5 min  
**Detector:** ECD @ 300 °C  
**Sample:** All compounds are 20 ppm

1. 1-Bromo-2-Nitrobenzene (IS)	9. δ-BHC	17. DDD
2. Tetrachloro-m-xylene (Surr)	10. Heptachlor Epoxide	18. Endosulfan II
3. α-BHC	11. Endosulfan I	19. DDT
4. Pentachloronitrobenzene (IS)	12. γ-Chlordane	20. Endrin Aldehyde
5. γ-BHC (Lindane)	13. α-Chlordane	21. Methoxychlor
6. Heptachlor	14. DDE	22. Endosulfan Sulfate
7. Aldrin	15. Dieldrin	23. Endrin Ketone
8. β-BHC	16. Endrin	24. Decachlorobiphenyl (Surr)


### Ordering Information

#### Zebron ZB-1701P GC Columns

ID(mm)	df(µm)	Temp. Limits °C	Part No.
30-Meter			
0.25	0.25	-20 to 280/300	<a href="#">7HG-G012-11</a>

Note: If you need a 5 in. cage, contact technical support at [www.phenomenex.com/chat](http://www.phenomenex.com/chat) or simply reach out to your Technical consultant. Conditions may apply. Agilent 6850, some SRI and process GC systems use only 5 in. cages.


### Column Profile



\*Thicker films (≥ 1.0 µm) are rated to 260/280 °C.

### Phase Chemistry

14% Cyanopropylphenyl



86% Dimethylpolysiloxane

### Recommended Applications

- Nitrogen Containing Pesticides
- Organochlorine Pesticides
- Organophosphorous Pesticides
- PCBs / Aroclors



**ZB-1701 Test Mix**  
**Part No.:** [AGO-5156](#)



Contact Phenomenex or your local Phenomenex distributor for additional GC products and applications.



Extend column lifetime. Add a Z-Guard™ to your next Zebron GC order.

## ZB-WAX

### PEG Versatility for Solvents, Acids, and Amines


- High polarity column with low bleed (MS certified) for improved results
- Highly stable, long lifetime
- Low activity for amines
- Bonded, solvent rinsible
- Excellent chromatography of polar complex mixtures
- Widely used for profiling and “fingerprinting”

Upgrade to Zebron from any polyethylene glycol phase:

Agilent®	Restek®	SGE®	Supelco®
<ul style="list-style-type: none"> <li>• DB®-WAXetr</li> <li>• HP-INNOWax</li> <li>• CP-Wax 57 CB</li> </ul>	<ul style="list-style-type: none"> <li>• Rtx®-WAX</li> <li>• Famewax</li> <li>• Stabilwax®-DB</li> </ul>	<ul style="list-style-type: none"> <li>• SolGel-WAX™</li> </ul>	<ul style="list-style-type: none"> <li>• Met-Wax</li> <li>• Omegawax</li> </ul>

### Performs for Industrial Chemicals

#### BTEX by GC-FID




<p><b>Column:</b> Zebron ZB-WAX  <b>Dimensions:</b> 30 meter x 0.32 mm x 0.50 µm  <b>Part No.:</b> <a href="#">7HM-G007-17</a>  <b>Injection:</b> Split 20:1 @ 250 °C, 0.2 µL  <b>Carrier Gas:</b> Helium @ 2 mL/min (constant flow)  <b>Oven Program:</b> 60 °C to 75 °C @ 15 °C/min to 90 °C @ 3 °C/min (hold 3 min)  <b>Detector:</b> FID @ 300 °C</p>	<p><b>Sample:</b></p> <table border="0"> <tr> <td>1. Pentane</td> <td>7. Ethylbenzene</td> </tr> <tr> <td>2. Heptane</td> <td>8. p-Xylene</td> </tr> <tr> <td>3. Solvent (methylene chloride)</td> <td>9. m-Xylene</td> </tr> <tr> <td>4. Benzene</td> <td>10. Dodecane</td> </tr> <tr> <td>5. Decane</td> <td>11. o-Xylene</td> </tr> <tr> <td>6. Toluene</td> <td></td> </tr> </table>	1. Pentane	7. Ethylbenzene	2. Heptane	8. p-Xylene	3. Solvent (methylene chloride)	9. m-Xylene	4. Benzene	10. Dodecane	5. Decane	11. o-Xylene	6. Toluene	
1. Pentane	7. Ethylbenzene												
2. Heptane	8. p-Xylene												
3. Solvent (methylene chloride)	9. m-Xylene												
4. Benzene	10. Dodecane												
5. Decane	11. o-Xylene												
6. Toluene													

#### Ordering Information


Zebron ZB-WAX GC Columns			
ID (mm)	df (µm)	Temp. Limits °C	Part No.
<b>10-Meter</b>			
0.10	0.10	20 to 250/260	<a href="#">7CB-G007-02</a>
<b>15-Meter</b>			
0.25	0.25	20 to 250/260	<a href="#">7EG-G007-11</a>
0.32	0.25	20 to 250/260	<a href="#">7EM-G007-11</a>
0.53	1.00	20 to 250/260	<a href="#">7EK-G007-22</a>
<b>20-Meter</b>			
0.18	0.18	20 to 250/260	<a href="#">7FD-G007-08</a>
<b>30-Meter</b>			
0.25	0.15	20 to 250/260	<a href="#">7HG-G007-05</a>
0.25	0.25	20 to 250/260	<a href="#">7HG-G007-11</a>
0.25	0.50	20 to 250/260	<a href="#">7HG-G007-17</a>
0.25	1.00	20 to 250/260	<a href="#">7HG-G007-22</a>
0.32	0.15	20 to 250/260	<a href="#">7HM-G007-05</a>
0.32	0.25	20 to 250/260	<a href="#">7HM-G007-11</a>
0.32	0.50	20 to 250/260	<a href="#">7HM-G007-17</a>
0.53	0.50	20 to 250/260	<a href="#">7HK-G007-17</a>
0.53	1.00	20 to 250/260	<a href="#">7HK-G007-22</a>
<b>60-Meter</b>			
0.25	0.25	20 to 250/260	<a href="#">7KG-G007-11</a>
0.25	0.50	20 to 250/260	<a href="#">7KG-G007-17</a>
0.32	0.25	20 to 250/260	<a href="#">7KM-G007-11</a>
0.32	0.50	20 to 250/260	<a href="#">7KM-G007-17</a>
0.53	1.00	20 to 250/260	<a href="#">7KK-G007-22</a>

Note: If you need a 5 in. cage, contact technical support at [www.phenomenex.com/chat](http://www.phenomenex.com/chat) or simply reach out to your Technical consultant. Conditions may apply. Agilent 6850, some SRI and process GC systems use only 5 in. cages.

#### Column Profile



#### Phase Chemistry



100 % Polyethylene Glycol

#### Recommended Applications

- |                        |                   |
|------------------------|-------------------|
| • Alcohols             | • Glycols         |
| • Aldehydes            | • Pharmaceuticals |
| • Aromatics            | • Solvents        |
| • Basic Compounds      | • Styrene         |
| • Essential Oils       | • Xylene Isomers  |
| • Flavors & Fragrances |                   |



**ZB-WAX Test Mix**  
**Part No.:** [AGO-5158](#)



Contact Phenomenex or your local Phenomenex distributor for additional GC products and applications.



Extend column lifetime. Add a Z-Guard™ to your next Zebron GC order.

## ZB-FFAP


### Improve Resolution for Free Fatty Acids

- High polarity column; excellent thermal and chemical stability
- Provides better peak shape for underivatized acids
- Especially suited for organic acids, free fatty acids, and alcohols
- Bonded, solvent rinsable FFAP phase

Upgrade to Zebron from any nitroterephthalic acid modified polyethylene glycol phase:

Agilent®	Restek®	SGE®	Supelco®	OV®
<ul style="list-style-type: none"> <li>• DB®-FFAP</li> <li>• HP-FFAP</li> <li>• CP-Wax 58 FFAP CB</li> <li>• CP-FFAP CB</li> </ul>	<ul style="list-style-type: none"> <li>• Stabilwax®-DA</li> </ul>	<ul style="list-style-type: none"> <li>• BP21</li> </ul>	<ul style="list-style-type: none"> <li>• Nukol</li> <li>• SPB®-1000</li> </ul>	<ul style="list-style-type: none"> <li>• OV-351</li> </ul>

### Unsaturated Free Fatty Acids by GC-FID



**Column:** Zebron ZB-FFAP  
**Dimensions:** 60 meter x 0.25 mm x 0.25 µm  
**Part No.:** [7KG-G009-11](#)  
**Injection:** Split 40:1 @ 220 °C, 0.2 µL  
**Carrier Gas:** Helium @ 2.4 mL/min (constant flow)  
**Oven Program:** 200 °C to 260 °C @ 2 °C/min for 30 min  
**Detector:** FID @ 250 °C

**Sample:**


1. Myristic Acid (C14:0)	9. Linolelaidic Acid (C18:2t)
2. Myristoleic Acid (C14:1c)	10. Linoleic Acid (C18:2c)
3. Palmitic Acid (C16:0)	11. Linolenic Acid (C18:3c)
4. Palmitoleic Acid (C16:1t)	12. Arachidic Acid (C20:0)
5. Palmitoleic Acid (C16:1c)	13. Gondoic Acid (C20:1c)
6. Stearic Acid (C18:0)	14. Behenic Acid (C22:0)
7. Elaidic Acid (C18:1t)	15. Erucic Acid (C22:1c)
8. Oleic Acid (C18:1c)	

### Ordering Information


Zebron ZB-FFAP GC Columns			
ID (mm)	df (µm)	Temp. Limits °C	Part No.
<b>15-Meter</b>			
0.25	0.25	40 to 250/260	<a href="#">7EG-G009-11</a>
0.32	0.25	40 to 250/260	<a href="#">7EM-G009-11</a>
0.53	1.00	40 to 250/260	<a href="#">7EK-G009-22</a>
<b>30-Meter</b>			
0.25	0.25	40 to 250/260	<a href="#">7HG-G009-11</a>
0.32	0.25	40 to 250/260	<a href="#">7HM-G009-11</a>
0.32	0.50	40 to 250/260	<a href="#">7HM-G009-17</a>
0.32	1.00	40 to 250/260	<a href="#">7HM-G009-22</a>
0.53	1.00	40 to 250/260	<a href="#">7HK-G009-22</a>
<b>50-Meter</b>			
0.32	0.50	40 to 250/260	<a href="#">7JM-G009-17</a>
<b>60-Meter</b>			
0.25	0.25	40 to 250/260	<a href="#">7KG-G009-11</a>

Note: If you need a 5 in. cage, contact technical support at [www.phenomenex.com/chat](http://www.phenomenex.com/chat) or simply reach out to your Technical consultant. Conditions may apply. Agilent 6850, some SRI and process GC systems use only 5 in. cages.

### Column Profile



### Phase Chemistry



Nitroterephthalic Acid Modified Polyethylene Glycol

### Recommended Applications

- |                    |                       |
|--------------------|-----------------------|
| • Acrylates        | • Ketones             |
| • Alcohols         | • Organic Acids       |
| • Aldehydes        | • Phenols             |
| • Free Fatty Acids | • Volatile Free Acids |



**ZB-FFAP Test Mix**  
**Part No.:** [AGO-5158](#)



Contact Phenomenex or your local Phenomenex distributor for additional GC products and applications.



Extend column lifetime. Add a Z-Guard™ to your next Zebron GC order.

## ZB-XLB

### Get Extremely Low Bleed

- Unique, low polarity si-arylene column
- Engineered specifically for use with bleed sensitive detectors such as MS
- Provides alternate selectivity to standard 5-type phases
- Often used for confirmation of pesticides, PCBs, or other environmental samples
- Good tool for sample screening to identify unknown contaminants

Upgrade to Zebron from these similar\* phases:

#### Agilent®

- DB®-XLB
- VF-XMS

#### Restek®

- Rtx®-XLB
- Rxi®-XLB


#### Supelco®

- MDN-12

\*not exact equivalent, selectivity may differ


### Better Performance for Chlorinated Pesticides EPA Method 8081A

#### Zebron ZB-XLB



VS.

#### Restek Rxi-XLB



#### Conditions for both columns:

**Columns:** As listed  
**Dimensions:** 30 meter x 0.25 mm x 0.25 µm  
**Injection:** Split 111:1 @ 250 °C, 1.5 µL  
**Carrier Gas:** Helium @ 0.9 mL/min (constant flow)  
**Oven Program:** 110 °C to 320 °C @ 15 °C/min and hold until last peak elutes  
**Detector:** ECD @ 350 °C

**Sample:**

1. α-BHC	11. 4,4'-DDE
2. γ-BHC	12. Dieldrin
3. β-BHC	13. Endrin
4. δ-BHC	14. 4,4'-DDD
5. Heptachlor	15. Endosulfan II
6. Aldrin	16. Endrin aldehyde
7. Heptachlor epoxide	17. 4,4'-DDT
8. γ-Chlordane	18. Endosulfan sulfate
9. α-Chlordane	19. Methoxychlor
10. Endosulfan I	20. Endrin ketone

Comparative separations may not be representative of all applications.

#### Ordering Information

##### Zebron ZB-XLB GC Columns

ID(mm)	df(µm)	Temp. Limits °C	Part No.
<b>20-Meter</b>			
0.18	0.18	30 to 340/360	<a href="#">7FD-G019-08</a>


#### Ordering Information

##### Zebron ZB-XLB GC Columns (cont'd)

ID(mm)	df(µm)	Temp. Limits °C	Part No.
<b>30-Meter</b>			
0.25	0.25	30 to 340/360	<a href="#">7HG-G019-11</a>
0.25	0.50	30 to 340/360	<a href="#">7HG-G019-17</a>
0.32	0.50	30 to 340/360	<a href="#">7HM-G019-17</a>
<b>60-Meter</b>			
0.25	0.25	30 to 340/360	<a href="#">7KG-G019-11</a>

Note: If you need a 5 in. cage, contact technical support at [www.phenomenex.com/chat](http://www.phenomenex.com/chat) or simply reach out to your Technical consultant. Conditions may apply. Agilent 6850, some SRI and process GC systems use only 5 in. cages.

#### Column Profile



\*Thicker films (≥ 1.0 µm) are rated to 320/340 °C.


#### Engineered Self Cross-linking™ (ESC)

#### Phase Chemistry

- Proprietary

#### Recommended Applications

- Herbicides / Insecticides
- PCBs
- Pesticides
- Unknown Samples



Contact Phenomenex or your local Phenomenex distributor for additional GC products and applications.



Extend column lifetime. Add a Z-Guard™ to your next Zebron GC order.




## Guardian Integrated Guard Columns

### Built-in Column Protection: No Leaks, No Worries!


#### Why Choose Zebron With Guardian?

Guardian columns have the 2 m, 5 m or 10 m guard built directly into the analytical column in one continuous length of tubing. Unlike traditional guard columns, which are known to be difficult to seal and prone to leaking after normal column maintenance, the Guardian system provides the same inert column protection, but eliminates the possibility of leaks.


- Eliminate the potential for leaks
- Extend column life
- Improve analyte focusing for low boiling compounds
- Aggressively tested to ensure deactivation



#### How It Works



The sample is introduced onto the Guardian section of the column.



As temperature increases (oven ramp program), the sample is vaporized and moves unretained through the Guardian section of the column. Non-volatile contaminants are deposited on the Guardian section, better preserving the stationary phase and making it easier to trim contaminants off the front of the column.



When the analytes reach the stationary phase (analytical portion of the column), they are refocused, resulting in a narrower initial peak width. This can help improve resolution.

#### Ordering Information


##### Guardian: Integrated Guard Columns

Zebron GC Column Phase	Dimensions	2m Guardian Part No.	5m Guardian Part No.	10m Guardian Part No.
ZB-1PLUS	30 meter x 0.25 mm x 0.25 µm	—	<a href="#">7HG-G031-11-GGA</a>	<a href="#">7HG-G031-11-GGC</a>
ZB-5ms	15 meter x 0.25 mm x 0.25 µm	—	—	<a href="#">7EG-G010-11-GGC</a>
ZB-5ms	30 meter x 0.25 mm x 0.25 µm	—	<a href="#">7HG-G010-11-GGA</a>	<a href="#">7HG-G010-11-GGC</a>
ZB-5ms	30 meter x 0.32 mm x 0.25 µm	—	<a href="#">7HM-G010-11-GGA</a>	—
ZB-5ms	30 meter x 0.32 mm x 1.00 µm	—	<a href="#">7HM-G010-22-GGA</a>	—
ZB-5MSPLUS™	30 meter x 0.25 mm x 0.25 µm	—	<a href="#">7HG-G030-11-GGA</a>	<a href="#">7HG-G030-11-GGC</a>
ZB-5	30 meter x 0.25 mm x 0.25 µm	—	<a href="#">7HG-G002-11-GGA</a>	<a href="#">7HG-G002-11-GGC</a>
ZB-5	30 meter x 0.25 mm x 0.50 µm	—	—	<a href="#">7HG-G002-17-GGC</a>
ZB-5HT Inferno	30 meter x 0.25 mm x 0.10 µm	—	<a href="#">7HG-G015-02-GGA</a>	—
ZB-5HT Inferno	30 meter x 0.25 mm x 0.25 µm	—	<a href="#">7HG-G015-11-GGA</a>	—
ZB-5PLUS™	20 meter x 0.18 mm x 0.18 µm	—	<a href="#">7FD-G032-08-GGA</a>	—
ZB-5PLUS	30 meter x 0.25 mm x 0.10 µm	—	<a href="#">7HG-G032-02-GGA</a>	—
ZB-5PLUS	30 meter x 0.25 mm x 0.25 µm	—	<a href="#">7HG-G032-11-GGA</a>	—
ZB-50	10 meter x 0.18 mm x 0.18 µm	<a href="#">7CD-G004-08-GGT</a>	—	—
ZB-MultiResidue™-1	30 meter x 0.25 mm x 0.25 µm	—	—	<a href="#">7HG-G016-11-GGC</a>
ZB-SemiVolatiles	30 meter x 0.25 mm x 0.25 µm	—	<a href="#">7HG-G027-11-GGA</a>	<a href="#">7HG-G027-11-GGC</a>
ZB-Dioxin	60 meter x 0.25 mm x 0.20 µm	—	<a href="#">7KG-G045-10-GGA</a>	—

## Z-Guard™ Columns

### Protect and Extend Column Lifetime

- Individually QC tested to ensure the highest level of quality
- Extend column lifetime by preventing stationary phase damage
- Improve separation and peak shapes (especially early elutors)
- Improve sensitivity and accuracy of quantitative results
- Available as individual guard columns or as complete kits with connectors



To ensure that all Z-Guards are the highest possible quality, we individually test each one! The columns are attached to a reference Zebron ZB-5 column and are tested using our specially designed QC mix. We carefully monitor activity, bleed, and stability. This way, we are able to say with confidence that Z-Guards will provide the low activity and high quality your methods require.

#### Ordering Information

Metal Z-Guard Column			
ID (mm)	Description	Part No.	
5-Meter			
0.53	Guard Column	<a href="#">7AK-G000-00-GMO</a>	

High Temperature Z-Guard Columns and Kits			
ID (mm)	Description	Part No.	Part No.
5-Meter			
10-Meter			
0.25	Guard Column	<a href="#">7AG-G000-00-GHO</a>	<a href="#">7CG-G000-00-GHO</a>
	Guard Column Kit	<a href="#">7AG-G000-00-GHK</a>	<a href="#">7CG-G000-00-GHK</a>
0.32	Guard Column	<a href="#">7AM-G000-00-GHO</a>	<a href="#">7CM-G000-00-GHO</a>
	Guard Column Kit	<a href="#">7AM-G000-00-GHK</a>	—
0.53	Guard Column	<a href="#">7AK-G000-00-GHO</a>	<a href="#">7CK-G000-00-GHO</a>
	Guard Column Kit	<a href="#">7AK-G000-00-GHK</a>	<a href="#">7CK-G000-00-GHK</a>

Standard Z-Guard Columns and Kits			
ID (mm)	Description	Part No.	Part No.
5-Meter			
10-Meter			
0.10	Guard Column	<a href="#">7AB-G000-00-GZO</a>	<a href="#">7CB-G000-00-GZO</a>
	Guard Column Kit	<a href="#">7AD-G000-00-GZK</a>	<a href="#">7CD-G000-00-GZK</a>
0.18	Guard Column	<a href="#">7AD-G000-00-GZO</a>	<a href="#">7CD-G000-00-GZO</a>
	Guard Column Kit	<a href="#">7AD-G000-00-GZK</a>	—
0.25	Guard Column	<a href="#">7AG-G000-00-GZO</a>	<a href="#">7CG-G000-00-GZO</a>
	Guard Column Kit	<a href="#">7AG-G000-00-GZK</a>	<a href="#">7CG-G000-00-GZK</a>
0.32	Guard Column	<a href="#">7AM-G000-00-GZO</a>	<a href="#">7CM-G000-00-GZO</a>
	Guard Column Kit	<a href="#">7AM-G000-00-GZK</a>	<a href="#">7CM-G000-00-GZK</a>
0.53	Guard Column	<a href="#">7AK-G000-00-GZO</a>	<a href="#">7CK-G000-00-GZO</a>
	Guard Column Kit	<a href="#">7AK-G000-00-GZK</a>	<a href="#">7CK-G000-00-GZK</a>

Bulk Z-Guard Columns			
ID (mm)	Description	Part No.	Unit
50-Meter			
0.25	Guard Column	<a href="#">7JG-G000-00-GZO</a>	ea
0.32	Guard Column	<a href="#">7JM-G000-00-GZO</a>	ea
0.53	Guard Column	<a href="#">7JK-G000-00-GZO</a>	ea
5-Meter			
0.53	Guard Column	<a href="#">7AK-G000-00-GZ1</a>	10/pk

ZB-5 Z-Guard Column Multi-Pak			
ID (mm)	Description	Part No.	Unit
2-Meter			
0.25	Zebron ZB-5 Z-Guard Column	<a href="#">KG0-7868</a>	25/pk

**i** Universal GC Guard Column. Designed for use with virtually any GC capillary column from virtually any manufacturer. Alternative to: Restek, Supelco, Agilent Technologies, and many more.

**i** Z-Guard Column Kits include 5 or 10 meters of deactivated fused silica tubing, 5 universal connectors and 0.5 mL of high-temperature polyimide resin.

Replacement Parts for Z-Guard Kits		
Description	Part No.	Unit
Universal Capillary Column Union, Borosilicate	<a href="#">AG0-4716</a>	5/pk
High Temperature Polyimide Resin, 0.5 mL	<a href="#">AG0-8514</a>	ea